




## TELEFLEX AIRWAYS

A solid tradition in innovation



## TELEFLEX AIRWAY MANAGEMENT – A SOLID TRADITION IN INNOVATION

Trusted brands make Teleflex a reliable and strong partner. Built on a solid tradition of innovation, Teleflex is a global leader in superior medical supplies designed to help providers minimise risk and maximise outcomes for their patients. Our understanding of the importance to our customers of a full range of products has led to the development of a unique line of medical devices, all of which complement one another.

Every article supplied by the Teleflex brands RÜSCH, HUDSON RCI, GIBECK and SHERIDAN is the result of many years of manufacturing experience, specialisation, as well as cooperation with well-known medical experts. Our primary concern is to offer you products of constant quality & unparalleled performance in the field of airway management, to help you focus on what's most important: your patient.

You can find further details and the technical specifications of our products within this brochure.

### TELEFLEX – HIGH QUALITY MEDICAL SUPPLIES FROM A SINGLE SOURCE



#### PAEDIATRIC SIZES AVAILABLE

To meet the special requirements of paediatric care, Teleflex supplies a large variety of tracheal and bronchial tubes. Frequently used models also come in smaller sizes.




#### SINGLE USE PRODUCT



#### DEHP-FREE

To increase patient safety, Teleflex supplies a range of products that does not contain DEHP.



|                                                     |    |
|-----------------------------------------------------|----|
| <b>ENDOTRACHEAL TUBES</b> .....                     | 2  |
| RÜSCH endotracheal tubes.....                       | 2  |
| SHERIDAN endotracheal tubes.....                    | 19 |
| <b>LARYNGEAL MASKS</b> .....                        | 25 |
| <b>BRONCHIAL TUBES AND BLOCKERS</b> .....           | 29 |
| RÜSCH bronchial tubes and blockers .....            | 29 |
| SHERIDAN bronchial tubes .....                      | 36 |
| <b>OROPHARYNGEAL / NASOPHARYNGEAL AIRWAYS</b> ..... | 37 |
| Oropharyngeal Airways.....                          | 37 |
| Nasopharyngeal Airways.....                         | 39 |
| <b>ANAESTHESIA MASKS / REBREATHING BAGS</b> .....   | 41 |
| Reusable face masks .....                           | 41 |
| Disposable face masks .....                         | 43 |
| Resuscitators .....                                 | 44 |
| Rebreathing bags .....                              | 45 |
| <b>STYLETS</b> .....                                | 46 |
| <b>ACCESSORIES</b> .....                            | 48 |
| <b>INDEX</b> .....                                  | 50 |

# RÜSCH ENDOTRACHEAL TUBES

The optimum tube for all applications. At RÜSCH, this demand is met thanks to an extraordinarily extensive and versatile range of tubes. All RÜSCH tubes feature an outstanding I.D. to O.D. ratio. Our high quality standards are reflected in the high-grade materials we use, which are tested according to the highest medical standards.

To meet the special requirements of paediatric patients, Teleflex's RÜSCH brand supplies a large variety of tracheal tubes in smaller sizes. Tracheal tubes for paediatric care are marked with black tip to ensure safe tracheal positioning.

| CHILD AGE | TUBE I.D.  | BLACK TIP LENGTH |
|-----------|------------|------------------|
| 1.        | 2.0–3.5 mm | 20 mm            |
| 2.–5.     | 4.0–5.0 mm | 30/40 mm         |
| 6.–14.    | 5.5–6.5 mm | 40 mm            |

## TRACHEAL TUBES:

1. Before using tracheal tubes with the connector pushed on halfway, it must be ensured that the connector is pushed into the tube shaft as far as it will go.
2. To improve the tight hold, it is advisable to clean with ethanol, both the contact surface of a connector that has been loosened once or removed and the corresponding tube shaft, before reinserting it.
3. This particularly applies to those cases in which the connector was completely removed first (e.g. for fiberoptic insertions), or was pushed on again after the tube had been shortened.

## SUPER SAFETYCLEAR



### CLEAR TRACHEAL TUBE MADE OF PVC, WITH LOW-PRESSURE CUFF

nasal/oral

sizes: I.D. 2.5–10.0 mm

semi-seated connector, valve for Luer and Luer-lock syringes, cupped atraumatic tip, continuous X-ray marker, black position indicator for correct tube placement, blue pilot balloon, graduated

- latex-free
- sterile



| SUPER SAFETYCLEAR     |                 |         |         |         | RÜSCH |
|-----------------------|-----------------|---------|---------|---------|-------|
| REF.                  | ORDER SIZE/I.D. | O.D.    | CUFF Ø  | LENGTH* | QTY   |
| 112480, Magill        | 2.5 mm          | 4.0 mm  | 8.0 mm  | 160 mm  | 10    |
|                       | 3.0 mm          | 5.0 mm  | 8.0 mm  | 170 mm  |       |
|                       | 3.5 mm          | 5.3 mm  | 8.0 mm  | 190 mm  |       |
|                       | 4.0 mm          | 6.0 mm  | 10.5 mm | 220 mm  |       |
|                       | 4.5 mm          | 6.3 mm  | 10.5 mm | 230 mm  |       |
| 112482, Murphy, 1 eye | 5.0 mm          | 6.7 mm  | 13.0 mm | 250 mm  |       |
|                       | 5.5 mm          | 7.3 mm  | 16.5 mm | 280 mm  |       |
|                       | 6.0 mm          | 8.0 mm  | 18.5 mm | 290 mm  |       |
|                       | 6.5 mm          | 8.7 mm  | 20.5 mm | 300 mm  |       |
|                       | 7.0 mm          | 9.3 mm  | 24.0 mm | 320 mm  |       |
|                       | 7.5 mm          | 10.0 mm | 26.0 mm | 330 mm  |       |
|                       | 8.0 mm          | 10.7 mm | 26.0 mm | 340 mm  |       |
|                       | 8.5 mm          | 11.3 mm | 28.0 mm | 345 mm  |       |
|                       | 9.0 mm          | 12.0 mm | 28.0 mm | 350 mm  |       |
|                       | 9.5 mm          | 12.7 mm | 29.0 mm | 350 mm  |       |
| 10.0 mm               | 13.3 mm         | 29.0 mm | 350 mm  |         |       |

\* Length without connector

## FLEXISET



### SUPER SAFETYCLEAR TRACHEAL TUBE WITH MURPHY EYE AND FLEXISLIP INTUBATION STYLET

sizes: I.D. 5.0–9.0 mm

To order each product separately

- Super SafetyClear tracheal tube only, Ref. 112482
- FlexiSlip intubation stylet only, Ref. 503700  
(please refer to page 46)

#### FLEXISET

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.    | CUFF Ø  | LENGTH* | QTY |
|--------|-----------------|---------|---------|---------|-----|
| 112481 | 5.0 mm          | 6.7 mm  | 13.0 mm | 250 mm  | 10  |
|        | 5.5 mm          | 7.3 mm  | 16.5 mm | 280 mm  |     |
|        | 6.0 mm          | 8.0 mm  | 18.5 mm | 290 mm  |     |
|        | 6.5 mm          | 8.7 mm  | 20.5 mm | 300 mm  |     |
|        | 7.0 mm          | 9.3 mm  | 24.0 mm | 320 mm  |     |
|        | 7.5 mm          | 10.0 mm | 26.0 mm | 330 mm  |     |
|        | 8.0 mm          | 10.7 mm | 26.0 mm | 340 mm  |     |
|        | 8.5 mm          | 11.3 mm | 28.0 mm | 345 mm  |     |
| 9.0 mm | 12.0 mm         | 28.0 mm | 350 mm  |         |     |

\* Length without connector

## SAFETYCLEAR PLUS



### CLEAR TRACHEAL TUBE MADE OF PVC WITH HIGH-VOLUME LOW-PRESSURE SUPER PLUS CUFF nasal/oral

sizes: I.D. 5.0–10.0 mm

semi-seated connector, valve for Luer and Luer-lock syringes, cupped atraumatic tip, continuous X-ray marker, black position indicator for correct tube placement, blue pilot balloon, graduated

- latex-free
- sterile



#### SAFETYCLEAR PLUS

**RÜSCH**

| REF.                  | ORDER SIZE/I.D. | O.D.    | CUFF Ø | LENGTH* | QTY |
|-----------------------|-----------------|---------|--------|---------|-----|
| 112080, Magill        | 5.0 mm          | 6.7 mm  | 20 mm  | 250 mm  | 10  |
|                       | 5.5 mm          | 7.3 mm  | 20 mm  | 280 mm  |     |
|                       | 6.0 mm          | 8.0 mm  | 24 mm  | 290 mm  |     |
|                       | 6.5 mm          | 8.7 mm  | 24 mm  | 300 mm  |     |
|                       | 7.0 mm          | 9.3 mm  | 28 mm  | 320 mm  |     |
|                       | 7.5 mm          | 10.0 mm | 28 mm  | 330 mm  |     |
|                       | 8.0 mm          | 10.7 mm | 30 mm  | 340 mm  |     |
|                       | 8.5 mm          | 11.3 mm | 30 mm  | 345 mm  |     |
|                       | 9.0 mm          | 12.0 mm | 34 mm  | 350 mm  |     |
|                       | 9.5 mm          | 12.7 mm | 34 mm  | 350 mm  |     |
| 112082, Murphy, 1 eye | 9.0 mm          | 12.0 mm | 34 mm  | 350 mm  |     |
|                       | 10.0 mm         | 13.3 mm | 34 mm  | 350 mm  |     |

\* Length without connector

# SAFETYCLEAR




## CLEAR TRACHEAL TUBE MADE OF PVC WITHOUT CUFF

nasal/oral

sizes: I.D. 2.0–10.0 mm

semi-seated connector, cupped atraumatic tip, continuous X-ray marker, graduated

- latex-free
- sterile



### SAFETYCLEAR

**RÜSCH**

| REF.                                               | ORDER SIZE/I.D. | O.D.    | LENGTH* | LENGTH BLACK TIP | QTY |
|----------------------------------------------------|-----------------|---------|---------|------------------|-----|
| 100380, Magill<br>sizes:<br>I.D. 2.0–10.0 mm       | 2.0 mm          | 3.0 mm  | 150 mm  | 20 mm            | 10  |
|                                                    | 2.5 mm          | 3.3 mm  | 160 mm  | 20 mm            |     |
|                                                    | 3.0 mm          | 4.0 mm  | 170 mm  | 20 mm            |     |
|                                                    | 3.5 mm          | 4.7 mm  | 190 mm  | 20 mm            |     |
|                                                    | 4.0 mm          | 5.3 mm  | 210 mm  | 30 mm            |     |
|                                                    | 4.5 mm          | 6.0 mm  | 230 mm  | 30 mm            |     |
| 100382, Murphy, 1 eye<br>sizes:<br>I.D. 2.0–7.0 mm | 5.0 mm          | 6.7 mm  | 250 mm  | 40 mm            |     |
|                                                    | 5.5 mm          | 7.3 mm  | 280 mm  | 40 mm            |     |
|                                                    | 6.0 mm          | 8.0 mm  | 290 mm  | 40 mm            |     |
|                                                    | 6.5 mm          | 8.7 mm  | 300 mm  | 40 mm            |     |
|                                                    | 7.0 mm          | 9.3 mm  | 320 mm  | –                |     |
|                                                    | 7.5 mm          | 10.0 mm | 330 mm  | –                |     |
|                                                    | 8.0 mm          | 10.7 mm | 340 mm  | –                |     |
|                                                    | 8.5 mm          | 11.3 mm | 345 mm  | –                |     |
|                                                    | 9.0 mm          | 12.0 mm | 350 mm  | –                |     |
|                                                    | 9.5 mm          | 12.7 mm | 350 mm  | –                |     |
|                                                    | 10.0 mm         | 13.3 mm | 350 mm  | –                |     |

\* Length without connector

# ISIS ETT

## ENDOTRACHEAL TUBE WITH INTEGRATED PORT FOR OPTIONAL SUBGLOTTIC SECRETION SUCTIONING

The ideal tube to provide the best of care to all your patients, up to 30 days.

nasal/oral 1 eye (Murphy)

Sizes: I.D. 6.0–9.0 mm

This convertible endotracheal tube with separate suction line allows for subglottic secretion removal on demand. This frees the clinician from the burden of deciding for an expensive tube at the time of intubation or to extubate/re-intubate a patient to provide a VAP prevention method of subglottic secretion suctioning. Thanks to the ISIS flexible design, subglottic suctioning can be provided to any patient whenever necessary.

## SUCTION LINE FOR ISIS ENDOTRACHEAL TUBE

When connected to the ISIS endotracheal tube, the optional suction line provides access from the top of the ETT cuff for easy & safe suctioning of the subglottic secretions that spill over the cuff. Available separately.

1 size fits all ISIS tubes




### TELEFLEX ISIS ETT

| REF.   | ORDER SIZE/I.D. | O.D.    | CUFF Ø  | LENGTH* | QTY |
|--------|-----------------|---------|---------|---------|-----|
| 112662 | 6.0 mm          | 9.2 mm  | 20.5 mm | 290 mm  | 10  |
|        | 6.5 mm          | 9.9 mm  | 24 mm   | 300 mm  |     |
|        | 7.0 mm          | 10.6 mm | 26 mm   | 320 mm  |     |
|        | 7.5 mm          | 11.3 mm | 26 mm   | 330 mm  |     |
|        | 8.0 mm          | 12.0 mm | 28 mm   | 340 mm  |     |
|        | 8.5 mm          | 12.6 mm | 28 mm   | 345 mm  |     |
|        | 9.0 mm          | 13.2 mm | 29 mm   | 350 mm  |     |

\* Length without connector


Endotracheal tube with integrated port for optional subglottic secretion suctioning.



### SUCTION LINE

| REF.   | QTY |
|--------|-----|
| 502700 | 5   |

# SUPER SAFETY SILK



**TRACHEAL TUBE  
MADE OF DEHP-FREE PVC  
WITH LOW PRESSURE  
HIGH VOLUME CUFF**

nasal/oral  
sizes: I.D. 3.0–10.0 mm

tissue-friendly thermosensitive material, especially advantageous for long-term intubation, soft atraumatic tip, continuous X-ray marker, semi-seated connector, low pressure sealing, Safety Silk surface, black ring for intubation depth control

- latex-free
- sterile



## SUPER SAFETY SILK

**RÜSCH**

| REF.                                                | ORDER SIZE/I.D. | O.D.    | CUFF Ø  | LENGTH | QTY |
|-----------------------------------------------------|-----------------|---------|---------|--------|-----|
| 112682, Murphy, 1 eye<br>sizes:<br>I.D. 3.0–10.0 mm | 3.0 mm          | 5.3 mm  | 8.0 mm  | 170 mm | 10  |
|                                                     | 3.5 mm          | 6.0 mm  | 10.5 mm | 190 mm |     |
|                                                     | 4.0 mm          | 6.3 mm  | 10.5 mm | 220 mm |     |
|                                                     | 4.5 mm          | 6.7 mm  | 13.0 mm | 230 mm |     |
| 112680, Magill<br>sizes:<br>I.D. 5.0–10.0 mm        | 5.0 mm          | 7.3 mm  | 16.5 mm | 250 mm |     |
|                                                     | 5.5 mm          | 8.0 mm  | 18.5 mm | 280 mm |     |
|                                                     | 6.0 mm          | 8.7 mm  | 20.5 mm | 290 mm |     |
|                                                     | 6.5 mm          | 9.3 mm  | 24.0 mm | 300 mm |     |
|                                                     | 7.0 mm          | 10.0 mm | 26.0 mm | 320 mm |     |
|                                                     | 7.5 mm          | 10.7 mm | 26.0 mm | 330 mm |     |
| 8.0 mm                                              | 11.6 mm         | 28.0 mm | 340 mm  |        |     |
| 8.5 mm                                              | 12.4 mm         | 28.0 mm | 345 mm  |        |     |
| 9.0 mm                                              | 13.0 mm         | 28.0 mm | 350 mm  |        |     |
| 9.5 mm                                              | 14.2 mm         | 29.0 mm | 350 mm  |        |     |
| 10.0 mm                                             | 15.0 mm         | 29.0 mm | 350 mm  |        |     |



## PEDI SAFETY SILK



### CLEAR TRACHEAL TUBE MADE OF SOFT PVC

part of the Safety Silk family,  
nasal/oral (Magill)

sizes: I.D. 2.0–5.5 mm

semi-seated connector, cupped atrau-  
matic black tip, continuous X-ray  
marker, soft version, graduated

- latex-free
- sterile




### PEDI SAFETY SILK

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.   | LENGTH* | BLACK TIP LENGTH | QTY |
|--------|-----------------|--------|---------|------------------|-----|
| 100480 | 2.0 mm          | 3.3 mm | 140 mm  | 20 mm            | 10  |
|        | 2.5 mm          | 4.0 mm | 160 mm  | 20 mm            |     |
|        | 3.0 mm          | 4.7 mm | 180 mm  | 20 mm            |     |
|        | 3.5 mm          | 5.3 mm | 200 mm  | 20 mm            |     |
|        | 4.0 mm          | 6.0 mm | 220 mm  | 30 mm            |     |
|        | 4.5 mm          | 6.7 mm | 240 mm  | 30 mm            |     |
|        | 5.0 mm          | 7.3 mm | 260 mm  | 40 mm            |     |
|        | 5.5 mm          | 8.0 mm | 280 mm  | 40 mm            |     |

\* Length without connector

## NASAL SAFETY SILK



### CLEAR SOFT ANATOMICALLY SHAPED TUBE

MADE OF DEHP-FREE PVC

nasal only, 1 eye (Murphy)

sizes: I.D. 3.0–8.0 mm

extremely soft & tissue-friendly PVC,  
bonded connector, kink-resistant,  
cupped atraumatic tip, continuous blue  
X-ray marker, graduated, with low-  
pressure cuff, blue pilot balloon, valve  
for Luer and Luer-lock syringes

- latex-free
- sterile




### NASAL SAFETY SILK

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.    | CUFF Ø  | LENGTH | QTY |
|--------|-----------------|---------|---------|--------|-----|
| 111782 | 3.0 mm          | 5.3 mm  | 11.0 mm | 240 mm | 10  |
|        | 3.5 mm          | 6.0 mm  | 11.0 mm | 271 mm |     |
|        | 4.0 mm          | 6.3 mm  | 13.0 mm | 303 mm |     |
|        | 4.5 mm          | 6.7 mm  | 13.0 mm | 338 mm |     |
|        | 5.0 mm          | 7.3 mm  | 16.0 mm | 358 mm |     |
|        | 5.5 mm          | 8.0 mm  | 18.0 mm | 368 mm |     |
|        | 6.0 mm          | 8.7 mm  | 20.0 mm | 415 mm |     |
|        | 6.5 mm          | 9.3 mm  | 24.0 mm | 420 mm |     |
|        | 7.0 mm          | 10.0 mm | 26.0 mm | 430 mm |     |
|        | 7.5 mm          | 10.7 mm | 26.0 mm | 445 mm |     |
| 8.0 mm | 11.3 mm         | 28.0 mm | 455 mm  |        |     |

# AGT TRACHEAL TUBE ORAL



**CLEAR ANATOMICALLY SHAPED  
AGT TUBE MADE OF PVC  
WITH LOW-PRESSURE CUFF**


oral, 1 eye  
sizes: I.D. 3.5–9.0 mm

semi-seated connector, cupped  
atraumatic tip, continuous X-ray  
marker, graduated, low-pressure cuff,  
blue pilot balloon, valve for Luer  
and Luer-lock syringes

- latex-free
- sterile

**AGT TRACHEAL TUBE, ORAL** **RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.    | CUFF Ø | LENGTH* | QTY |
|--------|-----------------|---------|--------|---------|-----|
| 111780 | 3.5 mm          | 5.3 mm  | 11 mm  | 120 mm  | 10  |
|        | 4.0 mm          | 6.0 mm  | 11 mm  | 130 mm  |     |
|        | 4.5 mm          | 6.3 mm  | 13 mm  | 145 mm  |     |
|        | 5.0 mm          | 6.7 mm  | 13 mm  | 155 mm  |     |
|        | 5.5 mm          | 7.3 mm  | 16 mm  | 165 mm  |     |
|        | 6.0 mm          | 8.0 mm  | 18 mm  | 180 mm  |     |
|        | 6.5 mm          | 8.7 mm  | 20 mm  | 190 mm  |     |
|        | 7.0 mm          | 9.3 mm  | 24 mm  | 200 mm  |     |
|        | 7.5 mm          | 10.0 mm | 26 mm  | 210 mm  |     |
|        | 8.0 mm          | 10.7 mm | 26 mm  | 220 mm  |     |
|        | 8.5 mm          | 11.3 mm | 28 mm  | 230 mm  |     |
| 9.0 mm | 12.0 mm         | 28 mm   | 240 mm |         |     |



**CLEAR ANATOMICALLY SHAPED  
AGT TUBE MADE OF PVC  
WITHOUT CUFF**

oral, 2 eyes  
sizes: I.D. 3.0–7.0 mm

semi-seated connector, cupped  
atraumatic tip, continuous X-ray  
marker, graduated

- latex-free
- sterile

**AGT TRACHEAL TUBE, ORAL** **RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.   | O.D. (CH.) | LENGTH* | RECOMMENDED SIZE FOR | QTY |
|--------|-----------------|--------|------------|---------|----------------------|-----|
| 100180 | 3.0 mm          | 4.0 mm | Ch. 12     | 105 mm  | newborn babies       | 10  |
|        | 3.5 mm          | 4.7 mm | Ch. 14     | 120 mm  | 1–6 months           |     |
|        | 4.0 mm          | 5.3 mm | Ch. 16     | 130 mm  | 0.5–1 year           |     |
|        | 4.5 mm          | 6.0 mm | Ch. 18     | 145 mm  | 1–2 years            |     |
|        | 5.0 mm          | 6.7 mm | Ch. 20     | 155 mm  | 2–4 years            |     |
|        | 5.5 mm          | 7.3 mm | Ch. 22     | 165 mm  | 4–6 years            |     |
|        | 6.0 mm          | 8.0 mm | Ch. 24     | 180 mm  | 6–8 years            |     |
|        | 6.5 mm          | 8.7 mm | Ch. 26     | 190 mm  | 8–9 years            |     |
| 7.0 mm | 9.3 mm          | Ch. 28 | 200 mm     | over 9  |                      |     |

\* For length see picture


## AGT TRACHEAL TUBE NASAL

### CLEAR ANATOMICALLY SHAPED AGT TUBE MADE OF PVC WITH LOW-PRESSURE CUFF

nasal only, 1 eye  
sizes: I.D. 3.5–8.0 mm

semi-seated connector, cupped atraumatic tip, continuous X-ray marker, graduated, low-pressure cuff, blue pilot balloon, valve for Luer and Luer-lock syringes


- latex-free
- sterile



### AGT TRACHEAL TUBE, NASAL

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.    | CUFF Ø | LENGTH* | QTY |
|--------|-----------------|---------|--------|---------|-----|
| 111781 | 3.5 mm          | 5.3 mm  | 11 mm  | 160 mm  | 10  |
|        | 4.0 mm          | 6.0 mm  | 11 mm  | 185 mm  |     |
|        | 4.5 mm          | 6.3 mm  | 13 mm  | 220 mm  |     |
|        | 5.0 mm          | 6.7 mm  | 13 mm  | 235 mm  |     |
|        | 5.5 mm          | 7.3 mm  | 16 mm  | 245 mm  |     |
|        | 6.0 mm          | 8.0 mm  | 18 mm  | 260 mm  |     |
|        | 6.5 mm          | 8.7 mm  | 20 mm  | 270 mm  |     |
|        | 7.0 mm          | 9.3 mm  | 24 mm  | 280 mm  |     |
|        | 7.5 mm          | 10.0 mm | 26 mm  | 290 mm  |     |
| 8.0 mm | 10.7 mm         | 26 mm   | 300 mm |         |     |



### CLEAR ANATOMICALLY SHAPED AGT TUBE MADE OF PVC WITHOUT CUFF

nasal only, 2 eyes  
sizes: I.D. 3.0–7.0 mm

semi-seated connector, cupped atraumatic tip, continuous X-ray marker, graduated

- latex-free
- sterile

### AGT TRACHEAL TUBE, NASAL

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.   | O.D. (CH.) | LENGTH* | RECOMMENDED SIZE FOR | QTY |
|--------|-----------------|--------|------------|---------|----------------------|-----|
| 100181 | 3.0 mm          | 4.0 mm | Ch. 12     | 115 mm  | newborn babies       | 10  |
|        | 3.5 mm          | 4.7 mm | Ch. 14     | 140 mm  | 1–6 months           |     |
|        | 4.0 mm          | 5.3 mm | Ch. 16     | 155 mm  | 0.5–1 year           |     |
|        | 4.5 mm          | 6.0 mm | Ch. 18     | 180 mm  | 1–2 years            |     |
|        | 5.0 mm          | 6.7 mm | Ch. 20     | 190 mm  | 2–4 years            |     |
|        | 5.5 mm          | 7.3 mm | Ch. 22     | 195 mm  | 4–6 years            |     |
|        | 6.0 mm          | 8.0 mm | Ch. 24     | 205 mm  | 6–8 years            |     |
|        | 6.5 mm          | 8.7 mm | Ch. 26     | 210 mm  | 8–9 years            |     |
| 7.0 mm | 9.3 mm          | Ch. 28 | 215 mm     | over 9  |                      |     |

\* For length see picture

# SILKOCLEAR



## CLEAR TRACHEAL TUBE MADE OF SILICONE WITH PREFORMED CUFF

nasal/oral, 1 eye (Murphy)  
sizes: I.D. 5.0–9.0 mm

fixed connector, valve for Luer and Luer-lock adapter syringes, cupped atraumatic tip, continuous X-ray marker, black position indicator for correct tube placement, transparent pilot balloon, fluted inner lumen, graduated

- latex-free
- sterile

### SILKOCLEAR

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.    | CUFF Ø  | LENGTH* | QTY |
|--------|-----------------|---------|---------|---------|-----|
| 105102 | 5.0 mm          | 7.3 mm  | 8.5 mm  | 240 mm  | 5   |
|        | 5.5 mm          | 8.0 mm  | 12.0 mm | 270 mm  |     |
|        | 6.0 mm          | 8.7 mm  | 12.0 mm | 280 mm  |     |
|        | 6.5 mm          | 9.3 mm  | 13.5 mm | 290 mm  |     |
|        | 7.0 mm          | 10.0 mm | 13.5 mm | 300 mm  |     |
|        | 7.5 mm          | 10.7 mm | 13.5 mm | 310 mm  |     |
|        | 8.0 mm          | 11.3 mm | 16.0 mm | 320 mm  |     |
|        | 8.5 mm          | 12.0 mm | 16.0 mm | 330 mm  |     |
|        | 9.0 mm          | 12.7 mm | 16.0 mm | 340 mm  |     |



## CLEAR TRACHEAL TUBE MADE OF SILICONE WITHOUT CUFF

nasal/oral (Magill)  
sizes: I.D. 2.0–4.5 mm

semi-seated connector, black tip,  
graduated

- latex-free
- sterile

### SILKOCLEAR

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.   | LENGTH* | BLACK TIP LENGTH | QTY |
|--------|-----------------|--------|---------|------------------|-----|
| 100280 | 2.0 mm          | 3.3 mm | 120 mm  | 20 mm            | 10  |
|        | 2.5 mm          | 4.0 mm | 140 mm  | 20 mm            |     |
|        | 3.0 mm          | 4.7 mm | 160 mm  | 20 mm            |     |
|        | 3.5 mm          | 5.3 mm | 180 mm  | 20 mm            |     |
|        | 4.0 mm          | 6.0 mm | 200 mm  | 30 mm            |     |
|        | 4.5 mm          | 6.7 mm | 220 mm  | 30 mm            |     |

\* Length without connector

**Please note:**

SILKOSPRAY is not suitable for devices made of silicone. If you require a lubricant, we recommend our TrachJell (page 49) or commercially available gels.

# RÜSCHFLEX

## ARMOURED PREFORMED TRACHEAL TUBE, CURVED, MADE OF PVC WITH LOW-PRESSURE CUFF

nasal/oral

sizes: I.D. 3.5–10.0 mm

fixed white connector, blue pilot balloon,  
valve for Luer- and Luer-lock syringes,  
radiopaque, graduated, also available with  
integrated FlexiSlip stylet (can be ordered  
separately, Ref. 503700. see page 46 for  
more details)

- latex-free
- sterile



*Magill, Ref. 104201  
Size: I.D. 3.5–10.0 mm*



*Magill, with FlexiSlip,  
Ref. 104203  
Size: I.D. 5.0–9.0 mm*



## RÜSCHFLEX

**RÜSCH**

| REF.                                                    | ORDER SIZE/I.D.                                         | O.D.    | LENGTH* MAGILL<br>104201/104203 | LENGTH* MURPHY<br>104202/104204 | QTY |        |
|---------------------------------------------------------|---------------------------------------------------------|---------|---------------------------------|---------------------------------|-----|--------|
| 104201, Magill<br>sizes: I.D. 3.5–10.0 mm               | 3.5 mm                                                  | 5.8 mm  | 228 mm                          | 229 mm                          | 5   |        |
|                                                         | 4.0 mm                                                  | 6.3 mm  | 235 mm                          | 237 mm                          |     |        |
| 104203, Magill with FlexiSlip<br>sizes: I.D. 5.0–9.0 mm | 4.5 mm                                                  | 6.8 mm  | 256 mm                          | 259 mm                          |     |        |
|                                                         | 5.0 mm                                                  | 7.3 mm  | 319 mm                          | 321 mm                          |     |        |
|                                                         | 5.5 mm                                                  | 7.8 mm  | 319 mm                          | 322 mm                          |     |        |
|                                                         | 6.0 mm                                                  | 8.4 mm  | 326 mm                          | 329 mm                          |     |        |
|                                                         | 6.5 mm                                                  | 9.0 mm  | 330 mm                          | 334 mm                          |     |        |
|                                                         | 7.0 mm                                                  | 9.6 mm  | 349 mm                          | 352 mm                          |     |        |
|                                                         | 7.5 mm                                                  | 10.2 mm | 353 mm                          | 356 mm                          |     |        |
|                                                         | 8.0 mm                                                  | 10.8 mm | 360 mm                          | 364 mm                          |     |        |
| 104202, Murphy, 1 eye<br>sizes: I.D. 3.5–10.0 mm        | 8.5 mm                                                  | 11.3 mm | 384 mm                          | 388 mm                          |     |        |
|                                                         | 9.0 mm                                                  | 11.9 mm | 393 mm                          | 398 mm                          |     |        |
|                                                         | 104204, Murphy with FlexiSlip<br>sizes: I.D. 5.0–9.0 mm | 9.5 mm  | 12.4 mm                         | 394 mm                          |     | 399 mm |
|                                                         |                                                         | 10.0 mm | 12.9 mm                         | 395 mm                          |     | 399 mm |

\* Length without connector

# ARMoured TRACHEAL TUBES

Ref. 104004



**ARMoured TRACHEAL TUBE, STRAIGHT, MADE OF PVC, WITH LOW-PRESSURE CUFF**  
nasal/oral  
sizes: I.D. 5.0–9.0 mm

valve for Luer and Luer-lock syringes, radiopaque, blue pilot balloon, graduated

- latex-free
- sterile



## ARMoured TRACHEAL TUBE

**RÜSCH**

| REF.                               | ORDER SIZE/I.D. | O.D.    | CUFF Ø  | LENGTH*  | QTY |
|------------------------------------|-----------------|---------|---------|----------|-----|
| 104004<br>fixed white<br>connector | 5.0 mm          | 7.8 mm  | 17.5 mm | 253.0 mm | 2   |
|                                    | 5.5 mm          | 8.3 mm  | 22.0 mm | 253.5 mm |     |
|                                    | 6.0 mm          | 8.8 mm  | 22.0 mm | 313.5 mm |     |
|                                    | 6.5 mm          | 9.3 mm  | 26.5 mm | 313.5 mm |     |
|                                    | 7.0 mm          | 9.8 mm  | 26.5 mm | 314.0 mm |     |
|                                    | 7.5 mm          | 10.3 mm | 28.5 mm | 314.0 mm |     |
|                                    | 8.0 mm          | 10.8 mm | 28.5 mm | 314.0 mm |     |
|                                    | 8.5 mm          | 11.3 mm | 29.5 mm | 354.0 mm |     |
|                                    | 9.0 mm          | 11.8 mm | 29.5 mm | 354.5 mm |     |



104000  
blue mounted  
connector  
**FO intubation**

**Tip:**  
Ref. 104000 without connector can be used for fiberoptic intubation along with a Mainz universal adaptor (see page 48).

**Attention:**  
We recommend spraying the inside of the tracheal tube and the outside of the suction catheter with SILKOSPRAY (Ref. 556000, see page 49) prior to use.

**PAEDIATRIC ARMoured TRACHEAL TUBE, STRAIGHT, MADE OF PVC, WITHOUT CUFF**  
nasal/oral  
sizes: I.D. 2.5–6.5 mm

mounted connector, radiopaque, graduated

- latex-free
- sterile



## PAEDIATRIC ARMoured TRACHEAL TUBE

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.   | LENGTH*  | QTY |
|--------|-----------------|--------|----------|-----|
| 103600 | 2.5 mm          | 4.1 mm | 184.0 mm | 2   |
|        | 3.0 mm          | 4.8 mm | 184.5 mm |     |
|        | 3.5 mm          | 5.5 mm | 184.5 mm |     |
|        | 4.0 mm          | 6.0 mm | 194.5 mm |     |
|        | 4.5 mm          | 6.5 mm | 198.0 mm |     |
|        | 5.0 mm          | 7.2 mm | 208.0 mm |     |
|        | 5.5 mm          | 7.7 mm | 208.0 mm |     |
|        | 6.0 mm          | 8.2 mm | 208.0 mm |     |
|        | 6.5 mm          | 8.9 mm | 268.5 mm |     |



\* Length without connector

## SILKOCLEAR FLEX



**ARMoured TRACHEAL TUBE  
MADE OF SILICONE, WITH HIGHLY  
FLEXIBLE THIN-WALLED CUFF**  
nasal/oral, 1 eye (Murphy)  
sizes: I.D. 3.5–10.0 mm

fixed connector, valve for Luer and  
Luer-lock syringes, black position  
indicator for correct tube placement,  
pilot balloon, smooth inner tube  
coating, graduated

- latex-free
- sterile

### SILKOCLEAR FLEX

**RÜSCH**

| REF.    | ORDER SIZE/I.D. | O.D.    | CUFF Ø   | LENGTH*  | QTY |
|---------|-----------------|---------|----------|----------|-----|
| 105702  | 3.5 mm          | 5.0 mm  | 7.0 mm   | 226.0 mm | 2   |
|         | 4.0 mm          | 5.7 mm  | 7.0 mm   | 246.0 mm |     |
|         | 4.5 mm          | 6.3 mm  | 7.0 mm   | 250.5 mm |     |
|         | 5.0 mm          | 6.7 mm  | 8.5 mm   | 268.0 mm |     |
|         | 5.5 mm          | 7.3 mm  | 12.0 mm  | 299.0 mm |     |
|         | 6.0 mm          | 8.0 mm  | 12.0 mm  | 330.5 mm |     |
|         | 6.5 mm          | 8.7 mm  | 12.0 mm  | 331.0 mm |     |
|         | 7.0 mm          | 9.3 mm  | 13.5 mm  | 332.0 mm |     |
|         | 7.5 mm          | 10.0 mm | 13.5 mm  | 342.5 mm |     |
|         | 8.0 mm          | 10.7 mm | 16.0 mm  | 345.5 mm |     |
|         | 8.5 mm          | 11.3 mm | 16.0 mm  | 376.0 mm |     |
|         | 9.0 mm          | 12.0 mm | 16.0 mm  | 375.0 mm |     |
|         | 9.5 mm          | 12.7 mm | 18.0 mm  | 375.5 mm |     |
| 10.0 mm | 13.3 mm         | 18.0 mm | 378.0 mm |          |     |



**PAEDIATRIC ARMoured  
TRACHEAL TUBE, STRAIGHT,  
MADE OF SILICONE,  
WITHOUT CUFF**  
nasal/oral  
sizes: I.D. 2.0–6.5 mm

fixed connector, cupped atraumatic  
tip, radiopaque, graduated

- latex-free
- sterile

### PAEDIATRIC ARMoured TRACHEAL TUBE

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.   | LENGTH*  | QTY |
|--------|-----------------|--------|----------|-----|
| 105400 | 2.0 mm          | 3.3 mm | 223.0 mm | 2   |
|        | 2.5 mm          | 3.9 mm | 224.0 mm |     |
|        | 3.0 mm          | 4.3 mm | 225.5 mm |     |
|        | 3.5 mm          | 5.0 mm | 226.0 mm |     |
|        | 4.0 mm          | 5.7 mm | 246.0 mm |     |
|        | 4.5 mm          | 6.3 mm | 250.5 mm |     |
|        | 5.0 mm          | 6.7 mm | 263.0 mm |     |
|        | 5.5 mm          | 7.3 mm | 294.5 mm |     |
|        | 6.0 mm          | 8.0 mm | 326.0 mm |     |
|        | 6.5 mm          | 8.7 mm | 327.5 mm |     |

\* Length without connector

**Please note:** SILKOSPRAY is not suitable for devices made of silicone.  
If you require a lubricant, we recommend commercially available gels.

### ARMoured TRACHEAL TUBE, STRAIGHT, MADE OF SILKOLATEX, WITH SILICONE CUFF

nasal/oral

sizes: I.D. 3.0–10.5 mm

O.D. in CH 18.0–42.0

mounted connector, radiopaque,  
graduated, with cuff, pilot balloon  
and universal injection port

- sterile

**Tip:** Ref. 103040 without connector can be used  
for fiberoptic intubation along with a Mainz  
universal adaptor (see page 48)

**Attention:**

These products contain natural rubber latex which  
may cause allergic reactions.



### ARMoured TRACHEAL TUBE

**RÜSCH**

| REF.   | ORDER SIZE | I.D.    | O.D.    | LENGTH*  | QTY |
|--------|------------|---------|---------|----------|-----|
| 103040 | Ch. 18     | 3.0 mm  | 6.0 mm  | 229.0 mm | 5   |
|        | Ch. 20     | 3.5 mm  | 6.7 mm  | 232.5 mm |     |
|        | Ch. 22     | 4.0 mm  | 7.3 mm  | 232.5 mm |     |
|        | Ch. 24     | 5.0 mm  | 8.0 mm  | 253.0 mm |     |
|        | Ch. 26     | 5.5 mm  | 8.7 mm  | 253.0 mm |     |
|        | Ch. 28     | 6.0 mm  | 9.3 mm  | 313.5 mm |     |
|        | Ch. 30     | 6.5 mm  | 10.0 mm | 313.5 mm |     |
|        | Ch. 32     | 7.0 mm  | 10.7 mm | 314.0 mm |     |
|        | Ch. 34     | 8.0 mm  | 11.3 mm | 354.0 mm |     |
|        | Ch. 36     | 8.5 mm  | 12.0 mm | 354.5 mm |     |
|        | Ch. 38     | 9.0 mm  | 12.7 mm | 354.5 mm |     |
|        | Ch. 40     | 9.5 mm  | 13.3 mm | 355.0 mm |     |
|        | Ch. 42     | 10.5 mm | 14.0 mm | 355.0 mm |     |

## BRONCHOFLEX-SET

### ARMoured TRACHEAL/ INSERTION TUBE MADE OF PVC FOR USE IN BRONCHOSCOPY

(with separate oxygen tubing)

sizes: I.D. 7.5 mm and 8.5 mm

mounted connector, radiopaque,  
approx. 33 cm long, double lumen,  
graduated, siliconised, with or  
without cuff

- latex-free
- sterile

**Attention:**


We recommend spraying the inside of the  
bronchial tube and the outside of the suction  
catheter with SILKOSPRAY (Ref. 556000,  
see p. 49) prior to use.



### BRONCHOFLEX-SET

**RÜSCH**

| REF.                                                                                                                                                                                                                                                     | ORDER SIZE/I.D. | TOTAL O.D. | LENGTH* | QTY |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|------------|---------|-----|
| 104100                                                                                                                                                                                                                                                   | 7.5 mm          | 13.8 mm    | 319 mm  | 1   |
| <b>SET COMPONENTS</b> <ul style="list-style-type: none"> <li>• armoured tracheal/insertion tube made of PVC</li> <li>• oxygen tubing made of PVC, approx. 110 cm long</li> <li>• plastic bite block with fixation flange</li> <li>• neck band</li> </ul> |                 |            |         |     |



### BRONCHOFLEX-SET

**RÜSCH**

| REF.                                                                                                                                                                                                                                                                            | ORDER SIZE/I.D. | TOTAL O.D. | LENGTH* | QTY |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|------------|---------|-----|
| 104103                                                                                                                                                                                                                                                                          | 8.5 mm          | 11.7 mm    | 359 mm  | 1   |
| <b>SET COMPONENTS</b> <ul style="list-style-type: none"> <li>• armoured tracheal/insertion tube made of PVC with high-volume low-pressure cuff</li> <li>• oxygen tubing made of PVC, approx. 110 cm long</li> <li>• valve for Luer-lock syringes</li> <li>• neckband</li> </ul> |                 |            |         |     |

\* Length without connector



## LASER TUBE



**PREFORMED TRACHEAL TUBE  
WITH DEFINED LASER-RESISTANCE  
MADE OF WHITE SOFT RUBBER  
FOR LASER INTERVENTION IN  
THE LARYNGEAL SPACE**

oral/nasal (Magill), 1 eye (Murphy)  
sizes: I.D. 4.0–8.0 mm

fixed white connector, two-way stop-cocks with Luer and Luer-lock adapter and locking mechanism, cupped atraumatic tip, two cuffs inside one another, two pilot balloons, non-kinking shaft, Laserguard film 17 cm long consisting of Merocel® foam and corrugated silver foil

- sterile

**Attention:**

*These products contain natural rubber latex which may cause allergic reactions.*

### LASER TUBE

**RÜSCH**

| REF.    | ORDER SIZE/I.D. | O.D.    | CUFF Ø OUTSIDE | CUFF Ø INSIDE | LENGTH* | QTY |
|---------|-----------------|---------|----------------|---------------|---------|-----|
| 102004* | 4.0 mm          | 8.0 mm  | 20 mm          | 20 mm         | 400 mm  | 2   |
|         | 5.0 mm          | 9.0 mm  | 20 mm          | 20 mm         | 400 mm  |     |
|         | 6.0 mm          | 10.8 mm | 24 mm          | 24 mm         | 400 mm  |     |
|         | 7.0 mm          | 12.3 mm | 28 mm          | 28 mm         | 400 mm  |     |
|         | 8.0 mm          | 13.4 mm | 28 mm          | 28 mm         | 400 mm  |     |

\* Length without connector



## CLEAR MICROLARYNGEAL TUBE



**MADE OF PVC WITH HIGH-VOLUME LOW-PRESSURE CUFF**

oral/nasal

sizes: I.D. 4.0–6.0 mm

thinner tube body designed for interventions in the laryngeal space

semi-seated connector, valve for Luer and Luer-lock syringes, continuous X-ray marker, blue pilot balloon, graduated

- latex-free
- sterile



**MICROLARYNGEAL TUBE**

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.   | CUFF Ø | LENGTH* | QTY |
|--------|-----------------|--------|--------|---------|-----|
| 112460 | 4.0 mm          | 6.0 mm | 31 mm  | 360 mm  | 10  |
|        | 5.0 mm          | 7.3 mm | 31 mm  | 360 mm  |     |
|        | 6.0 mm          | 8.7 mm | 31 mm  | 360 mm  |     |

\* Length without connector

## OXFORD TUBE



**CLEAR ANATOMICALLY SHAPED NON-KINKING TUBE MADE OF PVC WITH LOW-PRESSURE CUFF**

sizes: I.D. 7.0–10.0 mm

semi-seated connector, valve for Luer and Luer-lock syringes, atraumatic tip, continuous X-ray marker, blue pilot balloon, graduated

- latex-free
- sterile



**OXFORD TUBE**

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | O.D.    | CUFF Ø | LENGTH*  | QTY |
|--------|-----------------|---------|--------|----------|-----|
| 112880 | 7.0 mm          | 10.0 mm | 26 mm  | 229.5 mm | 10  |
|        | 8.0 mm          | 11.3 mm | 28 mm  | 237.0 mm |     |
|        | 9.0 mm          | 12.7 mm | 29 mm  | 243.5 mm |     |
|        | 10.0 mm         | 14.0 mm | 29 mm  | 250.5 mm |     |

## SOFT RUBBER TRACHEAL TUBE



### MADE OF SOFT RED RUBBER, WITH CUFF

Magill  
sizes: I.D. 2.5–11.0 mm

semi-coated with SILKOLATEX®,  
cuff and pilot balloon made of  
SILKOLATEX®, universal injection  
port, graduated

- suitable for MRT
- without metal

#### Attention:


*These products contain natural rubber latex  
which may cause allergic reactions*

### SOFT RUBBER TRACHEAL TUBE

**RÜSCH**

| REF.    | ORDER SIZE/I.D. | O.D.    | LENGTH* | QTY |
|---------|-----------------|---------|---------|-----|
| 102000  | 2.5 mm          | 4.0 mm  | 140 mm  | 10  |
|         | 3.0 mm          | 4.7 mm  | 160 mm  |     |
|         | 3.5 mm          | 5.3 mm  | 180 mm  |     |
|         | 4.0 mm          | 6.0 mm  | 200 mm  |     |
|         | 4.5 mm          | 6.7 mm  | 220 mm  |     |
|         | 5.0 mm          | 7.3 mm  | 240 mm  |     |
|         | 5.5 mm          | 8.0 mm  | 270 mm  |     |
|         | 6.0 mm          | 8.7 mm  | 280 mm  |     |
|         | 6.5 mm          | 9.3 mm  | 290 mm  |     |
|         | 7.0 mm          | 10.0 mm | 300 mm  |     |
|         | 7.5 mm          | 10.7 mm | 310 mm  |     |
|         | 8.0 mm          | 11.3 mm | 320 mm  |     |
|         | 8.5 mm          | 12.0 mm | 330 mm  |     |
|         | 9.0 mm          | 12.7 mm | 340 mm  |     |
|         | 9.5 mm          | 13.3 mm | 350 mm  |     |
|         | 10.0 mm         | 14.0 mm | 360 mm  |     |
| 10.5 mm | 14.7 mm         | 360 mm  |         |     |
| 11.0 mm | 15.3 mm         | 360 mm  |         |     |

\* Length without connector



### STANDARD CONNECTOR

**RÜSCH**

| REF.   | ORDER SIZE/I.D. | DESCRIPTION                                          | QTY |
|--------|-----------------|------------------------------------------------------|-----|
| 501003 | 1.5 mm–11.0 mm  | made of plastic, latex-free, 15 mm connection (O.D.) | 10  |

# AID INSTILLATION CATHETER

**FOR DRUG INSTILLATION DURING CPR AND SURFACTANT APPLICATION IN NEONATOLOGY**

with sealable Luer-lock adapter, multiple perforations over a length of 1 cm, opening 1 cm from distal end of the catheter, with intubation stylet, graduated

- latex-free
- sterile

# EDGAR TUBE

**CLEAR TRACHEAL TUBE MADE OF PVC FOR ENDO-BRONCHIAL DRUG & GAS APPLICATION DURING RESUSCITATION**  
nasal/oral (Magill)

semi-seated connector, additional instillation channel with sealable Luer-lock adapter, continuous X-ray marker, graduated

**Ref. 111480**

sizes: I.D. 6.5–10.0 mm with low-pressure cuff, valve for Luer and Luer-lock syringes and black position indicator for correct tube placement


**Ref. 111380**

sizes: I.D. 2.5–6.0 mm without cuff, with black tip

- latex-free
- sterile



| AID ADAPTER SET             |           |                                | RÜSCH |
|-----------------------------|-----------|--------------------------------|-------|
| REF.                        | SIZE/O.D. | DESCRIPTION                    | QTY   |
| 111100                      | Ch. 4     | catheter with angled connector | 1     |
| SET COMPONENTS              |           |                                |       |
| • AID instillation catheter |           |                                |       |
| • angled plastic connector  |           |                                |       |



| EDGAR TUBE |                 |         |        |         | RÜSCH |
|------------|-----------------|---------|--------|---------|-------|
| REF.       | ORDER SIZE/I.D. | O.D.    | CUFF Ø | LENGTH* | QTY   |
| 111480     | 6.5 mm          | 9.3 mm  | 24 mm  | 290 mm  | 10    |
|            | 7.0 mm          | 10.0 mm | 26 mm  | 300 mm  |       |
|            | 7.5 mm          | 10.7 mm | 26 mm  | 310 mm  |       |
|            | 8.0 mm          | 11.3 mm | 28 mm  | 320 mm  |       |
|            | 8.5 mm          | 12.0 mm | 28 mm  | 330 mm  |       |
|            | 9.0 mm          | 12.7 mm | 29 mm  | 340 mm  |       |
|            | 9.5 mm          | 13.3 mm | 29 mm  | 350 mm  |       |
|            | 10.0 mm         | 14.0 mm | 29 mm  | 360 mm  |       |
| 111380     | 2.5 mm          | 4.0 mm  | –      | 140 mm  | 10    |
|            | 3.0 mm          | 4.7 mm  | –      | 160 mm  |       |
|            | 3.5 mm          | 5.3 mm  | –      | 180 mm  |       |
|            | 4.0 mm          | 6.0 mm  | –      | 200 mm  |       |
|            | 4.5 mm          | 6.7 mm  | –      | 220 mm  |       |
|            | 5.0 mm          | 7.3 mm  | –      | 240 mm  |       |
|            | 5.5 mm          | 8.0 mm  | –      | 270 mm  |       |
|            | 6.0 mm          | 8.7 mm  | –      | 280 mm  |       |

\* Length without connector

# SHERIDAN ENDOTRACHEAL TUBES

Supplied under the SHERIDAN brand name, our endotracheal tubes offer all the features needed by today's health-care professionals, including a high level of tactile feed-back from the pilot balloon, thermosensitive PVC,

which allows the tubes to soften at body temperature, as well as a smooth tube tip that is moulded and gently cupped to minimise trauma during intubation.

## SHERIDAN ISIS HVT

**ENDOTRACHEAL TUBE WITH INTEGRATED PORT FOR OPTIONAL SUBGLOTTIC SECRETION SUCTIONING WITH HIGH VOLUME TAPERED CUFF**  
nasal/oral, 1 eye (Murphy)

The ideal tube to provide the best of care to all your patients, both over the short- & long-term.

The unique convertible nature of the Teleflex ISIS HVT frees clinicians from the burden of choosing which tube is best for the patient at the time of intubation. ISIS features an integrated suction port and separate suction line allowing for subglottic secretion removal whenever necessary. Both connection ports can be sealed upon disconnection, reducing the risk of cross-contamination when not in use. When needed, the suction tube attaches to the ISIS HVT via a secure locking connection.

- sterile



### SHERIDAN ISIS HVT

**SHERIDAN**

| REF.    | ORDER SIZE/I.D. | QTY |
|---------|-----------------|-----|
| 5-13012 | 6.0 mm          | 10  |
| 5-13013 | 6.5 mm          |     |
| 5-13014 | 7.0 mm          |     |
| 5-13015 | 7.5 mm          |     |
| 5-13016 | 8.0 mm          |     |
| 5-13017 | 8.5 mm          |     |
| 5-13018 | 9.0 mm          |     |



### SUCTION LINE

| REF.    | QTY |
|---------|-----|
| 5-23000 | 5   |

# SHERIDAN/HVT® TRACHEAL TUBE



SHERIDAN/HVT TRACHEAL TUBE **SHERIDAN**

| REF.    | ORDER SIZE/I.D. | QTY |
|---------|-----------------|-----|
| 5-10310 | 5.0 mm          | 10  |
| 5-10311 | 5.5 mm          |     |
| 5-10312 | 6.0 mm          |     |
| 5-10313 | 6.5 mm          |     |
| 5-10314 | 7.0 mm          |     |
| 5-10315 | 7.5 mm          |     |
| 5-10316 | 8.0 mm          |     |
| 5-10317 | 8.5 mm          |     |
| 5-10318 | 9.0 mm          |     |
| 5-10320 | 10.0 mm         |     |

**HIGH-VOLUME TAPERED, LOW-PRESSURE CUFF**  
nasal/oral, 1 eye (Murphy)

The high-volume, large diameter cuff provides a positive tracheal wallseal, while the low-pressure cuff minimises capillary restriction. The tapered cuff reduces the risk of micro-aspiration, increases ease of use & minimises the risk of herniation.

- sterile

# SHERIDAN/CF® TRACHEAL TUBE



SHERIDAN/CF TRACHEAL TUBE **SHERIDAN**

| REF.           | ORDER SIZE/I.D. | REF.          | ORDER SIZE/I.D. | QTY |
|----------------|-----------------|---------------|-----------------|-----|
| CF Magill Type |                 | CF Murphy Eye |                 | 10  |
| 5-10206        | 3.0 mm          | 5-10106       | 3.0 mm          |     |
|                | -               | 5-10107       | 3.5 mm          |     |
| 5-10208        | 4.0 mm          | 5-10108       | 4.0 mm          |     |
|                | -               | 5-10109       | 4.5 mm          |     |
| 5-10210        | 5.0 mm          | 5-10110       | 5.0 mm          |     |
|                | -               | 5-10111       | 5.5 mm          |     |
| 5-10212        | 6.0 mm          | 5-10112       | 6.0 mm          |     |
| 5-10213        | 6.5 mm          | 5-10113       | 6.5 mm          |     |
| 5-10214        | 7.0 mm          | 5-10114       | 7.0 mm          |     |
| 5-10215        | 7.5 mm          | 5-10115       | 7.5 mm          |     |
| 5-10216        | 8.0 mm          | 5-10116       | 8.0 mm          |     |
| 5-10217        | 8.5 mm          | 5-10117       | 8.5 mm          |     |
| 5-10218        | 9.0 mm          | 5-10118       | 9.0 mm          |     |
| 5-10219        | 9.5 mm          | 5-10119       | 9.5 mm          |     |
| 5-10220        | 10.0 mm         | 5-10120       | 10.0 mm         |     |

**CLOSE FITTING, LOW-PRESSURE CUFF**  
nasal/oral, 1 eye (Murphy or Magill)

By providing ease of intubation and increasing visualisation, the opaque, close fitting, low- pressure cuff is ideal for short-term surgical intubations.

- sterile

## SHERIDAN UNCUFFED™ TRACHEAL TUBE



### SHERIDAN UNCUFFED TRACHEAL TUBE

**SHERIDAN**

| REF.    | ORDER SIZE/I.D. | QTY |
|---------|-----------------|-----|
| 5-10404 | 2.0 mm          | 10  |
| 5-10405 | 2.5 mm          |     |
| 5-10406 | 3.0 mm          |     |
| 5-10407 | 3.5 mm          |     |
| 5-10408 | 4.0 mm          |     |
| 5-10409 | 4.5 mm          |     |
| 5-10410 | 5.0 mm          |     |
| 5-10411 | 5.5 mm          |     |
| 5-10412 | 6.0 mm          |     |
| 5-10413 | 6.5 mm          |     |
| 5-10414 | 7.0 mm          |     |

### MADE OF PVC, WITHOUT CUFF

nasal/oral, 1 eye (Murphy)  
cut on (\*) for oral use

Special depth lines at the distal tip indicate the depth of intubation below the chords. The smooth, moulded, gently cupped tip minimises trauma during intubation, while the tube softens at body temperature to conform to the airway anatomy.

- sterile

## SHERIDAN PED-SOFT™ UNCUFFED TRACHEAL TUBE



### SHERIDAN PED-SOFT UNCUFFED TRACHEAL TUBE

**SHERIDAN**

| REF.    | ORDER SIZE/I.D. | DISTAL TIP TO END OF MARK | QTY |
|---------|-----------------|---------------------------|-----|
| 5-30404 | 2.0 mm          | 20 mm                     | 10  |
| 5-30405 | 2.5 mm          | 22 mm                     |     |
| 5-30406 | 3.0 mm          | 24 mm                     |     |
| 5-30407 | 3.5 mm          | 26 mm                     |     |
| 5-30408 | 4.0 mm          | 28 mm                     |     |
| 5-30409 | 4.5 mm          | 30 mm                     |     |
| 5-30410 | 5.0 mm          | 32 mm                     |     |
| 5-30411 | 5.5 mm          | 34 mm                     |     |
| 5-30412 | 6.0 mm          | 36 mm                     |     |
| 5-30413 | 6.5 mm          | 39 mm                     |     |
| 5-30414 | 7.0 mm          | 41 mm                     |     |

### MADE OF PVC, WITHOUT CUFF

nasal/oral, 1 eye (Murphy)

Made of soft PVC, the tube conforms to a child's anatomy. Furthermore, a well defined mark at the distal tip serves as reference point during intubation. Additional depth marks assist in placement during nasal intubation, while the smooth moulded, gently cupped tip minimises trauma during intubation.

- sterile

# SHERIDAN PREFORMED™ TRACHEAL TUBE

## HIGH-VOLUME TAPERED, LOW-PRESSURE CUFF

nasal/oral, 1 eye (Murphy) with cuffed version  
 nasal/oral, 2 eyes (Murphy) with uncuffed version

Sheridan Preformed Tracheal Tubes for oral and maxillofacial surgery allow the anaesthesia circuit connection to be positioned out of the surgical field.

Features include: bold marks at the center of the bend with distance to distal tip indicated, high-volume, large diameter cuff for positive tracheal wall seal, low-pressure cuff to minimise capillary restrictions.

Uncuffed tracheal tubes have two Murphy eyes for enhanced patient safety.

- sterile



### SHERIDAN PREFORMED TRACHEAL TUBE, ORAL



| REF.          | ORDER SIZE/I.D. | REF.        | ORDER SIZE/I.D. | QTY |
|---------------|-----------------|-------------|-----------------|-----|
| Uncuffed oral |                 | Cuffed oral |                 | 10  |
| 5-22006       | 3.0 mm          | -           | -               |     |
| 5-22007       | 3.5 mm          | -           | -               |     |
| 5-22008       | 4.0 mm          | 5-22208     | 4.0 mm          |     |
| 5-22009       | 4.5 mm          | 5-22209     | 4.5 mm          |     |
| 5-22010       | 5.0 mm          | 5-22210     | 5.0 mm          |     |
| 5-22011       | 5.5 mm          | 5-22211     | 5.5 mm          |     |
| 5-22012       | 6.0 mm          | 5-22212     | 6.0 mm          |     |
| 5-22013       | 6.5 mm          | 5-22213     | 6.5 mm          |     |
| 5-22014       | 7.0 mm          | 5-22214     | 7.0 mm          |     |
|               |                 | 5-22215     | 7.5 mm          |     |
|               |                 | 5-22216     | 8.0 mm          |     |
|               |                 | 5-22217     | 8.5 mm          |     |
|               |                 | 5-22218     | 9.0 mm          |     |



### SHERIDAN PREFORMED TRACHEAL TUBE, NASAL



| REF.           | ORDER SIZE/I.D. | REF.         | ORDER SIZE/I.D. | QTY |
|----------------|-----------------|--------------|-----------------|-----|
| Uncuffed nasal |                 | Cuffed nasal |                 | 10  |
| 5-22106        | 3.0 mm          | -            | -               |     |
| 5-22107        | 3.5 mm          | -            | -               |     |
| 5-22108        | 4.0 mm          | -            | -               |     |
| 5-22109        | 4.5 mm          | -            | -               |     |
| 5-22110        | 5.0 mm          | -            | -               |     |
| 5-22111        | 5.5 mm          | -            | -               |     |
| 5-22112        | 6.0 mm          | 5-22312      | 6.0 mm          |     |
| 5-22113        | 6.5 mm          | 5-22313      | 6.5 mm          |     |
| 5-22114        | 7.0 mm          | 5-22314      | 7.0 mm          |     |
|                |                 | 5-22315      | 7.5 mm          |     |
|                |                 | 5-22316      | 8.0 mm          |     |



# SPIRAL-FLEX REINFORCED TRACHEAL TUBE

## HIGH-VOLUME TAPERED, LOW-PRESSURE CUFF

nasal/oral, with cuffed version

nasal/oral, 1 eye (Murphy) with uncuffed version



reduced risk of kinking due to stainless steel, spiral wound reinforcing wire within the tube wall; 15 mm connector is permanently bonded within the tube, blue cuff

- sterile

### CUFFED ORAL

cuffed, oral pre-cut length with reinforcing sleeve to minimise the risk of patient biting and occluding the tube; the oral length tube has a reinserted SHER-I-SLIP stylet.


| SPIRAL-FLEX REINFORCED TRACHEAL TUBE |                 |                   |                 |                     |                 | SHERIDAN |  |
|--------------------------------------|-----------------|-------------------|-----------------|---------------------|-----------------|----------|--|
| REF.                                 | ORDER SIZE/I.D. | REF.              | ORDER SIZE/I.D. | REF.                | ORDER SIZE/I.D. | QTY      |  |
| Cuffed oral                          |                 | Cuffed oral/nasal |                 | Uncuffed oral/nasal |                 | 5        |  |
|                                      | -               |                   | -               | 5-12906             | 3.0 mm          |          |  |
|                                      | -               |                   | -               | 5-12907             | 3.5 mm          |          |  |
|                                      | -               |                   | -               | 5-12908             | 4.0 mm          |          |  |
|                                      | -               |                   | -               | 5-12909             | 4.5 mm          |          |  |
| 5-12510                              | 5.0 mm          | 5-12610           | 5.0 mm          | 5-12910             | 5.0 mm          |          |  |
| 5-12511                              | 5.5 mm          | 5-12611           | 5.5 mm          |                     | -               |          |  |
| 5-12512                              | 6.0 mm          | 5-12612           | 6.0 mm          |                     | -               |          |  |
| 5-12513                              | 6.5 mm          | 5-12613           | 6.5 mm          |                     | -               |          |  |
| 5-12514                              | 7.0 mm          | 5-12614           | 7.0 mm          |                     | -               |          |  |
| 5-12515                              | 7.5 mm          | 5-12615           | 7.5 mm          |                     | -               |          |  |
| 5-12516                              | 8.0 mm          | 5-12616           | 8.0 mm          |                     | -               |          |  |
| 5-12517                              | 8.5 mm          | 5-12617           | 8.5 mm          |                     | -               |          |  |
| 5-12518                              | 9.0 mm          | 5-12618           | 9.0 mm          |                     | -               |          |  |

## LITA CUFFED TRACHEAL TUBE

**HIGH-VOLUME TAPERED, LOW-PRESSURE CUFF FOR LARYNGO-TRACHEAL INSTILLATION OF TOPICAL ANAESTHETIC**  
nasal/oral, 1 eye (Murphy)

allows for spray application of topical anaesthetic via specially designed lumen; patient tolerance of endotracheal tube is increased after administration of topical anaesthesia; black reference bar indicates the position of the uppermost opening

- sterile
- individually packed



### LITA CUFFED TRACHEAL TUBE *SHERIDAN*

| REF.    | ORDER SIZE/I.D. | QTY |
|---------|-----------------|-----|
| 5-20512 | 6.0 mm          | 10  |
| 5-20513 | 6.5 mm          |     |
| 5-20514 | 7.0 mm          |     |
| 5-20515 | 7.5 mm          |     |
| 5-20516 | 8.0 mm          |     |
| 5-20517 | 8.5 mm          |     |

## SHERIDAN LTS™ TRACHEAL TUBE FOR MICROLARYNGEAL TRACHEAL SURGERY

nasal/oral, 1 eye (Murphy)

small O.D. tube allows greater access to the surgical field, yet provides sufficient ventilation of the patient; same tube length and cuff diameter as an 8.0mm cuffed tracheal tube; high-volume, yellow cuff is more visible and provides a low-pressure tracheal wall seal

- sterile
- individually packed



### SHERIDAN LTS *SHERIDAN*

| REF.    | ORDER SIZE/I.D. | SIZE/O.D. | QTY |
|---------|-----------------|-----------|-----|
| 5-11108 | 4.0 mm          | 5.8 mm    | 10  |
| 5-11110 | 5.0 mm          | 7.1 mm    |     |
| 5-11112 | 6.0 mm          | 8.5 mm    |     |

# LARYNGEAL MASKS

The true benefits of laryngeal masks are firstly that they lead to less pain and coughing than endotracheal tubes, and secondly, that they can be inserted almost blindly. As such laryngeal masks are used for airway management in both anaesthesia and emergency medicine. Thanks to their soft, thin-walled cuff, they

efficiently help to prevent potential traumas caused by excessive pressure. To top it all off, the Teleflex Sure Seal laryngeal mask features an integrated unique cuff pressure device providing at-a-glance feedback on cuff pressure and alerting you instantly to changes before they affect patient safety.

## SURE SEAL & CUFF PILOT™

SURE SEAL combines all the features you would expect to find in a high-quality laryngeal mask. The mask itself is made of 100 % silicone and is finely shaped to perfectly contour the oropharyngeal space without trauma. The tube is available in both silicone and PVC according to preference.

A reinforced laryngeal mask complements the range. Our Sure Seal product line is further distinguished by a unique cuff pressure monitoring device that enables you to improve patient care and safety: the Cuff Pilot™ provides constant, at-a-glance, visual information on the pressure in the mask.



*Packaging colour coding for the disposable Sure Seal product*

| SIZE   | 1 | 1.5    | 2      | 2.5        | 3     | 4      | 5   | 6    |
|--------|---|--------|--------|------------|-------|--------|-----|------|
| COLOUR |   | Yellow | Purple | Light Blue | Green | Yellow | Red | Pink |

The world's first integrated indicator that constantly provides at-a-glance feedback on the pressure inside the cuff, alerting you instantly to changes before they affect patient safety.




**SINGLE USE LARYNGEAL MASK WITH CUFF PILOT™**

This laryngeal mask offers the advantages of a reusable mask (softness, flexibility), while avoiding the infection risks inherent in resterilization. Furthermore, it provides you with the comfort & safety of the Cuff Pilot.

- latex-free
- sterile


**TRACHJELL  
SINGLE DOSE WATER-SOLUBLE LUBRICANT**

To ensure successful intubation, a laryngeal mask needs to be lubricated. Teleflex now provides a water-soluble lubricant, sterile, single-dose, quick & easy to use.




**TRACHJELL**

| REF.   | DESCRIPTION                                                                                                           | QTY |
|--------|-----------------------------------------------------------------------------------------------------------------------|-----|
| 556100 | <ul style="list-style-type: none"> <li>• water soluble lubricating gel</li> <li>• 8.5 g</li> <li>• sterile</li> </ul> | 25  |



**SURE SEAL WITH CUFF PILOT™, SILICONE CUFF & PVC TUBE**

| REF.   | SIZE 1  | SIZE 1.5 | SIZE 2  | SIZE 2.5 | SIZE 3  | SIZE 4  | SIZE 5  | SIZE 6  |  |
|--------|---------|----------|---------|----------|---------|---------|---------|---------|--|
| 105200 | -000010 | -000015  | -000020 | -000025  | -000030 | -000040 | -000050 | -000060 |  |
| QTY    | 10      |          |         |          |         |         |         | 5       |  |



**SURE SEAL WITH CUFF PILOT™, 100% SILICONE**

| REF.   | SIZE 1  | SIZE 1.5 | SIZE 2  | SIZE 2.5 | SIZE 3  | SIZE 4  | SIZE 5  | SIZE 6  |  |
|--------|---------|----------|---------|----------|---------|---------|---------|---------|--|
| 105210 | -000010 | -000015  | -000020 | -000025  | -000030 | -000040 | -000050 | -000060 |  |
| QTY    | 10      |          |         |          |         |         |         | 5       |  |




**SURE SEAL WITH CUFF PILOT™, 100% SILICONE REINFORCED**

| REF.   | SIZE 1  | SIZE 1.5 | SIZE 2  | SIZE 2.5 | SIZE 3  | SIZE 4  | SIZE 5  | SIZE 6  |  |
|--------|---------|----------|---------|----------|---------|---------|---------|---------|--|
| 105220 | -000010 | -000015  | -000020 | -000025  | -000030 | -000040 | -000050 | -000060 |  |
| QTY    | 10      |          |         |          |         |         |         | 5       |  |

**SINGLE USE LARYNGEAL MASK WITH STANDARD PILOT BALLOON**


This mask offers you the quality and features you would expect from Teleflex's new range of laryngeal masks. Excellent materials and an outstanding design enable you to concentrate on your patients' breathing rather than on the device.

- latex-free
- sterile




**SURE SEAL, SILICONE CUFF & PVC TUBE**

| REF.   | SIZE 1  | SIZE 1.5 | SIZE 2  | SIZE 2.5 | SIZE 3  | SIZE 4  | SIZE 5  | SIZE 6  |  |
|--------|---------|----------|---------|----------|---------|---------|---------|---------|--|
| 105300 | -000010 | -000015  | -000020 | -000025  | -000030 | -000040 | -000050 | -000060 |  |
| QTY    | 10      |          |         |          |         |         |         | 5       |  |



**SURE SEAL, 100% SILICONE**

| REF.   | SIZE 1  | SIZE 1.5 | SIZE 2  | SIZE 2.5 | SIZE 3  | SIZE 4  | SIZE 5  | SIZE 6  |  |
|--------|---------|----------|---------|----------|---------|---------|---------|---------|--|
| 105310 | -000010 | -000015  | -000020 | -000025  | -000030 | -000040 | -000050 | -000060 |  |
| QTY    | 10      |          |         |          |         |         |         | 5       |  |



**SURE SEAL, 100% SILICONE REINFORCED**

| REF.   | SIZE 1  | SIZE 1.5 | SIZE 2  | SIZE 2.5 | SIZE 3  | SIZE 4  | SIZE 5  | SIZE 6  |  |
|--------|---------|----------|---------|----------|---------|---------|---------|---------|--|
| 105320 | -000010 | -000015  | -000020 | -000025  | -000030 | -000040 | -000050 | -000060 |  |
| QTY    | 10      |          |         |          |         |         |         | 5       |  |

**LARYNGEAL MASKS**

| SIZE | PATIENT WEIGHT | CUFF MAX VOL. |
|------|----------------|---------------|
| 1.0  | < 5 kg         | 5 ml          |
| 1.5  | 5–10 kg        | 7 ml          |
| 2.0  | 10–20 kg       | 10 ml         |
| 2.5  | 20–30 kg       | 15 ml         |
| 3.0  | 30–50 kg       | 20 ml         |
| 4.0  | 50–70 kg       | 35 ml         |
| 5.0  | 70–100 kg      | 48 ml         |
| 6.0  | > 100 kg       | 50 ml         |

### REUSABLE LARYNGEAL MASK WITH STANDARD PILOT BALLOON


Our reusable standard laryngeal mask incorporates everything you need for easy and safe insertion, as well as a reliable seal throughout the procedure.

- latex-free



#### SURE SEAL, REUSABLE, 100% SILICONE

| REF.   | SIZE 1  | SIZE 1.5 | SIZE 2  | SIZE 2.5 | SIZE 3  | SIZE 4  | SIZE 5  | SIZE 6  |
|--------|---------|----------|---------|----------|---------|---------|---------|---------|
| 105000 | -000010 | -000015  | -000020 | -000025  | -000030 | -000040 | -000050 | -000060 |
| QTY    | 2       |          |         |          |         |         |         |         |



#### SURE SEAL, REUSABLE, 100% SILICONE REINFORCED

| REF.   | SIZE 1  | SIZE 1.5 | SIZE 2  | SIZE 2.5 | SIZE 3  | SIZE 4  | SIZE 5  | SIZE 6  |
|--------|---------|----------|---------|----------|---------|---------|---------|---------|
| 105010 | -000010 | -000015  | -000020 | -000025  | -000030 | -000040 | -000050 | -000060 |
| QTY    | 2       |          |         |          |         |         |         |         |

## CRYSTAL AIRWAY MASK®

### SINGLE USE LARYNGEAL MASK

Our standard full PVC laryngeal mask

rotationally-moulded cuff with a strong yet thin soft reinforced posterior wall of the cuff, inflation channel integrated into the tube, colour-coded according to size, mask size & patient weight, cuff volume indicated on the tube, black line on the tube for better placement, 15 mm universal connector, 7 sizes

- latex-free
- sterile



#### CRYSTAL AIRWAY MASK

**RÜSCH**

| REF.   | SIZE | I.D.    | O.D.   | COLOUR CODE | QTY |
|--------|------|---------|--------|-------------|-----|
| 111000 | 1.0  | 6.2 mm  | 8.6 mm | clear       | 10  |
|        | 1.5  | 7.2 mm  | 10 mm  | blue        |     |
|        | 2.0  | 8.6 mm  | 12 mm  | green       |     |
|        | 2.5  | 10.4 mm | 14 mm  | purple      |     |
|        | 3.0  | 11.0 mm | 16 mm  | orange      |     |
|        | 4.0  | 12.4 mm | 18 mm  | pink        |     |
|        | 5.0  | 12.4 mm | 18 mm  | yellow      |     |

# RÜSCH BRONCHIAL TUBES AND BLOCKERS

All RÜSCH endobronchial tubes are characterised by high-quality workmanship. They are made of highly flexible and thermoplastic materials. Thanks to special features like the special shape of the bronchial cuff on right-sided models or the extremely soft Carina hook (Carlens and White type), RÜSCH endobronchial tubes are both safe and patient-friendly.

This is particularly true in case of our EZ-Blocker. Thanks to its y-shaped design, this innovative bronchial blocker anchors itself on the carina and will not dislocate following inflation of the isolated lung. After the procedure, when the EZ-Blocker has been removed, there is no need to reintubate the patient.

## SINGLE-LUMEN BRONCHIAL TUBES

### BRONCHIAL TUBE MADE OF PVC


semi-seated connector, valves for Luer and Luer-lock syringes, colourless tracheal low-pressure cuff with colourless pilot balloon, blue bronchial low-pressure cuff with blue pilot balloon, continuous X-ray marker and additional cuff markers, stylet, graduated

- latex-free
- sterile



#### ANGLED CONNECTOR


- made of transparent plastic
- double sealing cap made of silicone with small closable opening
- rotating connectors



### SINGLE-LUMEN BRONCHIAL TUBE SET

**RÜSCH**

| REF.                                            | ORDER SIZE/I.D.                     | MAX. O.D. | CUFF Ø BRONCH. | CUFF Ø TRACHEAL | LENGTH* | QTY |
|-------------------------------------------------|-------------------------------------|-----------|----------------|-----------------|---------|-----|
| 115900 (set)                                    | for left-sided bronchial intubation |           |                |                 |         | 1   |
|                                                 | 6.5 mm                              | 9.3 mm    | 16.5 mm        | 22±3.3 mm       | 450 mm  |     |
|                                                 | 8.0 mm                              | 11.3 mm   | 18 mm          | 28±4.2 mm       | 470 mm  |     |
| SET COMPONENTS                                  |                                     |           |                |                 |         |     |
| • single-lumen bronchial tube made of PVC, left |                                     |           |                |                 |         |     |
| • angled connector                              |                                     |           |                |                 |         |     |



### SINGLE-LUMEN BRONCHIAL TUBE SET

**RÜSCH**

| REF.                                             | ORDER SIZE/I.D.                      | MAX. O.D. | CUFF Ø BRONCH. | CUFF Ø TRACHEAL | LENGTH* | QTY |
|--------------------------------------------------|--------------------------------------|-----------|----------------|-----------------|---------|-----|
| 115901 (set)                                     | for right-sided bronchial intubation |           |                |                 |         | 1   |
|                                                  | 6.5 mm                               | 9.3 mm    | 16.5 mm        | 22±3.3 mm       | 440 mm  |     |
|                                                  | 8.0 mm                               | 11.3 mm   | 18 mm          | 28±4.2 mm       | 460 mm  |     |
| SET COMPONENTS                                   |                                      |           |                |                 |         |     |
| • single-lumen bronchial tube made of PVC, right |                                      |           |                |                 |         |     |
| • angled connector                               |                                      |           |                |                 |         |     |

\* Shaft length

# BRONCHOPART



## CLEAR DOUBLE-LUMEN BRONCHIAL TUBE MADE OF PVC



semi-seated connectors, valves for Luer and Luer-lock syringes, continuous X-ray marker and additional cuff markers, colourless tracheal low-pressure cuff with colourless pilot balloon, blue bronchial low-pressure cuff with blue pilot balloon, forming stylet with fixation, graduated

- latex-free
- sterile

### BRONCHOPART TUBE SET


**RÜSCH**

| REF.                         | ORDER SIZE                           | 2 X I.D. | MAX. O.D. | CUFF Ø BRONCH. | CUFF Ø TRACHEAL | LENGTH* | QTY |
|------------------------------|--------------------------------------|----------|-----------|----------------|-----------------|---------|-----|
| 116100 (set)                 | for left-sided bronchial intubation  |          |           |                |                 |         | 1   |
| 116162 (without accessories) | 26                                   | 3.40 mm  | 9.3 mm    | 15.0 mm        | 22±3.3 mm       | 280 mm  |     |
|                              | 28                                   | 3.71 mm  | 10.0 mm   | 15.0 mm        | 22±3.3 mm       | 280 mm  |     |
|                              | 35                                   | 4.26 mm  | 12.5 mm   | 16.5 mm        | 30±4.5 mm       | 300 mm  |     |
|                              | 37                                   | 4.52 mm  | 13.2 mm   | 16.5 mm        | 30±4.5 mm       | 310 mm  |     |
|                              | 39                                   | 4.76 mm  | 13.9 mm   | 18.0 mm        | 31±4.6 mm       | 330 mm  |     |
|                              | 41                                   | 5.04 mm  | 14.6 mm   | 18.0 mm        | 31±4.6 mm       | 340 mm  |     |
|                              | SET COMPONENTS                       |          |           |                |                 |         |     |
|                              | • double-lumen bronchial tube, left  |          |           |                |                 |         |     |
|                              | • 2 suction catheters                |          |           |                |                 |         |     |
|                              | • 2 angled connectors                |          |           |                |                 |         |     |
|                              | • Y connector                        |          |           |                |                 |         |     |
| 116200 (set)                 | for right-sided bronchial intubation |          |           |                |                 |         | 1   |
| 116262 (without accessories) | 26                                   | 3.40 mm  | 9.3 mm    | 15.0 mm        | 22±3.3 mm       | 280 mm  |     |
|                              | 28                                   | 3.71 mm  | 10.0 mm   | 15.0 mm        | 22±3.3 mm       | 280 mm  |     |
|                              | 35                                   | 4.26 mm  | 12.5 mm   | 16.5 mm        | 30±4.5 mm       | 300 mm  |     |
|                              | 37                                   | 4.52 mm  | 13.2 mm   | 16.5 mm        | 30±4.5 mm       | 310 mm  |     |
|                              | 39                                   | 4.76 mm  | 13.9 mm   | 18.0 mm        | 31±4.6 mm       | 340 mm  |     |
|                              | 41                                   | 5.04 mm  | 14.6 mm   | 18.0 mm        | 31±4.6 mm       | 350 mm  |     |
|                              | SET COMPONENTS                       |          |           |                |                 |         |     |
|                              | • double-lumen bronchial tube, right |          |           |                |                 |         |     |
|                              | • 2 suction catheters                |          |           |                |                 |         |     |
|                              | • 2 angled connectors                |          |           |                |                 |         |     |
|                              | • Y connector                        |          |           |                |                 |         |     |



\* Shaft length






**BRONCHOPART CARLENS TUBE SET**

**RÜSCH**


| REF.         | ORDER SIZE                                                                                                                                                                                               | 2 X I.D. | MAX. O.D. | CUFF Ø BRONCH. | CUFF Ø TRACHEAL | LENGTH* | QTY |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|-----------|----------------|-----------------|---------|-----|
| 116101 (set) | for left-sided bronchial intubation with Carina hook                                                                                                                                                     |          |           |                |                 |         | 1   |
|              | 35                                                                                                                                                                                                       | 4.26 mm  | 12.5 mm   | 16.5 mm        | 30±4.5 mm       | 300 mm  |     |
|              | 37                                                                                                                                                                                                       | 4.52 mm  | 13.2 mm   | 16.5 mm        | 30±4.5 mm       | 310 mm  |     |
|              | 39                                                                                                                                                                                                       | 4.76 mm  | 13.9 mm   | 18.0 mm        | 31±4.6 mm       | 330 mm  |     |
|              | 41                                                                                                                                                                                                       | 5.04 mm  | 14.6 mm   | 18.0 mm        | 31±4.6 mm       | 340 mm  |     |
|              | <b>SET COMPONENTS</b> <ul style="list-style-type: none"> <li>• Carlens double-lumen bronchial tube, left</li> <li>• 2 suction catheters</li> <li>• 2 angled connectors</li> <li>• Y connector</li> </ul> |          |           |                |                 |         |     |



**BRONCHOPART WHITE TUBE SET**

**RÜSCH**

| REF.         | ORDER SIZE                                                                                                                                                                                              | 2 X I.D. | MAX. O.D. | CUFF Ø BRONCH. | CUFF Ø TRACHEAL | LENGTH* | QTY |
|--------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|-----------|----------------|-----------------|---------|-----|
| 116201 (set) | for right-sided bronchial intubation with Carina hook                                                                                                                                                   |          |           |                |                 |         | 1   |
|              | 35                                                                                                                                                                                                      | 4.26 mm  | 12.5 mm   | 16.5 mm        | 30±4.5 mm       | 300 mm  |     |
|              | 37                                                                                                                                                                                                      | 4.52 mm  | 13.2 mm   | 16.5 mm        | 30±4.5 mm       | 310 mm  |     |
|              | 39                                                                                                                                                                                                      | 4.76 mm  | 13.9 mm   | 18.0 mm        | 31±4.6 mm       | 340 mm  |     |
|              | 41                                                                                                                                                                                                      | 5.04 mm  | 14.6 mm   | 18.0 mm        | 31±4.6 mm       | 350 mm  |     |
|              | <b>SET COMPONENTS</b> <ul style="list-style-type: none"> <li>• White double-lumen bronchial tube, right</li> <li>• 2 suction catheters</li> <li>• 2 angled connectors</li> <li>• Y connector</li> </ul> |          |           |                |                 |         |     |



**ACCESSORIES INCLUDED IN SET**



**SUCTION CATHETER**

- made of plastic
- transparent
- with suction control



**ANGLED CONNECTORS**

- with rotating connectors
- double sealing cap
- standard connectors



**Y CONNECTOR**

- made of plastic
- 2 x I.D. 15.0 mm
- 1 x O.D. 15.0 mm



**NECK BAND**

- adjustable (only Tracheopart Ref. 116400/116401)

\* Shaft length

# TRACHEOPART

## CLEAR DOUBLE-LUMEN TRACHEOSTOMY TUBE MADE OF PVC

semi-seated connectors, valves for Luer and Luer-lock syringes, continuous X-ray marker and additional cuff markers, colourless tracheal low-pressure cuff with colourless pilot balloon, blue bronchial low-pressure cuff with blue pilot balloon, fixation flange


- latex-free
- sterile




### TRACHEOPART TUBE SET

**RÜSCH**

| REF.         | ORDER SIZE                           | 2 X I.D. | O.D.    | CUFF Ø BRONCH. | CUFF Ø TRACH. | LENGTH* | INTRA. LENGTH** | QTY |
|--------------|--------------------------------------|----------|---------|----------------|---------------|---------|-----------------|-----|
| 116400 (set) | for left-sided bronchial intubation  |          |         |                |               |         |                 | 1   |
|              | 7.5                                  | 4.76 mm  | 13.9 mm | 18 mm          | 26 mm         | 180 mm  | 75 mm           |     |
|              | 8.5                                  | 4.76 mm  | 13.9 mm | 18 mm          | 31 ± 4.6 mm   | 190 mm  | 85 mm           |     |
|              | 9.5                                  | 4.76 mm  | 13.9 mm | 18 mm          | 31 ± 4.6 mm   | 200 mm  | 95 mm           |     |
| 116401 (set) | for right-sided bronchial intubation |          |         |                |               |         |                 | 1   |
|              | 7.5                                  | 4.76 mm  | 13.9 mm | 18 mm          | 26 mm         | 190 mm  | 75 mm           |     |
|              | 8.5                                  | 4.76 mm  | 13.9 mm | 18 mm          | 31 ± 4.6 mm   | 200 mm  | 85 mm           |     |
|              | 9.5                                  | 4.76 mm  | 13.9 mm | 18 mm          | 31 ± 4.6 mm   | 210 mm  | 95 mm           |     |



- SET COMPONENTS**
- double-lumen tracheostomy tube made of PVC, left
  - neck band
  - 2 suction catheters
  - 2 angled connectors
  - Y connector




- SET COMPONENTS**
- double-lumen tracheostomy tube made of PVC, right
  - neck band
  - 2 suction catheters
  - 2 angled connectors
  - Y connector

\* Shaft length      \*\* Intratracheal limb length


# SOFT RUBBER BRONCHIAL TUBES

## BRONCHIAL TUBE MADE OF RED SOFT RUBBER

universal injection ports, radiopaque markers, 2 cuffs and 2 pilot balloons made of SILKOLATEX® (red soft rubber coated with SILKOLATEX®)

| CARLENS TUBE SET                                                                                                                                                                           |                                     |          |           |         | RÜSCH |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------|----------|-----------|---------|-------|
| REF.                                                                                                                                                                                       | ORDER SIZE                          | 2 X I.D. | MAX. O.D. | LENGTH* | QTY   |
| 116000 (set)                                                                                                                                                                               | for left-sided bronchial intubation |          |           |         | 1     |
|                                                                                                                                                                                            | 35                                  | 4.0 mm   | 11.5 mm   | 350 mm  |       |
|                                                                                                                                                                                            | 37                                  | 4.7 mm   | 13.4 mm   | 365 mm  |       |
|                                                                                                                                                                                            | 39                                  | 4.8 mm   | 14.3 mm   | 385 mm  |       |
|                                                                                                                                                                                            | 41                                  | 5.5 mm   | 15.5 mm   | 388 mm  |       |
| SET COMPONENTS<br>• Carlens tube made of red soft rubber, left, with Carina hook<br>• suction catheter  |                                     |          |           |         |       |



| WHITE TUBE SET                                                                                                                                                                            |                                      |          |           |         | RÜSCH |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|----------|-----------|---------|-------|
| REF.                                                                                                                                                                                      | ORDER SIZE                           | 2 X I.D. | MAX. O.D. | LENGTH* | QTY   |
| 116300 (set)                                                                                                                                                                              | for right-sided bronchial intubation |          |           |         | 1     |
|                                                                                                                                                                                           | 35                                   | 4.0 mm   | 11.5 mm   | 350 mm  |       |
|                                                                                                                                                                                           | 37                                   | 4.7 mm   | 13.4 mm   | 365 mm  |       |
|                                                                                                                                                                                           | 39                                   | 4.8 mm   | 14.3 mm   | 385 mm  |       |
|                                                                                                                                                                                           | 41                                   | 5.5 mm   | 15.5 mm   | 388 mm  |       |
| SET COMPONENTS<br>• White tube made of red soft rubber, right, with Carina hook<br>• suction catheter  |                                      |          |           |         |       |



\* Shaft length

### Attention:

These products contain natural rubber latex which may cause allergic reactions.

# EZ-BLOCKER™ SET

## Y-SHAPED BRONCHUS BLOCKER FOR LUNGISOLATION AND ONE-LUNG VENTILATION

- for surgery with fiberoptic or video bronchoscopes
- 75 cm long, exact positioning via an endotracheal tube
- bifurcated distal extensions, both ends with an inflatable cuff and a central lumen, colour coded, minimal risk of dislocation during procedure
- radiopaque shaft with depth markers, cuffs made of polyurethane, providing excellent seal for up to 8 hours
- EZ-Multiport adapter to facilitate ventilation, EZ-Blocker placement and the introduction of fiberoptic or video bronchoscopes and suction catheters
- latex-free
- sterile



### EZ-BLOCKER SET

**RÜSCH**

| REF.         | SET COMPONENTS                                                                                                                                                         | LENGTH | SIZE | QTY |
|--------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|------|-----|
| MG-02770-002 | <ul style="list-style-type: none"> <li>• 1 EZ-Blocker™</li> <li>• 1 EZ-Multiport™ adaptor with scope lid</li> <li>• 1 oxygen adaptor</li> <li>• 2 dust caps</li> </ul> | 75 cm  | 7 Fr | 1   |



*EZ-Blocker combined with ETT (not included in set)*

# BRONCHUS BLOCKER

## MADE OF POLYURETHANE

detachable Luer-lock adapter with valve and stopcock, conical tip with central opening, radiopaque, balloon made of SILKOLATEX® or WIRUPREN®, double-lumen, approx. 170 cm long, graduated

- sterile

### Ref. 330600/330601

- for surgery with flexible fibre bronchoscopes with operating channel from I.D. 2.8 mm up
- ideal for use with BronchoFlex (Ref. 1041000, see page 14)

#### Attention:

These products contain natural rubber latex which may cause allergic reactions.

### Ref. 330602

curved tip for exact positioning in the left or right bronchus via the in situ tracheal tube, high-volume low-pressure cuff, premounted silicone cap for angled connector

- latex-free
- sterile



## BRONCHUS BLOCKER SET

**RÜSCH**

| REF.                                                                                                           | ORDER SIZE/O.D. | BALLOON CAPACITY | QTY |
|----------------------------------------------------------------------------------------------------------------|-----------------|------------------|-----|
| 330600 (set)                                                                                                   | Ch. 6           | max. 3.0 ml      | 1   |
| SET COMPONENTS<br>• bronchus blocker made of polyurethane<br>• single-use syringe 3.0 ml                       |                 |                  |     |
| 330601 (set)                                                                                                   | Ch. 6           | max. 5.0 ml      | 1   |
| SET COMPONENTS<br>• bronchus blocker made of polyurethane<br>• angled connector<br>• single-use syringe 5.0 ml |                 |                  |     |



## BRONCHUS BLOCKER SET

**RÜSCH**

| REF.                                                                                                           | ORDER SIZE/O.D. | BALLOON CAPACITY | QTY |
|----------------------------------------------------------------------------------------------------------------|-----------------|------------------|-----|
| 330602 (set)                                                                                                   | Ch. 6           | max. 5.0 ml      | 1   |
| SET COMPONENTS<br>• bronchus blocker made of polyurethane<br>• angled connector<br>• single-use syringe 5.0 ml |                 |                  |     |



#### ANGLED CONNECTOR

- made of plastic
- transparent
- rotating connectors
- double sealing cap with closeable 1.8 mm opening



#### SINGLE-USE SYRINGE

- 5.0 ml or
- 3.0 ml

# SHERIDAN BRONCHIAL TUBES

## SHER-I-BRONCH®

### LEFT-SIDED DOUBLE LUMEN ENDOBRONCHIAL TUBE FOR INTUBATION OF LEFT MAINSTEM BRONCHUS

colour coding distinguishes bronchial from tracheal airways; radiopaque cuff attachment points and tube tip, full-scale depth marks, blue X-ray line for ease of visualisation during fiberoptic scope insertion, tear-resistant tracheal cuff

- sterile, individually packed



SHER-I-BRONCH®

SHERIDAN

| REF.    | ORDER SIZE | 2 X I.D. | QTY |
|---------|------------|----------|-----|
| 5-16028 | 28         | 3.78     | 1   |
| 5-16035 | 35         | 4.90     |     |
| 5-16037 | 37         | 5.16     |     |
| 5-16039 | 39         | 5.51     |     |
| 5-16041 | 41         | 5.56     |     |

### RIGHT-SIDED DOUBLE LUMEN ENDOBRONCHIAL TUBE FOR INTUBATION OF RIGHT MAINSTEM BRONCHUS

1 large eye (for ventilating right upper lobe bronchial orifice)

unique double cuff design increases the margin of safety, when positioning the tube in the right mainstem bronchus; full-scale depth marks, blue X-ray line for ease of visualisation during fiberoptic scope insertion, tear-resistant tracheal cuff

- sterile, individually packed



SHER-I-BRONCH®

SHERIDAN

| REF.    | ORDER SIZE | 2 X I.D. | QTY |
|---------|------------|----------|-----|
| 5-16128 | 28         | 3.78     | 1   |
| 5-16135 | 35         | 4.90     |     |
| 5-16137 | 37         | 5.16     |     |
| 5-16139 | 39         | 5.51     |     |
| 5-16141 | 41         | 5.56     |     |



SHER-I-BRONCH® ACCESSORY PACK\*

SHERIDAN

| REF.    | DESCRIPTION                                                                                                                                        | QTY |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------|-----|
| 5-16142 | SET COMPONENTS:<br>SHER-I-SWIV/FO® double swivel connectors for use with fiberoptic scope, Y connector, and three special length suction catheters | 20  |

\* packaged separately with every SHER-I-BRONCH® endobronchial tube



# OROPHARYNGEAL/ NASOPHARYNGEAL AIRWAYS

For quick and safe airway patency without intubation the user can choose from an extensive range of RÜSCH oropharyngeal and nasopharyngeal airways.

All versions have been specifically developed to meet the customer's requirements and thus represent a safe and adequate solution for most indications. RÜSCH offers not only




standard versions made from a variety of materials but also special models for endoscopy that offer simultaneous fixation of the tracheal tube. Thanks to the broad selection, these RÜSCH products have become indispensable in many situations.

## OROPHARYNGEAL AIRWAYS

| OROPHARYNGEAL AIRWAY                                                                                                                                                                                                              |        | RÜSCH                                                                                                                                                                                                                       |                   |             |        |    |            |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|-------------|--------|----|------------|
| DESCRIPTION                                                                                                                                                                                                                       | REF.   | ORDER SIZE                                                                                                                                                                                                                  | ISO SIZE          | COLOUR CODE | QTY    |    |            |
|  <ul style="list-style-type: none"> <li>made of PVC</li> <li>Guedel airway</li> <li>colour-coded</li> <li>clear</li> <li>latex-free</li> </ul> | 124700 | 000                                                                                                                                                                                                                         | 3                 | colourless  | 10     |    |            |
|                                                                                                                                                                                                                                   |        | 00                                                                                                                                                                                                                          | 4                 | pink        |        |    |            |
|                                                                                                                                                                                                                                   |        | 0                                                                                                                                                                                                                           | 5                 | blue        |        |    |            |
|                                                                                                                                                                                                                                   |        | 1                                                                                                                                                                                                                           | 6                 | black       |        |    |            |
|                                                                                                                                                                                                                                   |        | 2                                                                                                                                                                                                                           | 7                 | white       |        |    |            |
|                                                                                                                                                                                                                                   |        | 3                                                                                                                                                                                                                           | 8                 | green       |        |    |            |
|                                                                                                                                                                                                                                   |        | 4                                                                                                                                                                                                                           | 9                 | yellow      |        |    |            |
|                                                                                                                                                                                                                                   |        | 5                                                                                                                                                                                                                           | 10                | red         |        |    |            |
|                                                                                                                                                                                                                                   |        | 6                                                                                                                                                                                                                           | 12                | purple      |        |    |            |
|                                                                                                                                                                                                                                   |        |  <ul style="list-style-type: none"> <li>made of polyethylene</li> <li>Guedel airway</li> <li>colour-coded</li> <li>latex-free</li> </ul> | 124900<br>sterile | 000         |        | 3  | colourless |
| 00                                                                                                                                                                                                                                | 4      |                                                                                                                                                                                                                             |                   | pink        |        |    |            |
| 0                                                                                                                                                                                                                                 | 5      |                                                                                                                                                                                                                             |                   | blue        |        |    |            |
| 1                                                                                                                                                                                                                                 | 6      |                                                                                                                                                                                                                             |                   | black       |        |    |            |
| 2                                                                                                                                                                                                                                 | 7      |                                                                                                                                                                                                                             |                   | white       |        |    |            |
| 3                                                                                                                                                                                                                                 | 8      |                                                                                                                                                                                                                             |                   | green       |        |    |            |
| 4                                                                                                                                                                                                                                 | 9      |                                                                                                                                                                                                                             |                   | yellow      |        |    |            |
| 5                                                                                                                                                                                                                                 | 10     |                                                                                                                                                                                                                             |                   | red         |        |    |            |
| 672901<br>non-sterile                                                                                                                                                                                                             | 00     |                                                                                                                                                                                                                             |                   | 4           | pink   | 50 |            |
|                                                                                                                                                                                                                                   | 0      |                                                                                                                                                                                                                             |                   | 5           | blue   |    |            |
|                                                                                                                                                                                                                                   | 1      |                                                                                                                                                                                                                             | 6                 | black       |        |    |            |
|                                                                                                                                                                                                                                   | 2      |                                                                                                                                                                                                                             | 7                 | white       |        |    |            |
|                                                                                                                                                                                                                                   | 3      |                                                                                                                                                                                                                             | 8                 | green       |        |    |            |
|                                                                                                                                                                                                                                   |        |                                                                                                                                                                                                                             | 4                 | 9           | yellow |    |            |
|                                                                                                                                                                                                                                   |        |                                                                                                                                                                                                                             | 5                 | 10          | red    |    |            |

OROPHARYNGEAL AIRWAYS




RUSCH

|                                                                                                                         | DESCRIPTION                                                                                                                                 | REF.             | ORDER SIZE | ISO SIZE |                | QTY |
|-------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|------------------|------------|----------|----------------|-----|
|                                        | made of soft rubber for fiberoptic intubation<br><br>• Guedel airway<br>• slit along upper surface                                          | 124400*          | 3          | 8        |                | 2   |
|                                                                                                                         |                                                                                                                                             |                  | 4          | 9        |                |     |
|                                                                                                                         |                                                                                                                                             |                  | 5          | 10       |                |     |
|  <p>Ref. 124500</p> <p>Ref. 124501</p> | made of soft rubber<br><br>• Guedel airway                                                                                                  | 124500*<br>black | 000        | 3        |                | 10  |
|                                                                                                                         |                                                                                                                                             |                  | 00         | 4        |                |     |
|                                                                                                                         |                                                                                                                                             |                  | 0          | 5        |                |     |
|                                                                                                                         |                                                                                                                                             | 124501*<br>red   | 1          | 6        |                |     |
|                                                                                                                         |                                                                                                                                             |                  | 2          | 6.5      |                |     |
|                                                                                                                         |                                                                                                                                             |                  | 3          | 8        |                |     |
|                                                                                                                         |                                                                                                                                             |                  | 4          | 9        |                |     |
|                                                                                                                         |                                                                                                                                             |                  | 5          | 10       |                |     |
|                                                                                                                         |                                                                                                                                             |                  | 6          | 12       |                |     |
|                                                                                                                         | DESCRIPTION                                                                                                                                 | REF.             | ORDER SIZE | I.D.     | LENGTH APPROX. | QTY |
|                                      | OPTOSAFE bite block made of PVC, for fiberoptic intubation<br><br>• Guedel airway<br>• circular inner diameter<br>• latex-free<br>• sterile | 124200           | 7          | 7 mm     | 67 mm          | 2   |
|                                                                                                                         |                                                                                                                                             |                  | 11         | 11 mm    | 80 mm          |     |
|                                                                                                                         |                                                                                                                                             |                  | 13         | 13 mm    | 80 mm          |     |
|                                                                                                                         |                                                                                                                                             |                  | 15         | 15 mm    | 95 mm          |     |

**\*Attention**  
 These products contain natural rubber latex which may cause allergic reactions.




CATH-GUIDE® AIRWAY

HUDSON RCI

|                                                                                     | DESCRIPTION                                                                                                                                    | REF.                                                                                                                                                                    | ORDER SIZE | I.D. | QTY    |    |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------|--------|----|
|  | • flexible vinyl with a rigid bite block<br>• Guedel style with three internal channels<br>• available in eight sizes<br>• individually packed |   | 41165      | 5    | 120 mm | 48 |
|                                                                                     |                                                                                                                                                |                                                                                                                                                                         | 41166      | 4    | 110 mm |    |
|                                                                                     |                                                                                                                                                |                                                                                                                                                                         | 41167      | 3    | 100 mm |    |
|                                                                                     |                                                                                                                                                |                                                                                                                                                                         | 41164      | 2    | 90 mm  |    |
|                                                                                     |                                                                                                                                                |                                                                                                                                                                         | 41168      | 1    | 80 mm  |    |
|                                                                                     |                                                                                                                                                |                                                                                                                                                                         | 41169      | 0    | 70 mm  |    |
|                                                                                     |                                                                                                                                                |                                                                                                                                                                         | 41170      | 00   | 60 mm  |    |
|                                                                                     |                                                                                                                                                |                                                                                                                                                                         | 41171      | 000  | 55 mm  |    |

BITEGARD® ORAL BITE BLOCK

HUDSON RCI

|                                                                                     | DESCRIPTION                                                                                                                                                                                      | REF.  | QTY |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|-----|
|  | • oral bite block<br><br>  | 41140 | 2   |



# NASOPHARYNGEAL AIRWAYS

## NASOPHARYNGEAL AIRWAY, WENDL PATTERN, MADE OF WIRUPREN®

- with adjustable flange
- latex-free
- sterile



## NASOPHARYNGEAL AIRWAY, WENDL PATTERN, MADE OF SOFT RUBBER

- with adjustable flange



**Attention:**  
These products contain natural rubber latex which may cause allergic reactions.




## NASOPHARYNGEAL AIRWAY MADE OF PVC

- transparent
- sterile



### NASOPHARYNGEAL AIRWAYS

**RÜSCH**

| REF.                                                                                          | ORDER SIZE | I.D.     | O.D.     | LENGTH | QTY |
|-----------------------------------------------------------------------------------------------|------------|----------|----------|--------|-----|
| 185420<br>   | 12         | 2.0 mm   | 4.00 mm  | 95 mm  | 5   |
|                                                                                               | 14         | 2.5 mm   | 4.67 mm  | 95 mm  |     |
|                                                                                               | 16         | 3.4 mm   | 5.33 mm  | 95 mm  |     |
|                                                                                               | 18         | 3.5 mm   | 6.00 mm  | 95 mm  |     |
|                                                                                               | 20         | 4.5 mm   | 6.67 mm  | 125 mm |     |
|                                                                                               | 22         | 5.0 mm   | 7.33 mm  | 170 mm |     |
|                                                                                               | 24         | 5.5 mm   | 8.00 mm  | 170 mm |     |
|                                                                                               | 26         | 6.0 mm   | 8.67 mm  | 170 mm |     |
|                                                                                               | 28         | 6.5 mm   | 9.33 mm  | 170 mm |     |
|                                                                                               | 30         | 7.0 mm   | 10.00 mm | 170 mm |     |
|                                                                                               | 32         | 7.5 mm   | 10.67 mm | 170 mm |     |
|                                                                                               | 34         | 8.0 mm   | 11.33 mm | 170 mm |     |
|                                                                                               | 36         | 8.5 mm   | 12.00 mm | 170 mm |     |
| 125200<br>  | 12         | 2.0 mm   | 4.0 mm   | 95 mm  | 5   |
|                                                                                               | 14         | 2.5 mm   | 4.7 mm   | 95 mm  |     |
|                                                                                               | 16         | 3.0 mm   | 5.3 mm   | 95 mm  |     |
|                                                                                               | 18         | 3.5 mm   | 6.0 mm   | 95 mm  |     |
|                                                                                               | 20         | 4.0 mm   | 6.7 mm   | 115 mm |     |
|                                                                                               | 22         | 4.5 mm   | 7.3 mm   | 125 mm |     |
|                                                                                               | 24         | 5.0 mm   | 8.0 mm   | 170 mm |     |
|                                                                                               | 26         | 5.5 mm   | 8.7 mm   | 170 mm |     |
|                                                                                               | 28         | 6.0 mm   | 9.3 mm   | 170 mm |     |
|                                                                                               | 30         | 6.5 mm   | 10.0 mm  | 170 mm |     |
|                                                                                               | 32         | 7.0 mm   | 10.7 mm  | 170 mm |     |
|                                                                                               | 34         | 7.5 mm   | 11.3 mm  | 170 mm |     |
|                                                                                               | 36         | 8.0 mm   | 12.0 mm  | 170 mm |     |
| 125410<br> | 12         | 2.5 mm   | 3.97 mm  | 60 mm  | 10  |
|                                                                                               | 14         | 3.0 mm   | 4.64 mm  | 72 mm  |     |
|                                                                                               | 16         | 3.5 mm   | 5.30 mm  | 85 mm  |     |
|                                                                                               | 18         | 4.0 mm   | 5.97 mm  | 105 mm |     |
|                                                                                               | 19         | 4.5 mm   | 6.29 mm  | 110 mm |     |
|                                                                                               | 20         | 5.0 mm   | 6.63 mm  | 115 mm |     |
|                                                                                               | 22         | 5.5 mm   | 7.29 mm  | 125 mm |     |
|                                                                                               | 24         | 6.0 mm   | 7.96 mm  | 130 mm |     |
|                                                                                               | 26         | 6.5 mm   | 8.63 mm  | 140 mm |     |
|                                                                                               | 28         | 7.0 mm   | 9.29 mm  | 155 mm |     |
|                                                                                               | 30         | 7.5 mm   | 9.96 mm  | 165 mm |     |
|                                                                                               | 32         | 8.0 mm   | 10.63 mm | 170 mm |     |
|                                                                                               | 34         | 8.5 mm   | 11.29 mm | 175 mm |     |
| 36                                                                                            | 9.0 mm     | 11.95 mm | 180 mm   |        |     |

**NASOPHARYNGEAL AIRWAY,  
MADE OF SILKOLATEX®  
RÜSCH GOLD**

- with adjustable flange and widened end
- sterile



**NASOPHARYNGEAL AIRWAY,  
ROBERTAZZI PATTERN,  
MADE OF WIRUPREN®**

- latex-free
- sterile



**NASOPHARYNGEAL AIRWAY,  
WENDL PATTERN, MADE OF  
SOFT RUBBER**

- double lumen
- with adjustable flange
- oxygen insufflation tube




**NASOPHARYNGEAL AIRWAYS**

**RÜSCH**

| REF.         | ORDER SIZE  | I.D.   | O.D.     | LENGTH  | QTY    |
|--------------|-------------|--------|----------|---------|--------|
| 185200*<br>⊗ | 20          | 4.5 mm | 6.7 mm   | 125 mm  | 10     |
|              | 22          | 5.0 mm | 7.3 mm   | 170 mm  |        |
|              | 24          | 5.5 mm | 8.0 mm   | 170 mm  |        |
|              | 26          | 6.0 mm | 8.7 mm   | 170 mm  |        |
|              | 28          | 6.5 mm | 9.3 mm   | 170 mm  |        |
|              | 30          | 7.0 mm | 10.0 mm  | 170 mm  |        |
|              | 32          | 7.5 mm | 10.7 mm  | 170 mm  |        |
|              | 34          | 8.0 mm | 11.3 mm  | 170 mm  |        |
|              | 36          | 8.5 mm | 12.0 mm  | 170 mm  |        |
|              | 185410<br>⊗ | 20     | 3.6 mm   | 6.67 mm |        |
| 22           |             | 4.6 mm | 7.33 mm  | 105 mm  |        |
| 24           |             | 5.2 mm | 8.00 mm  | 105 mm  |        |
| 26           |             | 5.8 mm | 8.67 mm  | 110 mm  |        |
| 28           |             | 6.5 mm | 9.33 mm  | 117 mm  |        |
| 30           |             | 7.0 mm | 10.00 mm | 118 mm  |        |
| 32           |             | 7.6 mm | 10.67 mm | 140 mm  |        |
| 34           |             | 7.9 mm | 11.33 mm | 155 mm  |        |
| 36           |             | 9.0 mm | 12.00 mm | 170 mm  |        |
| 125600*<br>⊗ |             | 24     | 5.5 mm   | 8.0 mm  | 160 mm |
|              | 26          | 6.0 mm | 8.7 mm   | 160 mm  |        |
|              | 28          | 6.5 mm | 9.3 mm   | 170 mm  |        |
|              | 30          | 7.0 mm | 10.0 mm  | 170 mm  |        |
|              | 32          | 7.6 mm | 10.7 mm  | 170 mm  |        |

**\*Attention**

*These products contain natural rubber latex which may cause allergic reactions.*



**STANDARD CONNECTOR**

**RÜSCH**

| REF.   | ORDER SIZE/I.D.  | DESCRIPTION                                          | QTY |
|--------|------------------|------------------------------------------------------|-----|
| 501003 | 1.5 mm – 11.0 mm | made of plastic, latex-free, 15 mm connection (O.D.) | 10  |

# ANAESTHESIA MASKS/ REBREATHING BAGS

We have a long tradition of serving the medical market with tried and tested quality products, and our RÜSCH and HUDSON RCI anaesthesia masks and rebreathing bags are classic examples of these. Whether they are made of sili-

cone, Wirupren or rubber latex, the selection of masks and bags at your disposal meets every requirement. You and your patients can rely on us.

## REUSABLE FACE MASKS

### SILICONE FACE MASK

- blue
- inflatable
- 22 mm connection
- latex-free



### SILICONE FACE MASK *RÜSCH*

| REF.   | ORDER SIZE | QTY |
|--------|------------|-----|
| 154200 | 1          | 1   |
|        | 2          |     |
|        | 3          |     |

### PAEDIATRIC SILICONE FACE MASK

- Rendall-Baker-Soucek pattern
- transparent
- 22 mm connection
- latex-free



### PAEDIATRIC SILICONE FACE MASK *RÜSCH*

| REF.   | ORDER SIZE | QTY |
|--------|------------|-----|
| 154300 | 0          | 2   |
|        | 1          |     |
|        | 2          |     |
|        | 3          |     |

### PAEDIATRIC BLACK RUBBER FACE MASK

- Rendall-Baker-Soucek pattern
- black
- 22 mm connection
- size 3 available only with harness ring



### PAEDIATRIC BLACK RUBBER FACE MASK *RÜSCH*

| REF.    | ORDER SIZE | QTY |
|---------|------------|-----|
| 154700* | 0          | 1   |
|         | 1          |     |
|         | 2          |     |
|         | 3          |     |

\* **Attention:** These products contain natural rubber latex which may cause allergic reactions.

**SILICONE FACE MASK**

- Size 0-2: round, transparent
- Size 3-5: anatomically shaped, non-inflatable
- 22 mm connection
- latex-free



**SILICONE FACE MASK** *RÜSCH*

| REF.   | ORDER SIZE | QTY |
|--------|------------|-----|
| 154900 | 0          | 2   |
|        | 1          |     |
|        | 2          |     |
|        | 3          |     |
|        | 4          |     |
|        | 5          |     |

**BLACK RUBBER FACE MASK**

- 154600 without hook ring
- 154601 with hook ring
- inflatable
- 22 mm connection



**BLACK RUBBER FACE MASK** *RÜSCH*

| REF.              | ORDER SIZE | QTY |
|-------------------|------------|-----|
| 154600*           | 0          | 1   |
| 154601* hook ring | 1          |     |
|                   | 2          |     |
|                   | 3          |     |
|                   | 4          |     |
|                   | 5          |     |
|                   | 6          |     |

\* **Attention:** These products contain natural rubber latex which may cause allergic reactions.

**SILICONE FACE MASK**

- with transparent visual field made of PSU
- inflatable pad
- latex-free



**SILICONE FACE MASK** *RÜSCH*

| REF.   | ORDER SIZE | QTY |
|--------|------------|-----|
| 672910 | 2          | 1   |
|        | 3          |     |
|        | 5          |     |

# DISPOSABLE FACE MASKS



## AIR CUSHION FACE MASK

- single use
- latex-free



## CLEAR COMFORT MASK

- sizes 1 to 3: scented
- sizes 3 to 6: 22 mm colour-coded retaining ring, removable for hand-held procedures
- 15 mm connector for sizes 1 & 2
- 22 mm connector for sizes 3 to 6
- single use



### AIR CUSHION FACE MASKS

**HUDSON RCI**

| REF.  | DESCRIPTION        | CONNECTOR            | RING COLOUR | QTY |        |
|-------|--------------------|----------------------|-------------|-----|--------|
| 41271 | neonate            | 15 mm                | white       | 20  |        |
| 41272 | infant             | 15 mm                | pink        |     |        |
| 41273 | paediatric/toddler | 22 mm                | yellow      |     |        |
| 41274 | small adult/youth  | 22 mm                | green       |     |        |
| 41275 | medium adult       | 22 mm                | red         |     |        |
| 41276 | large adult        | 22 mm                | blue        |     |        |
| 41277 | neonate            | with inflation valve | 15 mm       |     | white  |
| 41278 | infant             | with inflation valve | 15 mm       |     | pink   |
| 41279 | paediatric/toddler | with inflation valve | 22 mm       |     | yellow |
| 41280 | small adult/youth  | with inflation valve | 22 mm       |     | green  |
| 41281 | medium adult       | with inflation valve | 22 mm       |     | red    |
| 41282 | large adult        | with inflation valve | 22 mm       |     | blue   |



### CLEAR COMFORT MASK WITH INFLATION PORT

**RÜSCH**

| REF.   | ORDER SIZE | DESCRIPTION  | QTY |
|--------|------------|--------------|-----|
| 415800 | 1          | neonate      | 20  |
| 415802 | 2          | infant       |     |
| 415804 | 3          | child        |     |
| 415806 | 4          | small adult  |     |
| 415808 | 5          | medium adult |     |
| 415810 | 6          | large adult  |     |

## RUBBER BAND TO HOLD FACE MASKS

- black
- conductive

### RUBBER BAND

**RÜSCH**

| REF.    | QTY |
|---------|-----|
| 154500* | 10  |

**\*Attention**

These products contain natural rubber latex which may cause allergic reactions.



# RESUSCITATORS

## DISPOSABLE LIFESAVER™ RESUSCITATORS

Disposable lifesaver resuscitators from HUDSON RCI combine quality, features, and performance with the convenience of a disposable product. They are designed to meet or exceed ISO and ASTM standards. The unique right-angle oxygen tubing connector swivels 360°. The paediatric and neonate bags have a built-in pressure relief valve and a pressure monitoring port.

- with mask and flow diverter
- individually packed



| RESUSCITATORS |                    | <i>HUDSON RCI</i> |
|---------------|--------------------|-------------------|
| REF.          | DESCRIPTION        | QTY               |
| 45372         | adult              |                   |
| 45367         | paediatric (child) |                   |
| 45362         | neonate (infant)   |                   |





| PEEP VALVE |                                                                                                                                                                                                                                                               | <i>HUDSON RCI</i> |
|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|
| REF.       | DESCRIPTION                                                                                                                                                                                                                                                   | QTY               |
| 45385      |  <ul style="list-style-type: none"> <li>• adjustable from 5 to 20 cm H<sub>2</sub>O</li> <li>• ISO standard 30 mm I.D. connector</li> <li>• individually packed</li> </ul> | 10                |



# REBREATHING BAGS




## SILKOBAG

**RUSCH**

|                                                                                   | REF.   | ORDER SIZE / VOLUME |       | DESCRIPTION                                                                                                                         | QTY |
|-----------------------------------------------------------------------------------|--------|---------------------|-------|-------------------------------------------------------------------------------------------------------------------------------------|-----|
|  | 150700 | 0.5 l               | 2.0 l | <ul style="list-style-type: none"> <li>made of silicone</li> <li>with loop</li> <li>22 mm connection</li> <li>latex-free</li> </ul> | 2   |
|                                                                                   |        | 1.0 l               | 2.3 l |                                                                                                                                     |     |
|                                                                                   |        | 1.5 l               | 3.0 l |                                                                                                                                     |     |
|  | 150560 | 0.5 l               |       | <ul style="list-style-type: none"> <li>with loop</li> <li>with hole</li> <li>22 mm connection</li> <li>latex-free</li> </ul>        | 2   |
|                                                                                   |        | 1.5 l               |       |                                                                                                                                     |     |
|                                                                                   |        | 2.3 l               |       |                                                                                                                                     |     |


## LATEX REBREATHING BAG

**RUSCH**

|                                                                                     | REF.    | ORDER SIZE / VOLUME |       | DESCRIPTION                                                                                                                             | QTY |
|-------------------------------------------------------------------------------------|---------|---------------------|-------|-----------------------------------------------------------------------------------------------------------------------------------------|-----|
|  | 151000* | 0.5 l               | 2.3 l | <ul style="list-style-type: none"> <li>black</li> <li>with loop</li> <li>22 mm connection</li> </ul>                                    | 5   |
|                                                                                     |         | 1.0 l               | 3.0 l |                                                                                                                                         |     |
|                                                                                     |         | 1.5 l               | 3.5 l |                                                                                                                                         |     |
|                                                                                     |         | 2.0 l               | 5.0 l |                                                                                                                                         |     |
|  | 151300* | 0.5 l               | 2.3 l | <ul style="list-style-type: none"> <li>black</li> <li>with extended loop</li> <li>22 mm connection</li> </ul>                           | 5   |
|                                                                                     |         | 0.75 l              | 3.0 l |                                                                                                                                         |     |
|                                                                                     |         | 1.0 l               | 4.0 l |                                                                                                                                         |     |
|                                                                                     |         | 1.5 l               | 5.0 l |                                                                                                                                         |     |
|  | 153000* | 2.0 l               | 2.3 l | <ul style="list-style-type: none"> <li>black</li> <li>with extended closed tube</li> <li>with loop</li> <li>22 mm connection</li> </ul> | 5   |
|                                                                                     |         | 0.5 l               | 3.0 l |                                                                                                                                         |     |
|                                                                                     |         | 1.0 l               | 3.5 l |                                                                                                                                         |     |
|                                                                                     |         | 1.5 l               | 5.0 l |                                                                                                                                         |     |

## NEOPRENE REBREATHING BAG

**RUSCH**

|                                                                                     | REF.   | ORDER SIZE / VOLUME |  | DESCRIPTION                                                                                           | QTY |
|-------------------------------------------------------------------------------------|--------|---------------------|--|-------------------------------------------------------------------------------------------------------|-----|
|  | 191670 | 0.5 l               |  | <ul style="list-style-type: none"> <li>green</li> <li>22 mm connection</li> <li>latex-free</li> </ul> | 10  |
|                                                                                     |        | 1.0 l               |  |                                                                                                       |     |
|                                                                                     |        | 2.0 l               |  |                                                                                                       |     |
|                                                                                     |        | 3.0 l               |  |                                                                                                       |     |

**\*Attention**

These products contain natural rubber latex which may cause allergic reactions.

# STYLETS

## FLEXISLIP INTUBATION STYLET

flexible metal stylet with plastic coating, smooth surface, curved handle

- for tracheal tubes
- cupped, soft tip
- latex-free
- sterile



FLEXISLIP INTUBATION STYLET

| REF.   | ORDER SIZE | LENGTH CURVED | O.D.    | FOR I.D. | QTY |
|--------|------------|---------------|---------|----------|-----|
| 503700 | 6          | 220 mm        | 1.95 mm | 2.5 mm   | 20  |
|        | 10         | 335 mm        | 3.28 mm | 3.5 mm   |     |
|        | 12         | 366 mm        | 3.93 mm | 4.0 mm   |     |
|        | 14         | 366 mm        | 4.60 mm | 5.0 mm   |     |

## INTUBATION STYLET MADE OF PVC

- for tracheal tubes
- flexible
- cupped, soft tip
- latex-free



INTUBATION STYLET MADE OF PVC

RÜSCH

| REF.   | ORDER SIZE | LENGTH APPROX. | O.D.   | FOR I.D.   | QTY |
|--------|------------|----------------|--------|------------|-----|
| 503000 | 15         | 350 mm         | 1.5 mm | 2.0–2.5 mm | 10  |
|        | 20         | 400 mm         | 2.0 mm | 2.5–3.0 mm |     |
|        | 26         | 400 mm         | 2.6 mm | 3.0–3.5 mm |     |
|        | 33         | 450 mm         | 3.3 mm | 3.5–4.5 mm |     |
|        | 43         | 450 mm         | 4.3 mm | 4.5–6.0 mm |     |
|        | 56         | 500 mm         | 5.6 mm | ab 6.0 mm  |     |

- for non-kinking (Oxford) tube
- right-angled
- hard
- latex-free



INTUBATION STYLET MADE OF PVC

RÜSCH

| REF.   | ORDER SIZE | LENGTH APPROX. | O.D.   | FOR I.D.    | QTY |
|--------|------------|----------------|--------|-------------|-----|
| 502900 | 1          | 280 mm         | 2.6 mm | 3.5–5.0 mm  | 5   |
|        | 2          | 370 mm         | 4.3 mm | 5.5–6.0 mm  |     |
|        | 3          | 410 mm         | 5.6 mm | 7.0–12.0 mm |     |

## ENDOGUIDE TUBE EXCHANGER AND INSERTION AID IN ONE

permanently available large lumen for constant oxygenation of the patient; connection system for 15 mm and Luer-lock connectors, depth markers

- latex-free
- sterile



ENDOGUIDE-T

RÜSCH

| REF.   | ORDER SIZE | I.D.   | O.D.   | LENGTH* | QTY |
|--------|------------|--------|--------|---------|-----|
| 503100 | 2.5        | 1.4 mm | 2.6 mm | 700 mm  | 1   |
|        | 6.0        | 3.2 mm | 5.0 mm | 830 mm  |     |
| 503110 | 6.0        | 3.2 mm | 5.0 mm | 525 mm  |     |

\* Length without connector

### CONNECTOR SET


- Luer-lock connector
- silicone cone
- 15 mm standard connector





### SHER-I-SLIP® INTUBATION STYLET

dependable stiffness retains desired curve during oral intubation; blue colour improves visualisation; special material reduces stylet-to-tube friction; smooth moulded distal tip end



#### SHER-I-SLIP INTUBATION STYLET *HUDSON RCI*

| REF.    | LENGTH APPROX. | O.D.  | FOR TUBE SIZES I.D. | QTY |
|---------|----------------|-------|---------------------|-----|
| 5-15101 | 340 mm         | 10 Fr | 5.0 - 10.0 mm       | 25  |
| 5-15102 | 250 mm         | 5 Fr  | 2.0 - 4.5 mm        |     |
| 5-15103 | X-Long         | 10 Fr | 5.0 - 10.0 mm       |     |

### SOFT-TIP™ INTUBATION STYLET

atraumatic tip to minimise trauma; slippery outer finish makes for easy transition into endotracheal tube



#### SOFT-TIP INTUBATION STYLET *HUDSON RCI*

| REF.    | LENGTH APPROX. | O.D.  | FOR TUBE SIZES I.D. | QTY |
|---------|----------------|-------|---------------------|-----|
| 5-15120 | 270 mm         | 6 Fr  | 2.5 - 4.5 mm        | 20  |
| 5-15110 | 330 mm         | 10 Fr | 4.5 - 6.5 mm        |     |
| 5-15100 | 310 mm         | 14 Fr | 6.0 - 10.0 mm       |     |

# ACCESSORIES

## UNIVERSAL ADAPTER, MAINZ PATTERN

| REF.   | DESCRIPTION                                                                                                                                                                                                                                                                                                                                                                                                   | QTY |
|--------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
| 514800 | <p>sizes:<br/>O.D. 15.0 mm and 22.0 mm</p> <ul style="list-style-type: none"> <li>made of plastic, white</li> <li>for fiberoptic intubation during simultaneous mask ventilation</li> <li>connector for face masks</li> <li>lateral connector for anaesthetic circuit</li> <li>double sealing cap made of silicone</li> <li>latex-free</li> <li>sterile</li> </ul>                                            | 5   |
| 514805 | <p>sizes: O.D. 15.0 mm and 22.0 mm / I.D. 15 mm</p> <ul style="list-style-type: none"> <li>made of blue plastic</li> <li>for fiberoptic intubation during simultaneous mask ventilation</li> <li>connector for face masks, tracheal tubes and laryngeal masks</li> <li>lateral connector for anaesthetic circuit</li> <li>double sealing cap made of silicone</li> <li>latex-free</li> <li>sterile</li> </ul> | 5   |

## SHER-I-SWIV® DOUBLE SWIVEL CONNECTOR

| REF.    | DESCRIPTION                                                                                                                                                                                                                                                                                                                            | QTY |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
| 5-15301 | <ul style="list-style-type: none"> <li>for use with tracheostomy and tracheal tubes</li> <li>reduces breathing circuit drag and torque on the tube</li> <li>strap connector for anti-disconnect security</li> <li>special swivel joints reduce leakage potential</li> <li>easy suction catheter passage with secure closure</li> </ul> | 25  |

## SHER-I-SWIV/FO™ DOUBLE SWIVEL CONNECTOR

| REF.    | DESCRIPTION                                                                                                                                                                                                    | QTY |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
| 5-15401 | <ul style="list-style-type: none"> <li>for use with fiberoptics</li> <li>allows easy passage of fiberoptic scope while maintaining an air-tight seal around scope for positive pressure ventilation</li> </ul> | 25  |


## DOUBLE SEALING CAP MADE OF SILICONE

| REF.   | DESCRIPTION                                                                                                                                                                                                       | QTY |
|--------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
| 514801 | <p>sizes:<br/>I.D. 2.0 mm, blue<br/>I.D. 3.5 mm, transparent</p> <ul style="list-style-type: none"> <li>for universal adapter, Mainz pattern (Ref. 514800/514805)</li> <li>latex-free</li> <li>sterile</li> </ul> | 5   |


## M-CONNECTOR MADE OF PLASTIC

| REF.   | DESCRIPTION                                                                                                                                                            | QTY |
|--------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
| 514700 | <ul style="list-style-type: none"> <li>m-connector</li> <li>blue</li> <li>2 x O.D. 22 mm</li> <li>1 x O.D. 22 mm / I.D. 15 mm, rotating</li> <li>latex-free</li> </ul> | 5   |

**MOUTH GAG, MADE OF RED SOFT RUBBER**

|                                                                                                              | REF.    | DESCRIPTION           | QTY |
|--------------------------------------------------------------------------------------------------------------|---------|-----------------------|-----|
|  <p>a)      b)</p>           | 610000* | a) blunt-nosed        | 5   |
|                                                                                                              | 610001* | b) long-tapered       | 5   |
|                                                                                                              |         | size: 95 x 27 x 18 mm |     |
| <p><b>*Attention</b><br/>These products contain natural rubber latex which may cause allergic reactions.</p> |         |                       |     |


**CUFF PROTECTION TUBE**

|                                                                                                              | REF.    | DESCRIPTION                                                                                                                                                                                                                                          | QTY |
|--------------------------------------------------------------------------------------------------------------|---------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
|                           | 502501* | <ul style="list-style-type: none"> <li>• O.D. 4.0 mm</li> <li>• for nasal insertion of tracheal tubes</li> <li>• length of yellow protection sleeve approx. 7 cm</li> <li>• total length approx. 40 cm</li> <li>• made of PVC/SILKOLATEX®</li> </ul> | 5   |
| <p><b>*Attention</b><br/>These products contain natural rubber latex which may cause allergic reactions.</p> |         |                                                                                                                                                                                                                                                      |     |


**ENDOTRACHEAL TUBE HOLDER**

|                                                                                    | REF.  | DESCRIPTION                                                                                                                                                                                         | QTY |
|------------------------------------------------------------------------------------|-------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
|  | 41065 | <ul style="list-style-type: none"> <li>• keeps the endotracheal tube in place</li> <li>• supplied with cloth tie, safety-pin, and soft foam mouth cushion</li> <li>• built-in bite block</li> </ul> | 10  |


**SINGLE USE MAGILL FORCEPS**

|                                                                                    | REF.          | DESCRIPTION | QTY |
|------------------------------------------------------------------------------------|---------------|-------------|-----|
|  | 503300-000150 | infant      | 10  |
|                                                                                    | 503300-000200 | child       | 10  |
|                                                                                    | 503300-000245 | adult       | 10  |


**ENDOTEST**

|                                                                                   | REF.   | DESCRIPTION                                                                                                                                                                                                                                                                 | QTY |
|-----------------------------------------------------------------------------------|--------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
|  | 112700 | device for filling and monitoring the pressure of low-pressure cuffs on tracheal, tracheostomy, bronchial tubes and laryngeal masks <ul style="list-style-type: none"> <li>• with connecting tube</li> </ul>                                                                | 1   |
|                                                                                   | 230100 | connecting tube made of PVC for ENDOTEST cuff pressure monitoring device <ul style="list-style-type: none"> <li>• transparent</li> <li>• with Luer-connector (male and female)</li> <li>• approx. 100 cm long, 2 x 3 mm</li> <li>• latex-free</li> <li>• sterile</li> </ul> | 10  |

**SILKOSPRAY**

|                                                                                     | REF.   | DESCRIPTION                                                                                                                                                                                                                                                                                                                                                                                                      | QTY |
|-------------------------------------------------------------------------------------|--------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
|  | 556000 | <ul style="list-style-type: none"> <li>• universal silicone spray</li> <li>• prevents incrustation of rubber, latex and PVC devices and also prevents them from sticking to the mucosa</li> <li>• CFC free</li> <li>• latex-free</li> </ul> <p>Please note:<br/>SILKOSPRAY is not suitable for devices made of silicone.<br/>For silicone products we recommend our water-soluble lubricating gel TrachJell.</p> | 10  |

**TRACHJELL**


|                                                                                     | REF.   | DESCRIPTION                                                                                                                       | QTY |
|-------------------------------------------------------------------------------------|--------|-----------------------------------------------------------------------------------------------------------------------------------|-----|
|  | 556100 | <ul style="list-style-type: none"> <li>• single dose water-soluble lubricating gel</li> <li>• 8.5 g</li> <li>• sterile</li> </ul> | 25  |

## INDEX

| REF.    | PAGE | REF.    | PAGE | REF.    | PAGE | REF.    | PAGE |
|---------|------|---------|------|---------|------|---------|------|
| 5-10206 | 20   | 5-11110 | 24   | 5-16128 | 36   | 5-30404 | 21   |
| 5-10208 | 20   | 5-11112 | 24   | 5-16135 | 36   | 5-30405 | 21   |
| 5-10210 | 20   | 5-11714 | 24   | 5-16137 | 36   | 5-30406 | 21   |
| 5-10212 | 20   | 5-11716 | 24   | 5-16139 | 36   | 5-30407 | 21   |
| 5-10213 | 20   | 5-11718 | 24   | 5-16141 | 36   | 5-30408 | 21   |
| 5-10214 | 20   | 5-12510 | 23   | 5-16142 | 36   | 5-30409 | 21   |
| 5-10215 | 20   | 5-12511 | 23   | 5-20512 | 23   | 5-30410 | 21   |
| 5-10216 | 20   | 5-12512 | 23   | 5-20513 | 23   | 5-30411 | 21   |
| 5-10217 | 20   | 5-12513 | 23   | 5-20514 | 23   | 5-30412 | 21   |
| 5-10218 | 20   | 5-12514 | 23   | 5-20515 | 23   | 5-30413 | 21   |
| 5-10219 | 20   | 5-12515 | 23   | 5-20516 | 23   | 5-30414 | 21   |
| 5-10220 | 20   | 5-12516 | 23   | 5-20517 | 23   | 41065   | 49   |
| 5-10310 | 20   | 5-12517 | 23   | 5-20608 | 24   | 41140   | 38   |
| 5-10311 | 20   | 5-12518 | 23   | 5-20610 | 24   | 41164   | 38   |
| 5-10312 | 20   | 5-13012 | 19   | 5-20612 | 24   | 41165   | 38   |
| 5-10313 | 20   | 5-13013 | 19   | 5-22006 | 22   | 41166   | 38   |
| 5-10314 | 20   | 5-13014 | 19   | 5-22007 | 22   | 41167   | 38   |
| 5-10315 | 20   | 5-13015 | 19   | 5-22008 | 22   | 41168   | 38   |
| 5-10316 | 20   | 5-13016 | 19   | 5-22009 | 22   | 41169   | 38   |
| 5-10317 | 20   | 5-13017 | 19   | 5-22010 | 22   | 41170   | 38   |
| 5-10318 | 20   | 5-13018 | 19   | 5-22011 | 22   | 41171   | 38   |
| 5-10320 | 20   | 5-15100 | 47   | 5-22012 | 22   | 41271   | 43   |
| 5-10404 | 21   | 5-15101 | 47   | 5-22013 | 22   | 41272   | 43   |
| 5-10405 | 21   | 5-15102 | 47   | 5-22014 | 22   | 41273   | 43   |
| 5-10406 | 21   | 5-15103 | 47   | 5-22106 | 22   | 41274   | 43   |
| 5-10407 | 21   | 5-15110 | 47   | 5-22107 | 22   | 41275   | 43   |
| 5-10408 | 21   | 5-15120 | 47   | 5-22108 | 22   | 41276   | 43   |
| 5-10409 | 21   | 5-15301 | 48   | 5-22109 | 22   | 41277   | 43   |
| 5-10410 | 21   | 5-15401 | 48   | 5-22110 | 22   | 41278   | 43   |
| 5-10411 | 21   | 5-16028 | 36   | 5-22111 | 22   | 41279   | 43   |
| 5-10412 | 21   | 5-16035 | 36   | 5-22112 | 22   | 41280   | 43   |
| 5-10413 | 21   | 5-16037 | 36   | 5-22113 | 22   | 41281   | 43   |
| 5-10414 | 21   | 5-16039 | 36   | 5-22114 | 22   | 41282   | 43   |
| 5-11108 | 24   | 5-16041 | 36   | 5-23000 | 19   | 45362   | 44   |

| REF.   | PAGE | REF.   | PAGE | REF.   | PAGE | REF.          | PAGE |
|--------|------|--------|------|--------|------|---------------|------|
| 45367  | 44   | 111380 | 18   | 125200 | 39   | 503000        | 46   |
| 45372  | 44   | 111480 | 18   | 125410 | 39   | 503100        | 46   |
| 45385  | 44   | 111780 | 8    | 125600 | 40   | 503110        | 46   |
| 100180 | 8    | 111781 | 9    | 150560 | 45   | 503300-000150 | 49   |
| 100181 | 9    | 111782 | 7    | 150700 | 45   | 503300-000200 | 49   |
| 100280 | 10   | 112080 | 3    | 151000 | 45   | 503300-000245 | 49   |
| 100380 | 4    | 112082 | 3    | 151300 | 45   | 503700        | 46   |
| 100382 | 4    | 112460 | 16   | 153000 | 45   | 514700        | 48   |
| 100480 | 7    | 112480 | 2    | 154200 | 41   | 514800        | 48   |
| 102000 | 17   | 112481 | 3    | 154300 | 41   | 514801        | 48   |
| 102004 | 15   | 112482 | 2    | 154500 | 43   | 514805        | 48   |
| 103040 | 14   | 112662 | 5    | 154600 | 42   | 556000        | 49   |
| 103600 | 12   | 112680 | 6    | 154601 | 42   | 556100        | 26   |
| 104000 | 12   | 112682 | 6    | 154700 | 41   | 556100        | 49   |
| 104004 | 12   | 112700 | 49   | 154900 | 42   | 610000        | 49   |
| 104100 | 14   | 112880 | 16   | 185200 | 40   | 610001        | 49   |
| 104103 | 14   | 115900 | 29   | 185410 | 40   | 672901        | 37   |
| 104201 | 11   | 115901 | 29   | 185420 | 39   | 672910        | 42   |
| 104202 | 11   | 116000 | 33   | 191670 | 45   | MG-02770-002  | 35   |
| 104203 | 11   | 116100 | 30   | 230100 | 49   |               |      |
| 104204 | 11   | 116101 | 31   | 330600 | 35   |               |      |
| 105000 | 28   | 116162 | 30   | 330601 | 35   |               |      |
| 105010 | 28   | 116200 | 30   | 330602 | 35   |               |      |
| 105102 | 10   | 116201 | 31   | 415800 | 43   |               |      |
| 105200 | 26   | 116262 | 30   | 415802 | 43   |               |      |
| 105210 | 26   | 116300 | 33   | 415804 | 43   |               |      |
| 105220 | 26   | 116400 | 32   | 415806 | 43   |               |      |
| 105300 | 27   | 116401 | 32   | 415808 | 43   |               |      |
| 105310 | 27   | 124200 | 38   | 415810 | 43   |               |      |
| 105320 | 27   | 124400 | 38   | 501003 | 17   |               |      |
| 105400 | 13   | 124500 | 38   | 501003 | 40   |               |      |
| 105702 | 13   | 124501 | 38   | 502501 | 49   |               |      |
| 111000 | 28   | 124700 | 37   | 502700 | 5    |               |      |
| 111100 | 18   | 124900 | 37   | 502900 | 34   |               |      |





Central Venous Catheters  
Haemodialysis Access  
Haemodynamic Monitoring  
Trauma and Surgery  
Endotracheal Tubes & Double Lumen Tubes  
Oropharyngeal & Nasopharyngeal Airways  
Laryngoscopes  
Tracheostomy  
Regional Anaesthesia  
Aerosol & Oxygen Therapy  
Humidification  
Bacterial/Viral Filtration  
Surgical Instrumentation

ARROW | GIBECK | HUDSON  
KMEDIC | PILLING | RÜSCH | WECK

Central Venous Catheters  
PICC  
Haemodialysis Access  
Haemodynamic Monitoring  
Arterial Access  
Trauma and Surgery  
Regional Anaesthesia  
Standard Bladder Management  
Operative Urology and Endourology  
Intermittent Catheters  
Ligation  
Endo-Surgery/Laparoscopic Access  
Cardio-Thoracic Surgery & Drainage  
Surgical Instrumentation  
Intra-Aortic Balloon Pump & Catheters  
Right Heart Products  
Transradial Access

ARROW | DEKNATEL | KMEDIC  
MEDICAL SERVICE | PILLING  
PLEUR-EVAC | RÜSCH | TAUT | WECK

**READ UP ON OUR  
PRODUCT RANGES**

GO TO [WWW.TELEFLEX.COM/EMEA](http://WWW.TELEFLEX.COM/EMEA) AND CLICK ON  
SERVICES TO DOWNLOAD INFORMATION MATERIAL



Regional Anaesthesia  
Ligation  
Surgical Instrumentation

ARROW | KMEDIC | PILLING | WECK

Teleflex is a leading global provider of specialty medical devices used for diagnostic and therapeutic procedures in critical care, urology and surgery. Our mission is to provide solutions that enable healthcare providers to improve outcomes and enhance patient and provider safety.

We specialise in devices for general and regional anaesthesia, cardiac care, respiratory care, urology, vascular access and surgery and we serve healthcare providers in more than 150 countries. Teleflex also provides specialty products for medical device manufacturers.

Our well known brands include **ARROW**<sup>®</sup>, **DEKNATEL**<sup>®</sup>, **GIBECK**<sup>®</sup>, **HUDSON RCI**<sup>®</sup>, **KMEDIC**<sup>®</sup>, **LMA**<sup>™</sup>, **PILLING**<sup>®</sup>, **PLEUR-EVAC**<sup>®</sup>, **RÜSCH**<sup>®</sup>, **SHERIDAN**<sup>®</sup>, **TAUT**<sup>®</sup>, **TFX OEM**<sup>®</sup>, **VASONOVA**<sup>™</sup>, **VIDACARE**<sup>®</sup> and **WECK**<sup>®</sup>, all of which are trademarks or registered trademarks of Teleflex Incorporated.

Teleflex global operations: Australia, Austria, Belgium, Canada, China, Czech Republic, France, Germany, Greece, India, Ireland, Italy, Japan, Malaysia, Mexico, Netherlands, New Zealand, Portugal, Singapore, Slovak Republic, South Africa, Spain, Switzerland, United Kingdom, Uruguay and USA.

#### **YOUR INTERNATIONAL CONTACTS:**

##### **TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND**

Teleflex Medical Europe Ltd., IDA Business and Technology Park,  
Dublin Road, Athlone, Co Westmeath  
Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73  
orders.intl@teleflex.com

**AUSTRALIA/NEW ZEALAND** +61 (0)3 9081 0600

**AUSTRIA** +43 (0)1 402 47 72

**BELGIUM** +32 (0)2 333 24 60

**CHINA (SHANGHAI)** +86 (0)21 6163 0965

**CHINA (BEIJING)** +86 (0)10 6418 5699

**CZECH REPUBLIC** +420 (0)495 759 111

**FRANCE** +33 (0)5 62 18 79 40

**GERMANY** +49 (0)7151 406 0

**GREECE** +30 210 67 77 717

**INDIA** +91 (0)44-2836 5040

**ITALY** +39 0362 58 911

**JAPAN** +81 (0)3 3379 1511

**NETHERLANDS** +31 (0)88 00 215 00

**PORTUGAL** +351 22 541 90 85

**SINGAPORE** +65 6439 3000

**SLOVAK REPUBLIC** +421 (0)3377 254 28

**SOUTH AFRICA** +27 (0)11 807 4887

**SPAIN** +34 918 300 451

**SWITZERLAND** +41 (0)31 818 40 90

**UNITED KINGDOM** +44 (0)1494 53 27 61

For detailed information see [www.teleflex.com](http://www.teleflex.com)

The products in this catalogue may not be available in all countries.  
Please contact your local representative. All data current at time of printing (03/2014).  
Subject to technical changes without further notice.

94 02 02 - 00 00 11 · REV B · MC / WM · 03 14 01