





**The New Standard of Excellence,
Redefining Deliverability, Versatility, and Control**

 				
Product Code	Description	Balloon Diameter (mm)	Balloon Length (mm)	Distal Tip Length (mm)
1 balloon catheter, 1 introducer sheath, 1 shaping mandrel per box				
BC0410C	Compliant Occlusion Balloon 4.0mm x 10mm	4.0	10	5
BC0415C	Compliant Occlusion Balloon 4.0mm x 15mm	4.0	15	5
BC0420C	Compliant Occlusion Balloon 4.0mm x 20mm	4.0	20	5
 				
Product Code	Description	Balloon Diameter (mm)	Balloon Length (mm)	Distal Tip Length (mm)
1 balloon catheter, 1 introducer sheath, 1 shaping mandrel per box				
BC0411XC	X-tra Compliant Occlusion Balloon 4.0mm x 11mm	4.0	11	5

Features

14 Wire Compatible

Hydrophilic Coating on the Balloon

Soft & Long Tip

Low Profile

Distal Tip Marker

Benefits

Provides more choices of guidewire (Traxcess® docking wire compatible)
Improves trackability and stability with better performance of the wire

Significantly improves trackability
Enables use of longer balloons without compromising trackability

Improves balloon stability and trackability

Ensures compatibility with double catheter technique in 6F guide
Improves trackability

Positive identification of distal tip location



MicroVention, Inc.

Worldwide Headquarters
1311 Valencia Avenue
Tustin, CA 92780 USA

MicroVention UK Limited
MicroVention Europe, S.A.R.L.

MicroVention Deutschland GmbH
Web

PH +1.714.247.8000

PH +44 (0) 191 258 6777

PH +33 (1) 39 21 77 46

PH +49 211 210 798-0

microvention.com

