

Wound  
Management  
Product  
Catalogue

At **smith&nephew** our focus is on repairing and healing the human body. We are really proud of what we do.



# Product Catalogue

---

## Welcome

This comprehensive booklet gives full details of Smith & Nephew products currently available.

Comprehensive evidence is available upon request from either our head office in Hull or from your local sales representative.

Smith & Nephew also offer educational therapy information on topics such as wound bed preparation, leg ulcer care and pressure ulcer care.

The Smith & Nephew national study day programme offers places to over 5,000 nurses each year and focuses upon the latest advances in wound care.

These study days utilise eminent speakers from both within the UK and across the world. In addition, a local study day programme run by ENB18 trained staff is available.

The Smith & Nephew Formulary package has been developed to provide assistance locally to achieve clinically proven and cost-effective wound care. This package provides support on planning, organising, printing and training dependent on your specific needs and requirements.

For further details on products or any of the resources offered by Smith & Nephew, please contact your local representative or see back cover of this booklet for useful contact information.



# Product Catalogue

---

## Contents

- |    |   |    |  |
|----|---|----|--|
| 4  | ACTICOAT <sup>®</sup> /ACTICOAT 7               | 35 | ICHTHOPASTE <sup>®</sup>                     |
| 5  | ACTICOAT Absorbent                              | 36 | INTRASITE <sup>®</sup> Conformable           |
| 6  | ACTICOAT Flex 3/ACTICOAT Flex 7                 | 36 | INTRASITE Gel                                |
| 7  | ACTICOAT Moisture Control                       | 37 | IODOFLEX <sup>®</sup> /IODOSORB <sup>®</sup> |
| 8  | ACTICOAT Site                                   | 38 | IV3000 <sup>®</sup>                          |
| 9  | ALGISITE <sup>®</sup> Ag                        | 39 | JELONET <sup>®</sup>                         |
| 10 | ALGISITE M                                      | 40 | LEUKOSTRIP <sup>®</sup>                      |
| 11 | ALGOSTERIL <sup>®</sup>                         | 41 | MELOLIN <sup>®</sup>                         |
| 12 | ALLEVYN <sup>®</sup> Adhesive/Sacrum            | 42 | MELOLITE <sup>®</sup>                        |
| 13 | ALLEVYN Ag/Adhesive/Non<br>Adhesive/Heel/Sacrum | 43 | NO-STING SKIN-PREP <sup>®</sup>              |
| 14 | ALLEVYN Ag Gentle/<br>Gentle Border             | 44 | OPSITE <sup>®</sup> FLEXIFIX                 |
| 15 | ALLEVYN Cavity                                  | 45 | OPSITE FLEXIGRID <sup>®</sup>                |
| 16 | ALLEVYN Compression                             | 46 | OPSITE Incise                                |
| 17 | ALLEVYN Gentle/Gentle Border/<br>Border Heel    | 47 | OPSITE Post-Op/Plus                          |
| 18 | ALLEVYN Gentle Border Lite                      | 48 | OPSITE Post-Op Visible                       |
| 19 | ALLEVYN Heel/Ag Heel                            | 49 | OPSITE Spray                                 |
| 20 | ALLEVYN Lite                                    | 50 | PICO   |
| 21 | ALLEVYN Non-Adhesive                            | 51 | PRIMAFIX <sup>®</sup>                        |
| 22 | ALLEVYN Plus Adhesive/Sacrum                    | 52 | PRIMAPORE <sup>®</sup>                       |
| 23 | ALLEVYN Plus Cavity                             | 53 | PROFORE <sup>®</sup>                         |
| 24 | ALLEVYN Thin                                    | 54 | PROGUIDE <sup>®</sup>                        |
| 25 | ALLEVYN Tracheostomy                            | 55 | RENASYS <sup>®</sup> EZ PLUS                 |
| 26 | AMETOP <sup>®</sup> Gel                         | 56 | RENASYS GO                                   |
| 27 | BACTIGRAS <sup>®</sup>                          | 57 | TRICOTEX <sup>®</sup>                        |
| 27 | CADESORB <sup>®</sup>                           | 58 | TRIPLE CARE <sup>®</sup>                     |
| 28 | CARBONET <sup>®</sup>                           | 59 | VERSAJET <sup>®</sup>                        |
| 29 | CICA-CARE <sup>®</sup>                          | 60 | VISCOPASTE <sup>®</sup> PB7                  |
| 30 | CICA-PLAIE <sup>®</sup>                         | 61 | VISITRAK Depth                               |
| 31 | CUTINOVA <sup>®</sup> Hydro                     | 62 | VISITRAK Digital                             |
| 31 | CUTIPLAST <sup>®</sup> Steril                   | 62 | VISITRAK Grid                                |
| 32 | DURAFIBER <sup>®</sup>                          | 63 | ZIPZOC <sup>®</sup>                          |
| 33 | EXU-DRY <sup>®</sup>                            | 64 | ZOFF <sup>®</sup>                            |
| 34 | FLAMAZINE <sup>®</sup> Cream                    | 65 | Index  |
|    |   | 66 | Prescribing information                      |



# Product Catalogue

## ACTICOAT<sup>◇</sup> / ACTICOAT 7

### with SILCRYST<sup>◇</sup>

Silver-coated antimicrobial barrier dressing - sterile


#### Indications

ACTICOAT dressing is indicated as an antimicrobial barrier layer over partial and full thickness wounds such as pressure ulcers, venous leg ulcers, diabetic ulcers, burns, donor and recipient graft sites.

#### Contraindications

Do not use on patients with a known sensitivity to silver.

Do not use on patients undergoing Magnetic Resonance Imaging (MRI) examination.



S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>ACTICOAT</b>				
66000808	5cm x 5cm	5	ELY141	301-3224
66000789	10cm x 10cm	5	-	294-3942
66000791	10cm x 10cm	12	ELY071	295-5782
66000792	10cm x 20cm	12	ELY072	292-5774
66000793	20cm x 40cm	6	ELY073	301-3232
<b>ACTICOAT 7</b>				
66000809	5cm x 5cm	5	ELY143	301-3216
66000796	10cm x 12.5cm	5	ELY076	285-3372
66000797	15cm x 15cm	5	ELY077	285-3380

# ACTICOAT<sup>◇</sup> Absorbent with SILCRYST<sup>◇</sup>

Silver-coated antimicrobial absorbent dressing – sterile

## Indications

ACTICOAT Absorbent is indicated as an antimicrobial absorbent dressing for use over partial and full thickness wounds such as pressure ulcers, venous leg ulcers, diabetic ulcers, burns, cavity wounds, donor and recipient graft sites.



## Contraindications

Do not use on patients with a known sensitivity to silver or alginates.

Do not use on patients undergoing Magnetic Resonance Imaging (MRI) examination.

Do not use on dry or non-exuding wounds.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66001760	5cm x 5cm	5	ELY191	319-5021
66000841	10cm x 12.5cm	5	ELY171	304-2272
66000840	2cm x 30cm rope	5	ELY170	304-2264



ACTICOAT<sup>◇</sup> Flex 3/

## ACTICOAT Flex 7

with SILCRYST<sup>◇</sup>

Silver-coated antimicrobial barrier dressing - sterile

**Indications**

ACTICOAT Flex dressing is indicated as an antimicrobial barrier layer over partial and full thickness wounds such as pressure ulcers, venous leg ulcers, diabetic ulcers, burns, recipient graft sites and surgical sites.

**Contraindications**

Do not use on patients with a known sensitivity to silver.

Do not use on patients undergoing Magnetic Resonance Imaging (MRI) examination.

Prior to administering radiation therapy, remove ACTICOAT Flex 3 or 7. A new dressing can be applied following the treatment.

S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>ACTICOAT Flex 3</b>				
66800396	5cm x 5cm	5	ELY291	351-0781
66800398	10cm x 10cm	5	-	351-3256
66800399	10cm x 10cm	12	ELY292	351-0799
66800409	10cm x 20cm	12	ELY293	351-0807
66800419	20cm x 40cm	6	ELY295	351-0815
66800432	40cm x 40cm	6	ELY296	-
66800435	10cm x 120cm	6	ELY294	-
<b>ACTICOAT Flex 7</b>				
66800395	5cm x 5cm	5	ELY297	351-0823
66800397	10cm x 12.5cm	5	ELY298	351-0831
66800420	15cm x 15cm	5	ELY299	351-0849
66800400	20cm x 40cm	6	ELY300	-
66800401	40cm x 40cm	6	ELY301	-
66800545	2.5cm x 60cm	5	ELY370	364 - 7773

# Product Catalogue

## ACTICOAT<sup>◇</sup> Moisture Control with SILCRYST<sup>◇</sup>

Silver-coated antimicrobial barrier dressing - sterile

### Indications

ACTICOAT Moisture Control is indicated as an antimicrobial absorbent foam dressing over partial and full thickness wounds such as pressure ulcers, venous leg ulcers, diabetic ulcers, burns, epidermolysis bullosa, donor and recipient graft sites.

### Contraindications

Do not use on patients with a known sensitivity to silver.

Do not use on patients undergoing Magnetic Resonance Imaging (MRI) examination.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66001790	5cm x 5cm	10	ELY192	211-6671
66001791	10cm x 10cm	10	ELY193	211-6663
66000792	10cm x 20cm	10	ELY194	211-6655



# ACTICOAT<sup>◇</sup> Site with SILCRYST<sup>◇</sup>

Silver-coated absorbent antimicrobial dressing

## Indications

Silver-coated absorbent antimicrobial dressing for use around infected insertion sites.

## Contraindications

Do not use on patients with a known sensitivity to silver, undergoing MRI (Magnetic Resonance Imaging) examination or with a known sensitivity to polyurethane.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66001788	2.5cm disc with 4mm centre hole	10	ELY281	-



# Product Catalogue

## ALGISITE<sup>◇</sup> Ag

Calcium alginate wound dressing with silver

### Indications

ALGISITE Ag dressings are intended for use in wounds where there is a moderate to heavy level of exudate such as pressure sores, post-operative wounds, partial thickness and full thickness wounds, arterial, venous and diabetic leg ulcers, pressure sores, post-operative wounds and fungating lesions.

May be used on wounds which are infected or are likely to become infected.



### Contraindications

Do not use on patients with known sensitivity to ALGISITE Ag or any of its components. Do not use on patients undergoing Magnetic Resonance Imaging (MRI) examination.

S&N Code	Description	Pack Content	NHS Code	PIP Code
66800114	2g x 30cm Rope	5	ELA356	331-5090
66800111	5cm x 5cm	10	ELA353	331-5066
66800112	10cm x 10cm	10	ELA354	331-5082
66800113	10cm x 20cm	5	ELA355	331-5074



# Product Catalogue

---

## ALGISITE M

Calcium alginate dressing – sterile

### Indications

For the management of full and partial thickness wounds with moderate to heavy exudate and which may also be prone to minor bleeding such as pressure ulcers, leg ulcers, surgical wounds, diabetic foot ulcers, malignant wounds and burns.

### Contraindications

Do not use on patients with a known allergy to alginate dressings.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66000519	5cm x 5cm	10	ELS089	255-1117
66000520	10cm x 10cm	10	ELS090	255-1125
66000521	15cm x 20cm	10	ELS091	255-1133
66000522	2cm x 30cm (2g Rope)	5	ELS092	259-5049

# Product Catalogue

## ALGOSTERIL®

Calcium alginate dressing - sterile

### Indications

ALGOSTERIL is ideal for moderately to heavily exuding wounds including:-

- Bleeding wounds
- Infected wounds
- Skin graft donor sites
- Leg ulcers
- Pressure ulcers
- Amputation stumps
- Abscesses/pilonidal sinuses

### Contraindications

Do not use on patients with a known allergy to alginate dressings.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66076937	5cm x 5cm	10	ELS078	239-3692
66076934	10cm x 10cm	10	ELS077	239-3726
66076935	10cm x 20cm	10	ELS079	239-3700



# Product Catalogue

## ALLEVYN<sup>◇</sup> Adhesive/Sacrum

Adhesive polyurethane foam dressing – sterile

### Indications

Wound management by secondary intention on chronic and acute, full thickness or partial thickness or shallow, granulating, moderate to heavily exuding wounds such as pressure ulcers, leg ulcers, infected wounds, malignant wounds, first and second degree burns, donor sites, fungating ulcers, diabetic foot ulcers and surgical wounds.



S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>Adhesive</b>				
66150043	7.5cm x 7.5cm	10	ELA020	220-6720
66000599	10cm x 10cm	10	ELA113	290-9091
66150044	12.5cm x 12.5cm	10	ELA024	220-6738
66000744	12.5cm x 22.5cm	10	ELA046	276-6889
66150045	17.5cm x 17.5cm	10	ELA022	220-6746
66150046	22.5cm x 22.5cm	10	ELA018	220-6753
<b>Sacrum</b>				
66000700	17cm x 17cm (Small)	10	ELS046	273-4606
66000451	22cm x 22cm (Large)	10	ELA099	260-1672

# ALLEVYN<sup>◇</sup> Ag Range

Polyurethane antimicrobial foam dressing – sterile

## Indications

An absorbent antimicrobial dressing for the management by secondary intention of chronic and acute full thickness, partial thickness or shallow, granulating, exuding wounds such as pressure ulcers, venous ulcers, diabetic ulcers, burns, donor sites, fungating/malignant wounds and surgical wounds.



## Contraindications

Do not use on patients known to be hypersensitive to silver sulfadiazine or sulfonamides.

As sulfonamides are known to cause kernicterus, ALLEVYN Ag should not be used on females who are at, or near term pregnancy or lactating, on premature infants or on newborn infants during the first months of life.

S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>Ag Adhesive</b>				
66800073	7.5cm x 7.5cm	10	ELA344	333-1352
66800075	10cm x 10cm	10	ELA345	333-1147
66800078	12.5cm x 12.5cm	10	ELA347	333-1154
66800081	17.5cm x 17.5cm	10	ELA350	333-1162
<b>Ag Non-Adhesive</b>				
66800083	5cm x 5cm	10	ELA343	333-1170
66800086	10cm x 10cm	10	ELA346	333-1188
66800089	15cm x 15cm	10	ELA348	333-1196
66800091	20cm x 20cm	10	ELA351	333-1212
<b>Ag Heel</b>				
66800098	10.5cm x 13.5cm	5	ELA374	341-9397
<b>Ag Sacrum</b>				
66800094	15cm x 15cm (Small)	10	ELA349	333-1295
66800097	20cm x 20cm (Large)	10	ELA352	333-1394



# ALLEVYN<sup>◇</sup> Ag Gentle Range

Polyurethane antimicrobial foam dressing – sterile

## Indications

Wound management by secondary intention on chronic and acute, full thickness or partial thickness or shallow, granulating, exuding wounds such as pressure ulcers, leg ulcers, infected wounds, diabetic foot ulcers, surgical wounds. ALLEVYN Gentle and Gentle Border are suitable for use on fragile skin.

ALLEVYN Gentle and ALLEVYN Gentle Border are suitable for wounds with moderate to high levels of exudate.

ALLEVYN Gentle should be secured with a retention sheet, tape or bandage.

## Contraindications

Do not use on patients known to be hypersensitive to silver sulfadiazine or sulfonamides.

As sulfonamides are known to cause kernicterus, ALLEVYN Ag Gentle and ALLEVYN Ag Gentle Border should not be used on females who are at, or near term pregnancy or lactating, on premature infants or on newborn infants during the first months of life.



S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>Ag Gentle</b>				
66800464	5cm x 5cm	10	ELA393	346-7511
66800465	10cm x 10cm	10	ELA394	346-7529
66800466	10cm x 20cm	10	ELA395	346-7537
66800467	15cm x 15cm	10	ELA396	346-7545
66800468	20cm x 20cm	10	ELA397	346-7552
<b>Ag Gentle Border</b>				
66800460	7.5cm x 7.5cm	10	ELA398	346-7677
66800461	10cm x 10cm	10	ELA399	346-7685
66800462	12.5cm x 12.5cm	10	ELA400	346-7693
66800463	17.5cm x 17.5cm	10	ELA401	346-7701

# Product Catalogue

## ALLEVYN<sup>◇</sup> Cavity

Cavity wound dressing – sterile

### Indications

To manage wounds prior to delayed primary closure. For use on deep pressure sores and leg ulcers, surgical incisions and excisions e.g. pilonidal sinuses.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66007326	5cm Circular	10	ELA141	207-0852
66007327	10cm Circular	5	ELA143	207-0878
66007328	9cm x 2.5cm Tubular	10	ELA145	207-0910
66007329	12cm x 4cm Tubular	5	ELA147	207-0902



# Product Catalogue

## ALLEVYN<sup>◇</sup> Compression

Polyurethane dressing – sterile

### Indications

Indicated for moderate to heavily exuding wounds with sensitive peri-wound skin such as leg ulcers, diabetic foot ulcers, pressure ulcers, superficial and partial thickness burns.

### Contraindications

Do not use on wounds with exposed muscle, bone or tendon.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66047581	5cm x 6cm	10	ELA090	290-9174
66047583	10cm x 10cm	5	ELA091	290-9141
66045125	15cm x 15cm	3	ELA114	290-9158
66047584	15cm x 20cm	3	ELA092	290-9166



# Product Catalogue

## ALLEVYN<sup>◇</sup> Gentle

Soft gel adhesive hydrocellular foam dressing

## Gentle Border/Border Heel

Silicone gel adhesive hydrocellular foam dressing

## Gentle Border Sacrum

Silicone gel adhesive hydrocellular foam dressing

### Indications

Wound management by secondary intention on chronic and acute, full thickness or partial thickness or shallow, granulating, exuding wounds such as pressure ulcers, leg ulcers, infected wounds, diabetic foot ulcers, fungating ulcers, skin tears and surgical wounds.



S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>Gentle</b>				
66800247	5cm x 5cm	10	ELA364	334-2136
66800248	10cm x 10cm	10	ELA360	334-2235
66800249	10cm x 20cm	10	ELA365	334-2243
66800250	15cm x 15cm	10	ELA363	334-2276
66800251	20cm x 20cm	10	ELA357	334-2284
<b>Gentle Border</b>				
66800269	7.5cm x 7.5cm	10	ELA359	334-2052
66800270	10cm x 10cm	10	ELA362	334-2060
66800272	12.5cm x 12.5cm	10	ELA361	334-2102
66800273	17.5cm x 17.5cm	10	ELA358	334-2045
66800900	10cm x 20cm	10	ELA499	364-7765
<b>Gentle Border Heel</b>				
66800506	23cm x 23.2cm	5	ELA392	345-5243
<b>Gentle Border Sacrum</b>				
66800898	16.8cm x 17.1cm	10	ELA498	364-7781



# Product Catalogue

## ALLEVYN<sup>◇</sup> Gentle Border Lite

Silicone gel adhesive hydrocellular foam dressing

### Indications

Wound management by secondary intention on chronic and acute, full thickness or partial thickness or shallow, granulating, exuding wounds such as pressure ulcers, leg ulcers, infected wounds, diabetic foot ulcers, surgical wounds.

ALLEVYN Gentle Border Lite is suitable for use on fragile skin.

ALLEVYN Gentle Border Lite is suitable for wounds with low to moderate levels of exudate.

ALLEVYN Gentle Border Lite stays in place without the need for secondary retention.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66007630	5cm x 5cm	10	ELA467	359-1666
66800834	7.5cm x 7.5cm	10	ELA468	359-1674
66800835	10cm x 10cm	10	ELA469	359-1682
66800836	5.5cm x 12cm	10	ELA470	359-1690
66800837	8cm x 15cm	10	ELA471	359-1708
66800840	15cm x 15cm	10	ELA472	359-1716



# Product Catalogue

## ALLEVYN<sup>◇</sup> Heel/Ag Heel

Non-adhesive polyurethane dressing – sterile

### Indications

Wound healing by secondary intention for chronic, acute, full thickness, partial thickness, shallow, granulating and exuding wounds such as pressure ulcers, leg ulcers, infected wounds, diabetic foot ulcers and surgical wounds. ALLEVYN Heel and ALLEVYN Ag Heel are anatomically shaped dressings for use on heel wounds. ALLEVYN Ag Heel is indicated for wounds showing the signs and symptoms of infection. Where the product is used on infected wounds the wounds should be treated as per local protocol.

### Contraindications

ALLEVYN Ag Heel should not be used on patients known to be hypersensitive to silver sulfadiazine or sulfonamides.

As sulfonamides are known to cause kernicterus, ALLEVYN Ag should not be used on females who are at, or near term pregnancy or lactating, on premature infants or on newborn infants during the first months of life.



S&N Code	Description	Pack Content	NHS Code	PIP Code
ALLEVYN Heel 66007630	10.5cm x 13.5cm	5	ELA100	260-1714
ALLEVYN Ag Heel 66800098	10.5cm x 13.5cm	5	ELA374	341-9397



# Product Catalogue

---

## ALLEVYN<sup>◇</sup> Lite

Non-adhesive polyurethane dressing – sterile

### Indications

Indicated for low to moderately exuding wounds, healing by secondary intention.

ALLEVYN Lite should be secured with a retention sheet, tape or bandage.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66000690	5cm x 5cm	20	ELA030	272-8715
66000691	10cm x 10cm	20	ELA032	272-8723
66000692	10cm x 20cm	10	ELA033	272-8731
66000698	15cm x 20cm	10	ELA034	272-8749



# Product Catalogue

## ALLEVYN<sup>◇</sup> Non-Adhesive Polyurethane dressing – sterile

### Indications

Wound management by secondary intention on chronic and acute, full thickness or partial thickness or shallow, granulating, moderate to heavily exuding wounds such as pressure ulcers, leg ulcers, infected wounds, diabetic foot ulcers, fungating ulcers, donor sites and surgical wounds.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66157643	5cm x 5cm	10	ELA129	207-0787
66157637	10cm x 10cm	10	ELA131	207-0811
66157335	10cm x 20cm	10	ELA101	207-0829
66157638	20cm x 20cm	10	ELA133	207-0837



# ALLEVYN<sup>◇</sup> Plus Adhesive/ Sacrum

Adhesive polyurethane dressing – sterile

## Indications

Wound management by secondary intention on chronic and acute, full thickness or partial thickness or shallow granulating, heavy to extra heavy exuding wounds such as pressure ulcers, leg ulcers, infected wounds, diabetic foot ulcers, surgical wounds, malignant wounds, first and second degree burns, donor sites and fungating ulcers.



S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>Plus Adhesive</b>				
66000805	12.5cm x 12.5cm	10	ELA119	290-9182
66000807	12.5cm x 22.5cm	10	ELA122	290-9190
66000806	17.5cm x 17.5cm	10	ELA120	290-9208
<b>Plus Adhesive Sacrum</b>				
66000856	17cm x 17cm	10	ELA527	314-6974
66000856	22cm x 22cm	10	ELA258	314-6958

# Product Catalogue

## ALLEVYN<sup>◇</sup> Plus Cavity

Cavity wound dressing – sterile

### Indications

Management of deep secondary healing wounds, particularly stage III and IV, deep leg ulcers, cavity wounds, excisions and post-op wound dehiscence.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66047571	5cm x 6cm	10	ELA093	290-9240
66047573	10cm x 10cm	5	ELA094	290-9224
66047574	15cm x 20cm	3	ELA095	290-9232



# Product Catalogue

---

## ALLEVYN<sup>◇</sup> Thin

Polyurethane foam dressing – sterile

### Indications

Indicated for low to moderately exuding granulating wounds, such as pressure ulcers, post-operative wounds, donor sites, second degree burns, leg ulcers.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66047576	5cm x 6cm	10	ELA116	290-9133
66047578	10cm x 10cm	5	ELA117	290-9109
66045124	15cm x 15cm	3	ELA128	290-9117
66047579	15cm x 20cm	3	ELA118	290-9125



# Product Catalogue

---

## ALLEVYN Tracheostomy

Hydrophilic polyurethane (keyhole) dressing – sterile

### Indications

Absorbent apertured dressing for the stoma created by tracheostomy or for wound drainage sites.

ALLEVYN Tracheostomy should be secured with a retention sheet, tape or bandage.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66007640	9cm x 9cm Keyhole	10	ELA121	-



# Product Catalogue

## AMETOP<sup>◇</sup> Gel

4.0% w/w Tetracaine

### Indications:

Percutaneous local anaesthetic to produce anaesthesia of the skin prior to venepuncture or venous cannulation.

For prescribing information see page 67 of booklet.

Legal Category: P



S&N Code	Description	Pack Content	NHS Code	PIP Code
66000219	1.5g tube	12	-	224-3475
66000311	12 x 1.5g tubes with 24 occlusive dressings (Dispensing pack)	12	-	237-5350

# Product Catalogue


## BACTIGRAS<sup>◇</sup>

### Medicated tulle gras - sterile

#### Indications

Bactigras is a wound contact layer for wounds such as:

- Minor burns and scalds
- Lacerations, abrasions and other skin loss wounds
- Donor and recipient graft sites
- Where the wound is no more than 10% of the body area, where there is a risk of infection, or on already infected wounds in conjunction with systemic antibiotics.



#### Contraindications

It should not be used in known cases of sensitivity or allergy to chlorhexidine

S&N Code	Description	Pack Content	NHS Code	PIP Code
7456	5cm x 5cm	50	EKB011	013-7943
7457	10cm x 10cm	10	EKB503	017-9432
7461	15cm x 20cm	10	EKB015	000-9779
66007505	15cm x 1m roll	12	EKB002	226-8258

## CADESORB<sup>◇</sup>

### Protease Modulating ointment

#### Indications

Improving/managing the wound environment to promote healing in chronic wounds (healing by secondary intention).

#### Contraindications

Do not use on dry necrotic tissue (where no exudate is present) or in known cases of sensitivity to CADESORB Ointment or any of its ingredients.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66000056	10g tube	5	ELZ085	320-7040
66000055	20g tube	2	ELZ068	313-5472



# Product Catalogue

---

## CARBONET<sup>◇</sup>

### Odour Absorbing Charcoal Dressing

**Indications:**

To dress discharging, infected and malodorous wounds such as venous leg ulcers, fungating carcinomas, faecal fistulae and infected surgical wounds.



S&N Code	Description	Pack Content	NHS Code	PIP Code
7064	10cm x 10cm	10	ELV211	001-2807
7065	10cm x 20cm	10	ELV213	001-2997

# Product Catalogue

## CICA-CARE<sup>◇</sup>


### Silicone gel sheet

#### Indications

Treatment for red and raised, old and new scars.

CICA-CARE is designed for use:

- In the management of both existing and new red, raised scars (hypertrophic/keloid).
- As a preventative therapy on healed wounds to help prevent red and raised scars.



#### Contraindications

CICA-CARE is contra-indicated for patients with complicating medical factors which would make them unable to use the dressing properly and in patients with dermatological conditions which disrupt the integrity of the skin in areas of coverage.

Do not use on open or infected wounds or skin currently affected by acne.

S&N Code	Description	Pack Content	NHS Code	PIP Code
66000702	6cm x 12cm	3	ELY061	273-0752
66250706	15cm x 12cm	1	ELY055	273-0745
66250707	15cm x 12cm	10	-	-



# Product Catalogue

---

## CICA-PLAIE<sup>◇</sup>

Low allergy adhesive wound dressing

### Indications

To dress post-operative wounds; or cuts, lacerations or sutured wounds.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66660272	5cm x 7cm	100	-	-
66660273	10cm x 8cm	50	-	-
66660274	15cm x 10cm	50	-	-
66660275	20cm x 10cm	50	-	-
66660276	25cm x 10cm	25	-	-
66660277	30cm x 10cm	25	-	-
66660278	35cm x 10cm	25	-	-

# Product Catalogue

## CUTINOVA<sup>◇</sup> Hydro

### Self-adhesive hydro-selective dressing

#### Indications

For the treatment and management of exuding wounds that heal by secondary intention such as leg ulcers, pressure ulcers, superficial and partial thickness burns and diabetic foot ulcers and as a secondary dressing on cavity wounds where the surrounding skin is intact and in a normal condition.

#### Contraindications

Do not use on wounds with exposed muscle, bone or tendon where a risk of drying out occurs. Do not use on wounds where the surrounding skin is damaged or very sensitive.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66047441	5cm x 6cm	10	ELW751	099-4657
66047443	10cm x 10cm	10	ELW753	099-4988
66047445	15cm x 20cm	3	ELW755	099-5084

## CUTIPLAST<sup>◇</sup> Steril

### Absorbent dressing with adhesive border - sterile

#### Indications

- Post-operative wounds
- Cuts and lacerations
- Sutured wounds



S&N Code	Description	Pack Content	NHS Code	PIP Code
66001478	7.2cm x 5cm	100	-	036-0297
66001473	10cm x 8cm	50	-	036-0305
66001474	15cm x 8cm	50	-	036-0313
66001475	20cm x 10cm	50	-	036-0321
66001476	25cm x 10cm	50	-	291-3382
66001477	30cm x 10cm	50	-	291-3390



# Product Catalogue

## DURAFIBER<sup>◇</sup>

Highly absorbent gelling fibre dressing

### Indications

For the management of chronic and acute, full thickness, partial thickness, or shallow granulating exuding wounds. For example: leg ulcers; pressure ulcers; diabetic ulcers; surgical wounds; wounds left to heal by secondary intent; donor sites; tunnelling and fistulae wounds; partial thickness burns; traumatic wounds; and wounds that are prone to bleeding such as wounds that have been surgically or mechanically debrided.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66800559	5cm x 5cm	10	ELY330	356-4275
66800560	10cm x 10cm	10	ELY331	356-4283
66800561	15cm x 15cm	5	ELY332	356-4291
66800546	4cm x 10cm	5	ELY334	356-4317
66800547	4cm x 20cm	5	ELY335	356-4325
66800548	4cm x 30cm	5	ELY337	356-4333
66800563	2cm x 45cm	5	ELY333	356-4309



# Product Catalogue

## EXU-DRY<sup>◇</sup>

### One-piece anti-shear wound dressings

#### Indications

Under the care of a healthcare professional, EXU-DRY dressings may be used in the management of superficial to full thickness wounds such as:

- Ulcers (venous, arterial, diabetic)
- Pressure sores
- Donor Sites
- Surgical incisions and excisions
- Skin grafts
- Bio-engineered skin substitutes
- Burns
- Draining wounds
- Moist skin desquamation
- Fungating neoplasms
- Chemotherapy extravasation
- Kaposi's sarcoma
- Lymphoedema
- Graft Vs Host Disease
- Dermatological wounds
- Skin sloughing disorders



S&N Code	Description	Pack Content	NHS Code	PIP Code
5999004120	10cm x 15cm Dressing Full absorbency	10	ELS044	294-0872
5999006	15cm x 23cm Dressing Full absorbency	12	ELS045	294-0864
5999009	23cm x 38cm Dressing Full absorbency	30	ELS012	294-0880
5999018	38cm x 46cm Full absorbency	30	ELS013	-
5999024	38cm x 61cm Dressing Medium absorbency	30	ELS014	-
5999M36	61cm x 91cm Pad Non-permeable Heavy absorbency	15	ELS032	-
5999L72	91cm x 183cm Sheet Non-permeable	15	ELS020	-
5999LPA	Arm dressing	20	ELS021	-
5999LPL	Leg plus dressing Heavy absorbency	20	ELS029	-
5999LV1	Large burn vest	20	ELS030	-



# Product Catalogue

## FLAMAZINE<sup>◇</sup> Cream

1% w/w silver sulfadiazine

### Indications

Short-term treatment of infection in leg ulcers and pressure ulcers.

Conservative management of fingertip injuries.

Prophylaxis of infection in skin graft donor sites.

Prophylaxis of infection in extensive abrasions.

Prophylaxis and treatment of infection in burns wounds.

For prescribing information see page 68 of booklet.

Legal category POM



S&N Code	Description	Pack Content	NHS Code	PIP Code
66150257	20g tube POM	Single	-	267-7375
66150153	50g tube POM	Single	-	049-5432
66150154	250g jar POM	Single	-	001-7517
66150155	500g jar POM	Single	-	001-7525

# Product Catalogue

## ICHTHOPASTE<sup>◇</sup>

Zinc oxide paste and ichthammol bandage

### Indications

To assist the management of venous leg ulcers. Where venous insufficiency exists the paste bandage should be used under graduated compression bandaging. For the management of chronic eczema/dermatitis where occlusion is indicated.

### Contraindications

Do not use in known cases of sensitivity or allergy to any of the ingredients.



S&N Code	Description	Pack Content	NHS Code	PIP Code
4959	Bandage 7.5cm x 6m	12	EFA051	033-2668



# Product Catalogue

## INTRASITE<sup>◇</sup> Conformable

Non-woven material impregnated with INTRASITE Gel - sterile

### Indications

For use in shallow and deep open wounds healing by secondary intention.

Especially useful for the packing of deep cavity wounds to prevent premature closure.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66000324	10cm x 10cm (7.5g Gel loading)	10	ELG002	267-8043
66000325	10cm x 20cm (15g Gel loading)	10	ELG007	267-8084
66000326	10cm x 40cm (30g Gel loading)	10	ELG008	267-8100

## INTRASITE Gel

Amorphous hydrogel dressing - sterile

### Indications

For use in shallow and deep open wounds healing by secondary intention.



S&N Code	Description	Pack Content	NHS code	PIP Code
7308	APPLIPAK 8g	10	ELG015	210-2929
7311	APPLIPAK 15g	10	ELG000	004-2895
7313	APPLIPAK 25g	10	ELG011	006-3859

# Product Catalogue

## IODOFLEX<sup>◇</sup>

Cadexomer dressing with iodine

## IODOSORB<sup>◇</sup>

Cadexomer ointment/powder with iodine


### Indications

For the treatment of chronic exuding wounds.

IODOSORB is available as a dressing (IODOFLEX), ointment and powder which perform the same functions.

### Contraindications

Do not use on dry necrotic tissue or on patients with a known sensitivity to any of its ingredients. Do not use on children, pregnant or lactating women or people with thyroid disorders or renal impairment.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66001301	IODOFLEX 5g dressing	5	EKB007	073-1547
66001302	IODOFLEX 10g dressing	3	EKB008	014-9617
66001303	IODOFLEX 17g dressing	2	EKB009	216-3269
66001286	IODOSORB 3g Powder	7	EKB010	045-3696
66001298	IODOSORB 10g Ointment	4	EKB012	036-6658
66001297	IODOSORB 20g Ointment	2	EKB018	033-3906



# Product Catalogue

## IV3000<sup>◇</sup>

Vapour-permeable adhesive film IV dressing - sterile

### Indications

For dressing central line, peripheral, epidural and paediatric catheters.

### Contraindications

The dressings should be used on IV or epidural sites. If a film dressing for open wounds, such as pressure sores, burns and donor sites is required, OPSITE FLEXIGRID should be used.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66004011*	5cm x 6cm Ported paediatric catheter	100	ELW011	-
4006*	7cm x 9cm Ported peripheral catheters	100	ELW046	276-0783
66004009*	9cm x 12cm central/jugular catheters	50	ELW112	320-2827
4007*	6cm x 7cm peripheral catheters	100	ELW032	238-0269
66004005*	10cm x 12cm Central line catheters	10	-	305-7411
4008*	10cm x 12cm Central line catheters	50	ELW054	238-0285

\* 2 sterile securing strips and date label included

### Frame Delivery Application

59410082*	6cm x 7cm Peripheral catheters	100	ELW105	-
59410882*	10cm x 12cm Central catheters	50	ELW054	-

### IV3000 Orange Handle Application

4923	6cm x 8cm Ported peripheral catheters	100	ELW037	238-0251
4924	6cm x 8.5cm peripheral catheters	100	ELW039	320-2819
4925	10cm x 14cm Central line catheters	10	ELW041	238-0277
4649	10cm x 20cm Epidural catheters	50	ELW099	-

# Product Catalogue

## JELONET<sup>◇</sup>

Non-medicated tulle gras - sterile

### Indications

Minor burns and scalds, donor and recipient graft sites, skin loss wounds, lacerations, abrasions, leg ulcers.

### Contraindications

Do not use in cavity wounds.



S&N Code	Description	Pack Content	NHS Code	PIP Code
7403	5cm x 5cm BP	50	EKA009	217-0439
7404	10cm x 10cm BP	10	EKA501	027-8952
7409	10cm x 10cm BP	100	EKA011	217-0447
7459	10cm x 40cm BP	10	EKA013	004-8595
66007477	10cm x 7m BP tin	Single	EKA500	004-8710
7415	15cm x 2m roll	12	EKA014	217-0454



# Product Catalogue

## LEUKOSTRIP<sup>◇</sup>

Skin closure strips - sterile

### Indications

For primary and secondary closure of skin wounds. Used in the fixation of skin grafts.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66002903	4mm x 76mm x 4 strips	50 x 4	EIR510	046-1913
66002876	4mm x 38mm x 8 strips	50 x 8	EIR002	046-1731
66002878	6.4mm x 76mm x 3 strips	50 x 3	EIR512	046-2044
66002879	6.4mm x 102mm x 5 strips	50 x 5	EIR508	046-2051
66002880	13mm x 102mm x 6 strips	50 x 6	EIR509	046-2069
66002881	26mm x 102mm x 4 strips	25 x 4	EIR012	046-2135
66002952	6.4mm x 76mm x 3 strips	10 x 3	-	046-2317
66002953	13mm x 102mm x 6 strips	10 x 6	-	046-2325



# Product Catalogue

## MELOLIN<sup>◇</sup>

Low-adherent absorbent dressing - sterile

### Indications

A low-adherent absorbent dressing used for the management of a wide variety of light to moderately exuding wounds including clean sutured wounds, abrasions, lacerations and minor burns.



S&N Code	Description	Pack Contents	NHS Code	PIP Code
66974940	5cm x 5cm Sterile	100	EJE011	-
66974941	10cm x 10cm Sterile	100	EJE013	-
66974939	20cm x 10cm Sterile	100	EJE502	-
66974942	10cm x 10cm Bulk Non-Sterile	150	EJE033	-
66974943	10cm x 20cm Bulk Non-Sterile	75	EJE035	-
66974944	20cm x 30cm Bulk Non-Sterile	25	EJE037	-
66974930	50cm x 7m Non-Sterile	Single	EJE039	226-8241



# Product Catalogue

---

## MELOLITE<sup>◇</sup>

Low-adherent absorbent dressing - sterile

### Indications

For casualty and post-operative use over clean, sutured wounds and for lacerations and abrasions.



S&N Code	Description	Pack Contents	NHS Code	PIP Code
4811	5 x 7.5cm	100	EJE171	-
4812	10cm x 7.5cm	100	EJE173	-
4814	20cm x 2.5cm	100	EJE177	-



# Product Catalogue

## NO-STING SKIN-PREP<sup>◇</sup>

### Protective dressing

#### Indications

NO-STING SKIN-PREP is a sterile liquid film-forming product, containing a non-stinging solvent, which dries rapidly. NO-STING SKIN-PREP helps to protect intact or damaged skin from:

- Urinary incontinence
- Faecal incontinence
- Digestive secretions
- Fluids draining from wounds and drains
- Adhesives
- Frictional forces

#### Contraindications

Do not use on patients with a known allergy or sensitivity to any of the ingredients. This product is not recommended for use on sensitive areas, such as around the face and eyes unless under the guidance of a healthcare professional. Do not apply to infected areas.



S&N Code	Description	Pack Content
66800712	NO-STING SKIN-PREP wipes	50 wipes
66800789	NO-STING SKIN-PREP 1ml swabs	5 swabs
66800790	NO-STING SKIN-PREP 1ml swabs	25 swabs
66800787	NO-STING SKIN-PREP 3ml swabs	5 swabs
66800788	NO-STING SKIN-PREP 3ml swabs	25 swabs
66800711	NO-STING SKIN-PREP 3ml swabs	50 swabs
66800710	NO-STING SKIN-PREP spray	1 28ml/1oz bottle



# Product Catalogue

---

## OPSITE<sup>◇</sup> FLEXIFIX

Transparent Film Roll (non-sterile)

### Indications

To relieve pain in painful diabetic neuropathy, to retain primary dressings, Prophylaxis of pressure ulcers and to protect the skin around stoma and under leg bags.

### Contraindications

Do not apply to open wounds.



S&N Code	Description	Pack Contents	NHS Code	PIP Code
66000040	5cm x 10cm	Single	ELW101	276-0593
66000041	10cm x 10cm	Single	ELW102	276-0643

# Product Catalogue

## OPSITE<sup>◇</sup> FLEXIGRID<sup>◇</sup>

Moisture vapour permeable adhesive film dressing - sterile

### Indications

- Superficial pressure ulcers and their prophylaxis
- Abrasions and Lacerations
- Skin grafts and donor sites
- Post-operative wounds
- Minor burns and scalds
- Skin protection around stoma and leg bags
- Relief of pain in diabetic neuropathy



S&N Code	Description	Pack Content	NHS Code	PIP Code
4628	6cm x 7cm	100	ELW322	043-6170
4629	10cm x 12cm	10	-	-
4630	10cm x 12cm	50	ELW325	043-6386
66004964	12cm x 12cm	10	ELW118	034-4283
4632	12cm x 25cm	20	ELW329	043-6519
4631	15cm x 20cm	10	ELW327	043-6493



# Product Catalogue

## OPSITE<sup>◇</sup> Incise

### Adhesive incise drape

#### Indications

For use as an incise drape in major and minor surgery including:

- Orthopaedic
- Neurosurgery
- Major abdominal
- Plastic surgery
- Open heart and thoracic
- Paediatric



S&N Code	Description	Pack Content	NHS Code	PIP Code
4963	10cm x 14cm	10	-	217-0462
4975	10cm x 14cm	50	-	217-0470
4986	15cm x 28cm	10	-	034-5975
4987	30cm x 28cm	10	-	034-5983
4988	45cm x 28cm	10	-	004-8884
4995	42cm x 40cm	10	-	217-0488
4989	45cm x 55cm	10	-	217-0496
4994	56cm x 84cm	10	-	217-0504

# Product Catalogue

## OPSITE<sup>◇</sup> Post-Op / Plus

Adhesive film dressing with absorbent pad - sterile

### Indications

Low to moderately exuding acute wounds such as: Post-operative wounds, minor cuts, abrasions and lacerations.



S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>OPSITE Post-Op - for hospital use</b>				
66000708	6.5cm x 5cm	100	ELW052	-
66000709	9.5cm x 8.5cm	20	ELW051	-
66000710	12cm x 10cm	10	ELW090	-
66000712	15.5cm x 8.5cm	20	ELW050	-
66000713	20cm x 10cm	20	ELW091	-
66000714	25cm x 10cm	20	ELW092	-
66000715	30cm x 10cm	20	ELW045	-
66000716	35cm x 10cm	20	ELW049	-
<b>OPSITE Plus - for community use</b>				
66000760	5cm x 5cm	25	-	280-4904
66000761	8.5cm x 9.5cm	20	-	280-4912
66000763	10cm x 20cm	20	-	280-4839
66000762	12cm x 10cm	10	-	280-4821
66000764	35cm x 10cm	20	-	280-4896



# Product Catalogue

## OPSITE<sup>◇</sup> Post-Op Visible

Waterproof, bacteria-proof dressing with see-through absorbent pad

### Indications

Specifically designed for low to moderately exuding acute wounds, lacerations, minor cuts, and minor burns.



S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>OPSITE Post-Op Visible</b>				
66800136	10cm x 8cm	20	ELW307	-
66800137	15cm x 10cm	20	ELW227	-
66800138	20cm x 10cm	20	ELW228	-
66800139	25cm x 10cm	20	ELW229	-
66800140	30cm x 10cm	20	ELW307	-
66800141	35cm x 10cm	20	ELW308	-
<b>OPSITE Visible Drain Dressing</b>				
66800842	10cm x 9cm	20	ELW336	-



# Product Catalogue

## OPSITE<sup>◇</sup> Spray

Transparent film spray dressing

### Indications

For use on a variety of minor dry wounds, in particular:

- Minor cuts, abrasions and burns
- Over dry sutured wounds
- Sealing of external orthopaedic fixation, such as Steinman pins
- After suture removal
- Over vaccination sites
- Over intact blisters
- Skin graft fixation

### Contraindications

Do not use with plastic intravascular devices. Do not apply on children under the age of 18 months.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66004978	Protective spray Dressing 100ml	12	ELY031	217-0645
66004980	Protective spray Dressing 240ml	12	ELY034	217-0660



# Product Catalogue

## PICO<sup>◇</sup>

### Negative Pressure Wound Therapy

#### Indications

PICO is indicated for a range of wound types including chronic wounds, incision sites (preventative use of NPWT on closed wounds), acute wounds, traumatic wounds, sub-acute and dehisced wounds, partial-thickness burns, ulcers (such as diabetic or pressure), flaps and grafts, dermal substitutes.

#### Contraindications

The use of PICO is contraindicated in the presence of necrotic tissue with eschar, untreated osteomyelitis, malignancy in wound (except palliative care to enhance quality of life), exposed arteries, veins or organs, non-enteric and unexplored fistulae and anastomotic sites.

Not to be used on children under the age of 18 months.

#### Availability

PICO can be purchased as a single unit, with a variety of dressing sizes.



S&N Code Kits	Description	Pack Content	NHS Code	PIP Code
66800862*	10cm x 20cm	Unit	ELZ348	-
66800863	10cm x 30cm	Unit	ELZ349	366-0032
66800864*	10cm x 40cm	Unit	-	-
66800865	15cm x 15cm	Unit	ELZ350	366-0024
66800866	15cm x 20cm	Unit	ELZ351	366-0016
66800867	15cm x 30cm	Unit	-	-
66800868	20cm x 20cm	Unit	ELZ352	-
66800869	25cm x 25cm	Unit	-	-

\*Available in November 2011

# Product Catalogue

---


## PRIMAFIX<sup>◇</sup>

Hypoallergenic non-woven retention sheet

### Indications

Secure retention of primary dressings, particularly on joints and contoured areas of the body.

Reliable fixation of catheters and draining tubes.



S&N Code	Description	Pack Content
66800151	5cm x 10m	12 Rolls
66800152	10cm x 10m	12 Rolls
66800153	15cm x 10m	12 Rolls
66800154	20cm x 10m	12 Rolls



# Product Catalogue

---

## PRIMAPORE<sup>◇</sup>

Absorbent dressing with adhesive border - sterile

### Indications

To dress post-operative wounds or cuts, lacerations or sutured wounds.



S&N Code	Description	Pack Content	NHS Code	PIP Code
7135	8.3cm x 6cm	50	EUJ11	015-4088
66000317	10cm x 8cm	20	EUJ050	261-8932
66000318	15cm x 8cm	20	EUJ051	261-8965
66000319	20cm x 10cm	20	EUJ052	261-8973
66000320	25cm x 10cm	20	EUJ053	261-8981
66000321	30cm x 10cm	20	EUJ054	261-8999
66007140	35cm x 10cm	20	EUJ031	316-8192

# Product Catalogue

## PROFORE<sup>◇</sup>

### Multi-layer compression bandaging system

#### Indications

PROFORE multi-layer compression bandage system is suitable for venous leg ulcers, with a Doppler ABPI greater than 0.8.

PROFORE Lite is suitable for mixed aetiology leg ulcers (which require reduced compression), with a Doppler ABPI between 0.6 and 0.8.



S&N Code Kits	Description	Pack Content	NHS Code	PIP Code
66000015	Less than 18cm	4	EVN022	226-8076
66000016	18-25cm	8	EVN004	200-8506
66000017	25-30cm	8	EVN015	226-8134
66000018	Greater than 30cm	8	EVN018	226-8167
66000415	PROFORE Lite	8	EVN039	244-9338
<b>Components</b>				
66000701	PROFORE WCL	50	ELA127	274-9802
67000333	PROFORE #1	12	EBA053	246-1663
67000328	PROFORE #2	12	ECA029	258-5750
67000334	PROFORE #3	12	EBA050	258-5768
67000336	PROFORE #4	12	ECD007	246-1671
67000335	PROFORE +	12	EBA023	246-1689
<b>Kits Latex Free</b>				
66000770	18-25cm	8	ECA037	320-2793
66600771	PROFORE Lite	8	ECA055	320-2801
<b>Components Latex-free</b>				
66000701	PROFORE WCL	50	ELA127	274-9802
66000757	PROFORE #1	12	ECA131	322-2700
66000758	PROFORE #2	12	ECA132	322-2718
66000778	PROFORE #3	12	ECA133	322-2726
66000779	PROFORE #4	12	ECA134	322-2734
66000776	PROFORE +	12	ECA135	322-2742



# Product Catalogue

## PROGUIDE<sup>◇</sup>

VARI-STRETCH<sup>◇</sup> two layer bandaging system

### Indications

PROGUIDE is suitable for the management of venous leg ulcers and associated conditions.



S&N Code	Description	Pack Content	NHS Code	PIP Code
<b>Kits</b>				
66000780	18-22cm (Red kit)	8	ECA031	295-3990
66000781	22-28cm (Yellow kit)	8	ECA033	295-4006
66000782	28-32cm (Green kit)	8	ECA034	295-4014
<b>Components</b>				
66000783	WCL	10	EVN109	313-8047
66000784	PROGUIDE #1	12	EVN110	313-8054
66000798	PROGUIDE #2 (Red)	12	EVN111	313-8138
66000799	PROGUIDE #2 (Yellow)	12	EVN112	313-8112
66000800	PROGUIDE #2 (Green)	12	EVN113	313-8088

# Product Catalogue

## RENASYS<sup>◇</sup> EZ PLUS

### Negative Pressure Wound Therapy

#### Indications

RENASYS EZ PLUS is indicated for patients who would benefit from a vacuum device, particularly as the device may promote wound healing. RENASYS EZ PLUS is appropriate for post-operative via the removal of fluid and dehisced wounds, pressure ulcers, diabetic/neuropathic ulcers, traumatic wounds, skin flaps and grafts, venous insufficiency ulcers and explored fistulae.

#### Contraindications

The use of RENASYS EZ PLUS is contraindicated in the presence of necrotic tissue with eschar, untreated osteomyelitis, malignancy in wound (except palliative care to enhance quality of life), exposed arteries, veins or organs, non-enteric and unexplored fistulae and anastomotic sites.

#### Availability

RENASYS EZ PLUS devices are available as and when required through a range of flexible packages to suit your needs. Packages include ad hoc rental, monthly rental, rent to buy and outright purchase.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66800697	RENASYS EZ PLUS Device	unit	-	-
<b>Consumables</b>				
66800391	Antimicrobial Large Gauze Rolls (pack of 5)	Single	ELZ251	-
66800394	Large Drape 20cm x 30cm (pack of 10)	Single	ELZ249	-
<b>Dressing Kits</b>				
66800491	RENASYS-G Small Round Kit (10Fr)	Single	ELZ238	346-7941
66800492	RENASYS-G Small Flat Kit	Single	ELZ239	346-7958
66800493	RENASYS-G Medium Channel Kit	Single	ELZ242	346-7966
66800494	RENASYS-G Medium Flat Kit	Single	ELZ243	346-7974
66800495	RENASYS-G Large Round Kit (19Fr)	Single	ELZ241	346-7933
66800496	RENASYS-G Large Flat Kit	Single	ELZ240	346-7925
66800932	RENASYS-G Large Irrigation/Aspiration Kit	Single	ELZ170	-
66800133	RENASYS-G Sterile Kit	Single	ELZ236	-
66800639	RENASYS-F/P Small Kit	Single	ELZ274	355-8954
66800640	RENASYS-F/P Medium Kit	Single	ELZ275	355-8962
66800641	RENASYS-F/P Large Kit	Single	ELZ276	355-8970
66800694	RENASYS Port Kit	Single	ELZ282	357-5339
66800713	RENASYS-F/AB Abdominal Dressing Kit	Single	ELZ347	-
66800882	RENASYS-G/P Small Kit	Single	ELZ327	-
66800883	RENASYS-G/P Medium Kit	Single	ELZ328	-
66800884	RENASYS-G/P Large Kit	Single	ELZ329	-
<b>Canisters</b>				
66800946	800ml Sealed Canister Kit	Single	ELZ165	343-9882
<b>Accessories</b>				
66800504	RENASYS Y Connector	Single	ELZ250	-



# Product Catalogue

## RENASYS<sup>◇</sup> GO

### Negative Pressure Wound Therapy

#### Indications

RENASYS GO is indicated for a range of wound types including chronic, acute, traumatic, sub-acute, dehisced, partial thickness burns, flaps, grafts, pressure ulcers and diabetic ulcers.

#### Contraindications

The use of RENASYS GO is contraindicated in the presence of necrotic tissue with eschar, untreated osteomyelitis, malignancy in wound (except palliative care to enhance quality of life), exposed arteries, veins or organs, non-enteric and unexplored fistulae and anastomotic sites.

#### Availability

RENASYS GO devices are available as and when required through a range of flexible packages to suit your needs. Packages include ad hoc rental, monthly rental, rent to buy and outright purchase.



S&N Code	Description	Pack Content Unit	NHS Code	PIP Code
66800164	RENASYS GO Device	Unit	ELZ246	-
<b>Consumables</b>				
66800391	Antimicrobial Large Gauze Rolls (pack of 5)	Single	ELZ251	-
66800394	Large Drape 20cm x 30cm (pack of 10)	Single	ELZ249	-
<b>Dressing Kits</b>				
66800491	RENASYS-G Small Round Kit (10Fr)	Single	ELZ238	346-7941
66800492	RENASYS-G Small Flat Kit	Single	ELZ239	346-7958
66800493	RENASYS-G Medium Channel Kit	Single	ELZ242	346-7966
66800494	RENASYS-G Medium Flat Kit	Single	ELZ243	346-7974
66800495	RENASYS-G Large Round Kit (19Fr)	Single	ELZ241	346-7933
66800496	RENASYS-G Large Flat Kit	Single	ELZ240	346-7925
66800932	RENASYS-G Large Irrigation/Aspiration Kit	Single	ELZ170	-
66800133	RENASYS-G Sterile Kit	Single	ELZ236	-
66800639	RENASYS-F/P Small Kit	Single	ELZ274	355-8954
66800640	RENASYS-F/P Medium Kit	Single	ELZ275	355-8962
66800641	RENASYS-F/P Large Kit	Single	ELZ276	355-8970
66800694	RENASYS Port Kit	Single	ELZ282	357-5339
66800713	RENASYS-F/AB Abdominal Dressing Kit	Single	ELZ347	-
66800882	RENASYS-G/P Small Kit	Single	ELZ327	-
66800883	RENASYS-G/P Medium Kit	Single	ELZ328	-
66800884	RENASYS-G/P Large Kit	Single	ELZ329	-
<b>Canisters</b>				
66800941	RENASYS GO 300ml Canister	Single	ELZ237	347-3782
66800943	RENASYS GO Large Canister	Single	ELZ279	358-2657
<b>Accessories</b>				
66800162	RENASYS GO Carrying Case	Single	ELZ244	-
66800163	RENASYS GO Carrying Case Strap	Single	ELZ245	-
66800504	RENASYS Y Connector	1	ELZ250	-



# Product Catalogue


---

## TRICOTEX<sup>◇</sup>

Low adherent wound contact dressing - sterile

### Indications

Intended for use as an interface layer between ulcerating or granulating wounds and conventional absorptive dressings.



S&N Code	Description	Pack Content	NHS Code	PIP Code
7111	9.5cm x 9.5cm	50	EKG11	028-0438



# Product Catalogue

## TRIPLE CARE<sup>◇</sup>

### Skin protection creams

#### Indications

Moisturises dry skin and is a protectant that aids in reducing exposure to faeces or urine. Triple Care Cream contains Clove Oil which helps to control odour.



S&N Code	Description	Pack Content	NHS Code	PIP Code
59447525	Cream 78g	Single	-	208-6288
59446825	Cleanser 237ml	Single	-	208-6270
59450425	Extra protection Barrier cream 92g	Single	-	208-6304

# Product Catalogue

## VERSAJET<sup>◇</sup>

### Hydrosurgery system

The VERSAJET Hydrosurgery System is intended to cut, ablate and remove tissue from wounds, and to resect and remove material in various surgical procedures including wound debridement.

#### Indications

- Pressure ulcers
- Diabetic Foot ulcers
- Burns and plastic surgery
- Leg ulcer
- Dehisced wounds
- Non-healing wounds



S&N Code	Description	Pack Content	NHS Code	PIP Code
50750	Power Console	Single	-	-
50635	Disposable Handpiece 15°/14mm	Single	-	-
50636	Disposable Handpiece 45°/14mm	Single	-	-
50637	Disposable Handpiece 45°/8mm	Single	-	-

Smith & Nephew recommend that only VERSAJET Handpiece codes 50636 and 50637 are used out of theatre.



# VISCOPASTE<sup>◇</sup> PB7

## Zinc oxide paste bandage

### Indications

To assist the management of venous leg ulcers. Where venous insufficiency exists the paste bandage should be used under graduated compression bandaging. For the management of chronic eczema/dermatitis where occlusion is indicated.

### Contraindications

Do not use in known cases of sensitivity or allergy to any of the ingredients.



S&N Code	Description	Pack Content	NHS Code	PIP Code
4948	Zinc Oxide 7.5cm x 6m	12	EFA011	033-2734

# Product Catalogue

---

## VISITRAK<sup>◇</sup> Depth Wound Measurement System

### Indications

For measurement of wound depth; can be used on all wound types. VISITRAK Depth can be used to measure wounds up to 15cm in depth, and is indicated for use with VISITRAK Grid and VISITRAK Digital.

### Contraindications

Do not bend VISITRAK Depth. Only insert the foam tipped end into the wound. Rest the tip on the wound bed. Care should be taken not to depress the wound bed.

S&N Code	Description	Pack Content	NHS Code	PIP Code
66000852	-	50	-	-



# Product Catalogue

## VISITRAK<sup>◇</sup> Digital Wound Measurement System

### Indications

Wound measurement can be used on all wound types. Indicated for use only with VISITRAK Grid and Depth.



S&N Code	Description	Pack Content	NHS Code	PIP Code
66000850	-	Device	-	-

## VISITRAK Grid Wound Measurement System

### Indications

For wound area measurement; can be used on all wound types. Indicated for use with VISITRAK Depth and VISITRAK Digital.

S&N Code	Description	Pack Content	NHS Code	PIP Code
66000851	-	50	-	-



# Product Catalogue

## ZIPZOC<sup>◇</sup>

20% w/w Zinc oxide impregnated medicated stocking

### Indications

To be used as an aid in the treatment of chronic leg ulcers. Where chronic venous insufficiency exists the stocking can be used as a primary contact layer under compression bandaging or hosiery.

For prescribing information see page 69 of booklet.

Legal category P



S&N Code	Description	Pack Content
66151550	Zinc Oxide Impregnated Medicated Stocking	4
66151540	Zinc Oxide Impregnated Medicated Stocking	10



# Product Catalogue

---

## ZOFF<sup>◇</sup>

### Adhesive remover wipes

#### Indications

Indicated for the removal of adhesive devices and residue from skin. Softens adhesive on tapes, films and appliances thereby minimising trauma caused by adhesive removal.



S&N Code	Description	Pack Content	NHS Code	PIP Code
5000	Adhesive Remover Wipes	20	-	-





# Product Catalogue

## Index

<b>Absorbent Dressing Containing Activated Charcoal Cloth</b>			
CARBONET®	28	ALLEVYN® Adhesive/Sacrum	12
<b>Adhesive remover</b>		ALLEVYN Ag Gentle/ Gentle Border/Gentle Border Heel	14
ZOFF®	64	ALLEVYN Compression	16
<b>Antimicrobials</b>		ALLEVYN Non-Adhesive	21
ALGISITE® Ag	9	ALLEVYN Gentle	17
ALLEVYN Ag	13	ALLEVYN Gentle Border	17
FLAMAZINE® Cream	34	ALLEVYN Gentle Border Heel	17
IODOFLEX®	37	ALLEVYN Gentle Border Lite	18
IODOSORB®	37	ALLEVYN Heel/Ag Heel	19
<b>Antimicrobials barrier dressings</b>		ALLEVYN Lite	20
ACTICOAT®	4	ALLEVYN Thin	24
ACTICOAT 7	4	ALLEVYN Plus Adhesive/Sacrum	22
ACTICOAT Flex 3/ACTICOAT Flex 7	6	ALLEVYN Tracheostomy	25
ACTICOAT Moisture control	7	ALGISITE® M	10
ACTICOAT Site	8	ALGOSTERIL	11
<b>Cavity dressings</b>		CUTINOVA® Hydro	31
ALGOSTERIL®	11	OPSITE® FLEXIGRID®	45
ALLEVYN Cavity	15	<b>Paste bandages/stockings</b>	
ALLEVYN Plus Cavity	23	ICHTHOPASTE®	35
INTRASITE® Conformable	36	VISCOPLASTE® PB7	60
<b>Compression bandaging</b>		ZIPZOC®	63
PROFORE®	53	<b>Post-operative dressings</b>	
PROGUIDE®	54	OPSITE Post-Op/Plus	47
<b>Fixation</b>		<b>Skin care</b>	
IV3000®	38	CADESORB®	27
LEUKOSTRIP®	40	NO-STING SKIN-PREP®	43
PRIMAFIX®	51	TRIPLE CARE®	58
<b>Gelling fibre dressing</b>		<b>Scar therapy</b>	
DURAFIBER®	32	CICA-CARE®	29
<b>Hydrogels</b>		<b>Topical anaesthetic</b>	
INTRASITE Conformable	36	AMETOP® Gel	26
INTRASITE Gel	36	<b>Transparent film dressings</b>	
<b>Incise drapes</b>		OPSITE FLEXIFIX	44
OPSITE Incise	46	OPSITE Post-Op Visible	48
<b>Low adherent dressings</b>		<b>Transparent film spray dressing</b>	
BACTIGRAS®	27	OPSITE Spray	49
CICA-PLAIE®	30	<b>Hydrosurgery and Negative Pressure Wound Systems</b>	
CUTIPLAST® Steril	31	PICO®	50
EXU-DRY®	33	RENASYS® EZ PLUS	55
JELONET®	39	RENASYS GO	56
MELOLIN®	41	VERSAJET®	59
MELOLITE®	42	<b>Wound measurement</b>	
PRIMAPORE®	52	VISITRAK® Depth	61
TRICOTEX®	57	VISITRAK Digital	62
<b>Moist wound healing dressings</b>		VISITRAK Grid	62



# Product Catalogue

## UK Prescribing information

### AMETOP® 40mg/g Gel

(Tetracaine base Ph. Eur. 4.0% w/w)

#### ABBREVIATED PRESCRIBING INFORMATION:

Please refer to Summary of Product Characteristics (SmPC) before prescribing.

#### PRESENTATIONS:

Please refer to Summary of Product Characteristics (SmPC) before prescribing.

#### INDICATIONS:

Percutaneous local anaesthetic to produce anaesthesia of the skin prior to venepuncture or venous cannulation.

#### DOSAGE AND ADMINISTRATION:

Apply contents of tube to centre of area to be anaesthetised; cover with an occlusive dressing. Contents from 1 tube can cover and anaesthetise an area of up to 30 sq.cm. (6cm x 5cm). Smaller areas of anaesthetised skin may be adequate in infants and small children. Adequate anaesthesia usually achieved after application time of 30 minutes for venepuncture and 45 minutes for venous cannulation after which gel should be removed with a gauze swab and the site prepared with an antiseptic wipe in the normal manner. Anaesthesia remains for 4-6 hours in most patients after a single application. If anaesthesia of more than one site is required: in adults and children > 5 years, 5 tubes; children between 1 month and 5 years no more than one tube may be applied at separate sites at the same time. Dose can be repeated after 5 hours if necessary. Max. cumulative dose in a 24 hour period ≤ 7 tubes (adults) and ≤ 2 tubes (children).

#### CONTRAINDICATIONS, PRECAUTIONS AND

##### WARNINGS:

Known hypersensitivity to any of the ingredients or to local anaesthetics of the ester type. Not for infants under 1 month of age. For premature babies it is recommended that it is not used until 1 month after their expected delivery date (ie 44 weeks gestation). Do not apply to broken skin, mucous membranes, eyes or ears. Only apply to intact, normal skin. Not to be taken internally. Do not instil into the middle ear or use for procedures which might involve penetration into the middle ear. Repeated exposure to Ametop Gel may increase the risk of sensitisation reactions to tetracaine. The systemic availability of tetracaine by percutaneous absorption of AMETOP Gel is low but exercise caution in patients with epilepsy. Not recommended for use in breast feeding mothers.

#### UNDESIRABLE EFFECTS:

Erythema, oedema, itching, blistering. Consult SmPC for further information about adverse events.

#### LEGAL CATEGORY: P

PRODUCT AUTHORISATION NO.: PL 14038/0001

COST: Carton of 12: £12.90; Carton of 12 with 24 occlusive dressings: £18.54

#### MARKETING AUTHORISATION HOLDER:

Smith & Nephew Pharmaceuticals Ltd.,  
101 Hessle Road, Hull, HU3 2BN, England  
Further information is available on request from the  
Marketing Authorisation Holder

DATE OF PREPARATION: August 2011

ASSET NO: 30650

Adverse events should be reported.  
Reporting forms and information can be found at  
[www.yellowcard.gov.uk](http://www.yellowcard.gov.uk).

Adverse events should also be reported to  
[Advice@smith-nephew.com](mailto:Advice@smith-nephew.com)



# UK Prescribing information

## FLAMAZINE® Cream

1% w/w (Silver sulfadiazine 1% w/w)

### ABBREVIATED PRESCRIBING INFORMATION:

Please refer to Summary of Product Characteristics (SmPC) before prescribing.

### PRESENTATIONS:

Tube containing 20g, or 50g and jars containing 250g or 500g silver sulfadiazine 1% w/w.

### INDICATIONS:

Prophylaxis and treatment of infection in burn wounds, leg ulcers and pressure sores. As an aid to the prophylaxis of infection in skin graft donor sites and extensive abrasions. Conservative management of finger-tip injuries where pulp, nail loss and/or partial loss of the distal phalanx has occurred.

### DOSAGE AND ADMINISTRATION:

To be applied topically. Burns and Leg Ulcers /Pressure Sores: Apply a layer approx. 3-5mm thick to the affected area using a sterile glove or spatula. Then cover area with an absorbent gauze dressing and support bandage where necessary. Dressing should be changed and FLAMAZINE cream applied at least every 24 hours in burn treatment, or at least three times weekly otherwise, and debridement carried out as necessary.

### FINGER-TIP INJURIES:

Achieve haemostasis of injury prior to application of a 3-5mm layer of FLAMAZINE cream. A conventional finger dressing may be used. Change dressings every 2-3 days.

### CONTRAINDICATIONS:

Should not be used at, or near term pregnancy, on premature infants or on newborn infants during the first months of life, in patients known to be hypersensitive to sulfonamides, silver sulfadiazine or other components of the cream.

### PRECAUTIONS AND WARNINGS:

Use with caution in nursing mothers and in the presence of significant hepatic or renal impairment and in individuals known to have G6P deficiency. Not recommended for use in leg or pressure ulcers that are very exudative. Use of FLAMAZINE cream may delay separation of burn eschar and may alter the appearance of burn wounds. Interactions: Silver may inactivate enzymatic debriding agents. With use in large area burns, effects of oral hypoglycaemic agents and phenytoin may be potentiated.

### UNDESIRABLE EFFECTS:

Application site rash, pruritus, burning, leucopenia, renal failure, argyria. Consult SmPC for further information about adverse events.

### LEGAL CATEGORY: POM

PRODUCT AUTHORISATION NO.: PL 13374/0006

COST: 20g tube: £2.91; 50g tube: £3.85; 250g jar: £10.32; 500g jar: £18.27

### MARKETING AUTHORISATION HOLDER:

Smith & Nephew Pharmaceuticals Ltd.,  
101 Hessle Road, Hull, HU3 2BN, England  
Further information is available on request from the Marketing Authorisation Holder

DATE OF PREPARATION: June 2011

ASSET NO: 30651

### Adverse events should be reported.

Reporting forms and information can be found at [www.yellowcard.gov.uk](http://www.yellowcard.gov.uk).

Adverse events should also be reported to Advice@smith-nephew.com



# UK Prescribing information

## **ZIPZOC® Medicated Stocking**

(Zinc oxide impregnated medicated stocking)

### **ABBREVIATED PRESCRIBING INFORMATION**

Please refer to Summary of Product Characteristics (SmPC) before prescribing.

#### **PRESENTATIONS:**

Sterile rayon stocking (80cm. x 7cm.) impregnated with 41.5g of ointment containing Zinc Oxide Ph Eur 20% w/w.

#### **INDICATIONS:**

To be used as an aid in the treatment of chronic leg ulcers. Where chronic venous insufficiency exists, the stocking can be used as a primary contact layer under compression bandaging or hosiery.

#### **DOSAGE AND ADMINISTRATION:**

For Cutaneous use. Zipzoc should be applied to cover the lower leg from the base of the toes to below the knee. All folds should be smoothed out. To protect clothes from Zipzoc a suitable outer bandage should be worn. If chronic venous insufficiency exists, Zipzoc can be used as a primary contact layer under suitable compression bandaging. The duration and frequency of application will depend on the clinical circumstances, but Zipzoc must be changed no less frequently than weekly.

#### **CONTRAINDICATIONS:**

Arterial leg ulcers. Known hypersensitivity to zinc oxide or excipients of the ointment.

#### **PRECAUTIONS AND WARNINGS:**

If the outer bandage is inappropriate or the stocking is not changed with sufficient frequency it can cause the stocking to dry out.

#### **UNDESIRABLE EFFECTS:**

The skin of leg ulcer patients is easily sensitised to some topical medications. Reactions reported include, rash, erythema, itching and maceration of the wound edge. Consult SmPC for further information about adverse events.

#### **LEGAL CATEGORY: P**

**PRODUCT AUTHORISATION NO.:** PL14038/0010

**COST:** Carton of 4: £12.50; Carton of 10: £31.26

#### **MARKETING AUTHORISATION HOLDER:**

Smith & Nephew Healthcare Ltd., Healthcare House, 101 Hessle Road, Hull HU3 2BN

**DATE OF PREPARATION:** August 2011

**ASSET NO:** 30648

**Adverse events should be reported.**

**Reporting forms and information can be found at [www.yellowcard.gov.uk](http://www.yellowcard.gov.uk).**

Adverse events should also be reported to

[Advice@smith-nephew.com](mailto:Advice@smith-nephew.com)







**Wound Management**

Smith & Nephew  
Healthcare Ltd  
Healthcare House  
101 Hessle Road  
Hull HU3 2BN

T 01482 222200  
F 01482 222211

[advice@smith-nephew.com](mailto:advice@smith-nephew.com)  
[www.smith-nephew.com](http://www.smith-nephew.com)

®Trademark of Smith & Nephew  
© Smith & Nephew January 2012  
33753