

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

For further information please call your local Smiths Medical distributor or Smiths Medical on +44 (0)1303 260551

**smiths**


**Smiths Medical International Ltd**

Hythe, Kent CT21 6JL UK

Tel: +44 (0)1303 260551 Fax: +44 (0)1303 266761

[www.smiths-medical.com](http://www.smiths-medical.com)

Smiths, OPTIVA, Jelco, PROTECTIV, PROTECTIV ACUVANCE, PROTECTIV Plus, Kids Cath, LOGICATH, LOGICATH *AgTive*, PORT-A-CATH, PORT-A-CATH II, PolyFlow, SECURE SITE, Fluoro-Free, CATH-SHIELD, ProPort, P.A.S PORT, P.A.S PORT Elite, POINT-LOK, WIDESITE, GRIPPER and GRIPPER PLUS, are registered trademarks of Smiths Group plc. Statlock and CV PLUS are trademarks of VENETEC, Inc. Literature No. LIT/VA2614  
©2006 Smiths Medical Family of Companies. All rights reserved.



**smiths**

Smiths Medical - a part of Smiths Group plc

smiths

## Vascular Access

product catalogue




**VASCULARACCESS**




<b>Peripheral IV Catheters</b>	<b>4-9</b>
- Standard IV Catheters	
- Safety IV Catheters	
<b>Central Venous Catheters</b>	<b>10-25</b>
- Kids Cath™	
- LOGICATH™	
- LOGICATH™ <i>AgTive</i>	
<b>Implantable Ports</b>	<b>26-33</b>
- PORT-A-CATH®	
- PORT-A-CATH® II	
- ProPort®	
- P.A.S PORT®	
<b>Huber Needles</b>	<b>34-35</b>
- GRIPPER®	
- GRIPPER PLUS®	
- Traditional Huber Needles	

# STANDARD IV CATHETERS



The Thin Walled Design



Precision Bevel

With Jelco® and OPTIVA® Smiths Medical offers a comprehensive range of conventional IV catheters. The catheters are available in straight hub, winged, and ported versions in a variety of lengths and gauge sizes to help meet patients' and therapeutic requirements.

OPTIVA® catheters are available in clear yet radio opaque polyurethane. The polyurethane used in OPTIVA® IV catheters is firm for controlled insertion but softens whilst indwelling to conform to the contour of the vein.

Jelco® IV catheters are available with a FEP polymer construction which also enables firm, controlled insertion.

Catheters feature a thin walled design that optimizes flow rate potential for each gauge size, allowing the clinician to choose smaller catheters without compromising on flow rates.

The precision beveled catheter tip and optimally engineered needle ensures smooth insertion, resulting in less patient trauma.

## OPTIVA® IV Catheter Straight Hub

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.2mm	50mm	315	4215	<b>5068</b>
● 16G	1.8mm	32mm	225		<b>5061</b>
● 16G	1.8mm	50mm	210	4217	<b>5062</b>
● 18G	1.3mm	45mm	100	4218	<b>5064</b>
● 18G	1.3mm	32mm	110	4219	<b>5065</b>
● 20G	1.1mm	45mm	60		<b>5069</b>
● 20G	1.1mm	32mm	65	4220	<b>5066</b>
● 22G	0.9mm	25mm	38	4222	<b>5060</b>
● 24G	0.7mm	19mm	24	4224	<b>5063</b>



OPTIVA®  
I.V. Catheters

## OPTIVA® W IV Catheter Winged

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.2mm	45mm	310	4514	<b>1114</b>
● 16G	1.8mm	45mm	200	4516	<b>1116</b>
● 18G	1.3mm	45mm	95	4518	<b>1118</b>
● 18G	1.3mm	32mm	95	4519	<b>1119</b>
● 20G	1.1mm	32mm	60	4520	<b>1120</b>
● 22G	0.9mm	25mm	35	4522	<b>1122</b>
● 24G	0.7mm	19mm	20	4524	<b>1124</b>



OPTIVA® W  
I.V. Catheters  
(with fixation wings)

## OPTIVA® 2 IV Catheter Injection Port

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.2mm	45mm	315	4514	<b>1014</b>
● 16G	1.8mm	45mm	210	4516	<b>1016</b>
● 18G	1.3mm	45mm	100	4518	<b>1018</b>
● 18G	1.3mm	32mm	110	4519	<b>1019</b>
● 20G	1.1mm	32mm	65	4520	<b>1020</b>
● 22G	0.9mm	25mm	38	4522	<b>1022</b>



OPTIVA® 2 I.V. Catheters  
(with injection valve)

# STANDARD IV CATHETERS

## JELCO® IV Catheter Straight Hub Radiopaque

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.1mm	50mm	315	4215	<b>4038</b>
● 16G	1.7mm	50mm	210	4217	<b>4032</b>
● 18G	1.3mm	45mm	100	4218	<b>4034</b>
● 18G	1.3mm	32mm	110	4219	<b>4035</b>
● 20G	1.1mm	45mm	60		<b>4039</b>
● 20G	1.1mm	32mm	65	4220	<b>4036</b>
● 22G	0.9mm	25mm	38	4222	<b>4030</b>
● 24G	0.7mm	19mm	24	4224	<b>4033</b>



JELCO®  
I.V. Catheters  
(radiopaque)

## JELCO® IV Catheter Straight Hub Clear Non Radiopaque

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.1mm	50mm	315	4215	<b>4018</b>
● 16G	1.7mm	50mm	210	4217	<b>4012</b>
● 18G	1.3mm	45mm	100	4218	<b>4014</b>
● 20G	1.1mm	32mm	65	4220	<b>4016</b>
● 22G	0.9mm	25mm	38	4222	<b>4010</b>
● 24G	0.7mm	19mm	24	4224	<b>4013</b>



JELCO®  
I.V. Catheters Clear  
(non radiopaque)

## JELCO® IV Catheter Winged Radiopaque

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 20G	1.1mm	32mm	65	4220	<b>4076</b>
● 22G	0.9mm	25mm	37	4222	<b>4070</b>
● 24G	0.7mm	16mm	24	4224	<b>4073</b>



JELCO® W  
I.V. Catheters with  
fixation wings  
(radiopaque)

## JELCO® IV Catheter Injection Port Radiopaque


Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.1mm	45mm	315	4514	<b>MR14</b>
● 16G	1.7mm	45mm	210	4516	<b>MR16</b>
● 18G	1.3mm	45mm	100	4518	<b>MR18</b>
● 18G	1.3mm	32mm	110	4519	<b>MR19</b>
● 20G	1.1mm	32mm	65	4520	<b>MR20</b>
● 22G	0.9mm	25mm	38	4522	<b>MR22</b>



JELCO® 2  
I.V. Catheters  
with injection valve  
(radiopaque)



GOES IN SHARP...  
...COMES OUT BLUNT



Needlestick still pose the greatest occupational risk for transmitting HIV and other blood borne pathogens. I.V. catheter needles (hollow bore, blood filled needles) are associated with the greatest incidence of high – risk needlestick injuries.<sup>1,2</sup>

Smiths Medical offers 2 safety IV catheters, PROTECTIV® ACUVANCE® and PROTECTIV®, designed to reduce percutaneous injuries. Both catheters feature polyurethane tubing, a thin wall design and optimally engineered needle and precision bevel catheter.

PROTECTIV® ACUVANCE® insertion requires minimal change in technique and no extra step.

Once the catheter hub is threaded off the needle, the unique safety mechanism automatically extends the blunt beyond the needle tip prior to removal from the patient, thereby significantly reducing the risk of needlestick injury.

PROTECTIV® IV catheter offers a user friendly design, allowing one-handed catheter insertion technique and controlled catheter placement.

#### References

1. Ippolito G et al. Prevention Management and Chemoprophylaxis of Occupational Exposure to HIV International Health Care Worker Safety Center and Advances in Exposure Prevention, University of Virginia 1997.
2. Advances in Exposure Prevention 1998; 3(4): 37-48.



# SAFETY IV CATHETERS

## ACUVANCE® Plus IV Safety Catheter Straight Hub PUR

	Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
●	14G	2.2mm	52mm	314		<b>3358</b>
●	16G	1.8mm	52mm	208		<b>3352</b>
●	16G	1.8mm	32mm	215		<b>3342</b>
●	18G	1.3mm	45mm	100	4218	<b>3354</b>
●	18G	1.3mm	32mm	110	4619	<b>3355</b>
●	20G	1.1mm	32mm	63	4220	<b>3356</b>
●	20G	1.1mm	25mm	65		<b>3357</b>
●	22G	0.9mm	25mm	38	4622	<b>3350</b>
●	24G	0.8mm	16mm	24	4524	<b>3353</b>



ACUVANCE® Plus

## ACUVANCE® Plus IV Safety Catheter Winged PUR

	Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
●	18G	1.4mm	32mm	110		<b>3375</b>
●	20G	1.2mm	32mm	63		<b>3376</b>
●	22G	1.0mm	25mm	38		<b>3370</b>
●	24G	0.8mm	16mm	24		<b>3373</b>

## ACUVANCE® 2 IV Safety Catheter Injection Port PUR

	Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
●	14G	2.2mm	45mm	315	4514	<b>1714</b>
●	16G	1.8mm	45mm	210	4516	<b>1716</b>
●	18G	1.3mm	45mm	100	4518	<b>1718</b>
●	18G	1.3mm	32mm	105	4519	<b>1719</b>
●	20G	1.1mm	32mm	65	4520	<b>1720</b>
●	22G	0.9mm	25mm	38	4522	<b>1722</b>



ACUVANCE® 2

## PROTECTIV® IV Safety Catheter Straight Hub FEP

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.2mm	32mm	325	P214	<b>3048</b>
● 16G	1.8mm	32mm	215	P216	<b>3042</b>
● 18G	1.3mm	32mm	110	4219	<b>3055</b>
● 20G	1.1mm	32mm	63	4220	<b>3056</b>
● 22G	0.9mm	25mm	38	4222	<b>3050</b>
● 24G	0.7mm	19mm	24	4224	<b>3053</b>



PROTECTIV® Plus

## PROTECTIV® W IV Safety Catheter Winged FEP

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.2mm	30mm	325	P214	<b>3078</b>
● 16G	1.8mm	30mm	215	P216	<b>3072</b>
● 18G	1.3mm	30mm	110	4219	<b>3075</b>
● 20G	1.1mm	30mm	63	4220	<b>3076</b>
● 22G	0.9mm	24mm	38	4222	<b>3070</b>
● 24G	0.7mm	15mm	24	4224	<b>3073</b>

## PROTECTIV® Plus IV Safety Catheter Straight Hub PUR

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.2mm	32mm	325	P214	<b>3068</b>
● 16G	1.8mm	32mm	215	P216	<b>3062</b>
● 18G	1.3mm	32mm	110	4219	<b>3065</b>
● 20G	1.1mm	32mm	63	4220	<b>3066</b>
● 22G	0.9mm	25mm	38	4222	<b>3060</b>
● 24G	0.7mm	19mm	24	4224	<b>3063</b>

## PROTECTIV® Plus W IV Safety Catheter Winged PUR

Gauge	O.D	Length	Flowrate ml/min	Obturator	Product Codes
● 14G	2.2mm	30mm	325	P214	<b>3088</b>
● 16G	1.8mm	30mm	215	P216	<b>3082</b>
● 18G	1.3mm	30mm	110	4219	<b>3085</b>
● 20G	1.1mm	30mm	63	4220	<b>3086</b>
● 22G	0.9mm	24mm	38	4222	<b>3080</b>
● 24G	0.7mm	15mm	24	4224	<b>3083</b>



Kids Cath™ Single

Therapeutic regimes of neonatal and paediatric patients demand the use of specialised catheters. Smiths Medical has focused on that problem and offers a range of catheters to meet the special demands on both the material and design.

All catheters within the Kids Cath™ range are made of polyurethane and are therefore DEHP free.

#### **Kids Cath™ Neo**

A 27G single lumen catheter with integrated guidewire specially developed for premature infants with a birth weight less than 1000g. This catheter is available in different length and can be inserted with a 24G IV catheter.

#### **Kids Cath™ Single**

For infants with a birth weight more than 1000g. This catheter is provided with a 22G sheath introducer and is available in different length.

#### **Kids Cath™ Umbilical**

A range of single, double and triple lumen venous umbilical catheters.

## Premature infants birth weight less 1000g 1 LUMEN Kids Cath™ Neo

Description	I.D / O.D	Lumen	Introducer	Case Qty	Product Codes
KidsCathNeo 1 Lumen 27G 12cm	0.18x0.4mm			10	<a href="#">MXC0712</a>
KidsCathNeo 1 Lumen 27G 20cm	0.18x0.4mm			10	<a href="#">MXC0720</a>
KidsCathNeo 1 Lumen 27G 30cm	0.18x0.4mm			10	<a href="#">MXC0730</a>



Kids Cath™ Neo

## Infants birth weight more than 1000g 1 LUMEN Kids Cath™ Single

Description	I.D / O.D	Lumen	Introducer	Case Qty	Product Codes
KidsCath Single1 Lumen 23G 20cm	0.4x0.6mm		22G (0.71mm) 3cm	10	<a href="#">MXC1020-2</a>
KidsCath Single1 Lumen 23G 33cm	0.4x0.6mm		22G (0.71mm) 3cm	10	<a href="#">MXC1033-2</a>
KidsCath Single1 Lumen 23G 50cm	0.4x0.6mm		22G (0.71mm) 3cm	10	<a href="#">MXC1050-2</a>

Kids Cath™ Single



## 1 LUMEN Kids Cath™ Umbilical

Description	O.D	Lumen	Introducer	Case Qty	Product Codes
KidsCath Umbilical 1 Lumen 4F 20cm	1.33mm			10	<a href="#">MXC2120-4</a>

Kids Cath™ Umbilical

## 2 LUMEN Kids Cath™ Umbilical

Description	O.D	Lumen	Introducer	Case Qty	Product Codes
KidsCath Umbilical 2 Lumen 4F 20cm	1.33mm	24G/24G		10	<a href="#">MXC2220-4</a>
KidsCath Umbilical 2 Lumen 4F 30cm	1.33mm	24G/24G		10	<a href="#">MXC2230-4</a>




## 3 LUMEN Kids Cath™ Umbilical

Description	O.D	Lumen	Introducer	Case Qty	Product Codes
KidsCath Umbilical 3 Lumen 4F 30cm	1.33mm	24G/25G/25G		10	<a href="#">MXC2320-4</a>

## Flowrate and Priming Volume Kids Cath™

Birth Weight	Priming Volume (Distal)	Priming Volume Blue (Prox)	Priming Volume Yellow (Middle)	Flowrate White ml/h (Distal)	Flowrate Blue ml/h (Prox)	Flowrate Yellow ml/h (Middle)	Product Codes
<1000g	0,05ml			7,7			<a href="#">MXC0712</a>
<1000g	0,06ml			5,0			<a href="#">MXC0720</a>
<1000g	0,08ml			3,8			<a href="#">MXC0730</a>
>1kg	0,1ml			138,5			<a href="#">MXC1020-2</a>
>1kg	0,11ml			92,3			<a href="#">MXC1033-2</a>
>1kg	0,13ml			65,8			<a href="#">MXC1050-2</a>
	0,19ml			689			<a href="#">MXC2120-4</a>
	0,2ml	0,17ml		520,9	167,2		<a href="#">MXC2220-4</a>
	0,22ml	0,2ml		373,9	95,2		<a href="#">MXC2230-4</a>
	0,22ml	0,17ml	0,18ml	503	128,7	193	<a href="#">MXC2320-4</a>

Gravity flow 100cm head height normal saline solution.



LOGICATH™

Smiths Medical offers a comprehensive range of acute central venous catheters. LOGICATH™, our conventional catheter and LOGICATH™ *AgTive* our antimicrobial catheter are available in different lengths, outside diameters, lumen and set configurations to meet both healthcare worker and patient demands.

The benefits of LOGICATH™ and LOGICATH™ *AgTive* central venous catheters are achieved using a range of features:

#### **Polyurethane Catheter Material**

The family of LOGICATH™ Central Venous Catheters are all made of polyurethane. At room temperature the catheter is stiff and difficult to kink and so easy to introduce. Once in the body it softens and reduces the potential for vessel trauma.

#### **Soft Catheter Tip**

The soft radio opaque conical catheter tip, which helps to reduce vessel irritation, forms a seamless junction with the catheter body. Most of the LOGICATH™ family central venous catheters are provided with a soft tip.

#### **Clearly Marked Extensions and Colour Coded Hubs**

The clearly marked extensions allow catheter manipulation to be performed well away from the insertion site, thus reducing the risk of infection. Different lengths of extensions and colour coded hubs allow easy and reliable lumen identification.

#### **Dilator**

The specially formed dilator tip provides a smooth transition with the guidewire to reduce the risk of skin irritation and dilator tip breakage. The sharp dilator tip reduces the need for a scalpel to enlarge the access site. Dilators are tailored to each catheter size.



Sets and Kits

### Y-Site Adaptor

The special Y-site adaptor fits between the syringe and introducer needle during catheter introduction. The guidewire is introduced through the side arm of the Y, thus eliminating the need to remove the syringe. This prevents blood loss and exposure to blood borne pathogens.

### Sets and Kits

LOGICATH™ and LOGICATH™ *AgTive* central venous catheters are available in different set and kit configurations. The standard set provides all components necessary for the placement of a central venous catheter. The kit version features additional components for the insertion of a central venous catheter to maximise infection prevention by a more procedural approach.

### Standard Set Content:

- LOGICATH™ conventional or antimicrobial catheter.
- Dilator.
- Introducer needle.
- Guidewire.
- Syringe.
- Catheter fastener and clamp.
- Injection caps.
- Y-site adaptor (for 14G, 16G, 7F, 8.5F and 12F catheters only).
- Scalpel.
- Traceability label.

These standard sets are also available as cost effective yet high quality MINI sets without a scalpel and Y-site adaptor for most of our conventional LOGICATH™ catheters.

Standard catheter sets without a Y-site adaptor are also available for a majority of catheters.

### Standard Kit Content:

- LOGICATH™ conventional or antimicrobial catheter.
- Dilator.
- Introducer needle.
- Guidewire.
- Syringe.
- Catheter fastener and clamp.
- Injection caps.
- Y-site adaptor (for 14G, 16G, 7F, 8.5F and 12F catheters only).
- Scalpel.
- Traceability label.
- 90x90cm fenestrated CVC drape.
- Gauze pads.
- Initial 24Hr gauze dressing.
- Pre-mounted 4-way blue lipid resistant stopcocks.
- Needle container
- STATLOCK® CV PLUS.

## LOGICATH™ Adult Catheters

### 1 LUMEN LOGICATH™ 16G without extension

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 1 Lumen 16G 16cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x16x16</a>
LOGICATH™ Set 1 Lumen 16G 16cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x16x16G</a>
LOGICATH™ Set 1 Lumen 16G 20cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x20x16</a>
LOGICATH™ Mini Set 1 Lumen 16G 20cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x20x16M</a>
LOGICATH™ Set 1 Lumen 16G 20cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x20x16G</a>
LOGICATH™ Set 1 Lumen 16G 30cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x30x16</a>

### with extension

LOGICATH™ Set 1 Lumen 16G 16cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x16x16V</a>
LOGICATH™ Set 1 Lumen 16G 20cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x20x16V</a>
LOGICATH™ Mini Set 1 Lumen 16G	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x20x16VM</a>
LOGICATH™ Set 1 Lumen 16G 30cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x30x16V</a>

### 1 LUMEN LOGICATH™ Kit 16G

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 1 Lumen 16G 20cm	1.1x1.7mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x20x16V-E</a>

### 1 LUMEN LOGICATH™ 14G without extension

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 1 Lumen 14G 16cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x16x14</a>
LOGICATH™ Set 1 Lumen 14G 20cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x20x14</a>
LOGICATH™ Set 1 Lumen 14G 20cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x20x14G</a>
LOGICATH™ Set 1 Lumen 14G 30cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x30x14</a>

### with extension

LOGICATH™ Set 1 Lumen 14G 16cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x16x14V</a>
LOGICATH™ Set 1 Lumen 14G 20cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x20x14V</a>
LOGICATH™ Set 1 Lumen 14G 20cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x20x14VG</a>
LOGICATH™ Mini Set 1 Lumen 14G 20cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x20x14VM</a>
LOGICATH™ Set 1 Lumen 14G 30cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x30x14V</a>

### 1 LUMEN LOGICATH™ Kit 14G

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 1 Lumen 14G 20cm	1.4x2.1mm		18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	7.5F (2.4mm) 8cm	10	<a href="#">MXA231x20x14V-E</a>

2 LUMEN LOGICATH™ 7F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 2 Lumen 7F 16cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x16x70</a>
LOGICATH™ Set 2 Lumen 7F 16cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x16x70G</a>
LOGICATH™ Set 2 Lumen 7F 16cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x16x70M</a>
LOGICATH™ Set 2 Lumen 7F 20cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x20x70</a>
LOGICATH™ Set 2 Lumen 7F 20cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x20x70G</a>
LOGICATH™ Mini Set 2 Lumen 7F 20cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x20x70M</a>
LOGICATH™ Set 2 Lumen 7F 30cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x30x70</a>
LOGICATH™ Set 2 Lumen 7F 30cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x30x70G</a>
LOGICATH™ Mini Set 2 Lumen 7F 30cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x30x70M</a>
LOGICATH™ Set 2 Lumen 7F 16G/16G	2.34mm	16G/16G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x16x75</a>
LOGICATH™ Set 2 Lumen 7F 16G/16G	2.34mm	16G/16G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x20x75</a>
LOGICATH™ Set 2 Lumen 7F 16G/16G	2.34mm	16G/16G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x30x75</a>

2 LUMEN LOGICATH™ Kit 7F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 2 Lumen 7F 16cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x16x70-E</a>
LOGICATH™ Kit 2 Lumen 7F 20cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x20x70-E</a>
LOGICATH™ Kit 2 Lumen 7F 30cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA232x30x70-E</a>

2 LUMEN LOGICATH™ 8.5F High Flow

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set Lumen 8.5F 16cm	2.77mm	14G/14G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	<a href="#">MXA232x16x85</a>
LOGICATH™ Set 2 Lumen 8.5F 20cm	2.77mm	14G/14G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	<a href="#">MXA232x20x85</a>
LOGICATH™ Set 2 Lumen 8.5F 20cm	2.77mm	14G/14G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	<a href="#">MXA232x20x85G</a>
LOGICATH™ Set 2 Lumen 8.5F 25cm	2.77mm	14G/14G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	<a href="#">MXA232x25x85</a>

2 LUMEN LOGICATH™ 12F High Flow

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 2 Lumen 12F 16cm	4.00mm	12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F ( 4mm) 8cm	10	<a href="#">MXA232x16x12</a>
LOGICATH™ Set 2 Lumen 12F 20cm	4.00mm	12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F ( 4mm) 8cm	10	<a href="#">MXA232x20x12</a>

2 LUMEN LOGICATH™ 12F Hemo Dialysis

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 2 Lumen 12F 16cm	4.00mm	12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F ( 4mm) 8cm	10	<a href="#">MXA232x16x12S</a>
LOGICATH™ Set 2 Lumen 12F 20cm	4.00mm	12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F ( 4mm) 8cm	10	<a href="#">MXA232x20x12S</a>



## 3 LUMEN LOGICATH™ 7F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 3 Lumen 7F 16cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x16x70</a>
LOGICATH™ Set 3 Lumen 7F 16cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x16x70G</a>
LOGICATH™ mini Set 3 Lumen 7F 16cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x16x70M</a>
LOGICATH™ Set 3 Lumen 7F 20cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x20x70</a>
LOGICATH™ Set 3 Lumen 7F 20cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x20x70G</a>
LOGICATH™ Mini Set 3 Lumen 7F 20	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x20x70M</a>
LOGICATH™ Set 3 Lumen 7F 30cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x30x70</a>
LOGICATH™ Set 3 Lumen 7F 30cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x30x70G</a>
LOGICATH™ Mini Set 3 Lumen 7F 30cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x30x70M</a>

## 3 LUMEN LOGICATH™ Kit 7F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 3 Lumen 7F 16cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x16x70-E</a>
LOGICATH™ Kit 3 Lumen 7F 20cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x20x70-E</a>
LOGICATH™ Kit 3 Lumen 7F 30cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA233x30x70-E</a>

## 3 LUMEN LOGICATH™ 8.5F High Flow

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 3 Lumen 8.5F 16cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA233x16x85</a>
LOGICATH™ Set 3 Lumen 8.5F 16cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA233x16x85G</a>
LOGICATH™ Set 3 Lumen 8.5F 20cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA233x20x85</a>
LOGICATH™ Set 3 Lumen 8.5F 20cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA233x20x85G</a>
LOGICATH™ Set 3 Lumen 8.5F 25cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA233x25x85</a>
LOGICATH™ Set 3 Lumen 8.5F 30cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA233x30x85</a>

## 3 LUMEN LOGICATH™ Kit 8.5F High Flow

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 3 Lumen 8.5F 16cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA233x16x85-E</a>
LOGICATH™ Kit 3 Lumen 8.5F 20cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA233x20x85-E</a>

## 3 LUMEN LOGICATH™ 12F High Flow

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 3 Lumen 12F 16cm	4.00mm	16G/12G/12G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	12F (4mm) 8cm	10	<a href="#">MXA233x16x12</a>
LOGICATH™ Set 3 Lumen 12F 20cm	4.00mm	16G/12G/12G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	12F (4mm) 8cm	10	<a href="#">MXA233x20x12</a>
LOGICATH™ Set 3 Lumen 12F 25cm	4.00mm	16G/12G/12G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	12F (4mm) 8cm	10	<a href="#">MXA233x25x12</a>

3 LUMEN LOGICATH™ 12F Hemo Dialysis

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 3 Lumen 12F 16cm	4.00mm	16G/12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F (4mm) 8cm	10	MXA233x16x12S
LOGICATH™ Set 3 Lumen 12F 20cm	4.00mm	16G/12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F (4mm) 8cm	10	MXA233x20x12S

3 LUMEN LOGICATH™ Kit 12F High Flow

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 3 Lumen 12F 16cm HF	4.00mm	16G/12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F (4mm) 8cm	10	MXA233x16x12-E
LOGICATH™ Kit 3 Lumen 12F 20cm HF	4.00mm	16G/12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F (4mm) 8cm	10	MXA233x20x12-E

4 LUMEN LOGICATH™ 8.5F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 4 Lumen 8.5F 16cm	2.77mm	16G/18G/ 18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA234x16x85
LOGICATH™ Set 4 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA234x20x85G
LOGICATH™ Set 4 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA234x20x85
LOGICATH™ Set 4 Lumen 8.5F 25cm	2.77mm	16G/18G/ 18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA234x25x85
LOGICATH™ Set 4 Lumen 8.5F 30cm	2.77mm	16G/18G/ 18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA234x30x85

4 LUMEN LOGICATH™ Kit 8.5F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 4 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA234x20x85-E
LOGICATH™ Kit 4 Lumen 8.5F 30cm	2.77mm	16G/18G/ 18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA234x30x85-E

5 LUMEN LOGICATH™ 8.5F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 5 Lumen 8.5F 16cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x16x85
LOGICATH™ Set 5 5 Lumen 8.5F 16cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x16x85G
LOGICATH™ Set 5 5 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x20x85
LOGICATH™ Set 5 5 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x20x85G
LOGICATH™ Set 5 5 Lumen 8.5F 25cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x25x85
LOGICATH™ Set 5 5 Lumen 8.5F 30cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x30x85

5 LUMEN LOGICATH™ KIT 8.5F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 5 5 Lumen 8.5F 16cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x16x85-E
LOGICATH™ Kit 5 5 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x20x85-E
LOGICATH™ Kit 5 5 Lumen 8.5F 30cm	2.77mm	16G/18G/ 18G/18G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	MXA235x30x85-E

## LOGICATH™ Paediatric Catheters

### 1 LUMEN LOGICATH™ with extension

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 1 Lumen 24G 10cm	0.5x0.7mm		20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	4.5F (1.5mm) 6cm	10	<a href="#">MXA231x10x24V</a>
LOGICATH™ Set 1 Lumen 24G 16cm	0.5x0.7mm		20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	4.5F (1.5mm) 6cm	10	<a href="#">MXA231x16x24V</a>
LOGICATH™ Set 1 Lumen 24G 20cm	0.5x0.7mm		20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	4.5F (1.5mm) 6cm	10	<a href="#">MXA231x20x24V</a>
LOGICATH™ Set 1 Lumen 20G 10cm	0.6x0.9mm		20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	4.5F (1.5mm) 6cm	10	<a href="#">MXA231x10x20V</a>
LOGICATH™ Set 1 Lumen 20G 16cm	0.6x0.9mm		20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	4.5F (1.5mm) 6cm	10	<a href="#">MXA231x16x20V</a>
LOGICATH™ Set 1 Lumen 20G 20cm	0.6x0.9mm		20G (0.89mm) 4cm	.018" (.46mm) 60cm straight	4.5F (1.5mm) 6cm	10	<a href="#">MXA231x20x20V</a>
LOGICATH™ Set 1 1 Lumen 20G 16cm	0.6x0.9mm		peelaway cannula			10	<a href="#">MXA231x16x20VS</a>
LOGICATH™ Set 1 Lumen 20G 20cm	0.6x0.9mm		peelaway cannula			10	<a href="#">MXA231x20x20VS</a>

### without extension

LOGICATH™ Set 1 Lumen 18G 16cm	0.8x1.3mm		20G (0.89mm) 4cm	.025" (.64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x16x18</a>
LOGICATH™ Set 1 Lumen 18G 20cm	0.8x1.3mm		20G (0.89mm) 4cm	.025" (.64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x20x18</a>
LOGICATH™ Set 1 Lumen 18G 30cm	0.8x1.3mm		20G (0.89mm) 4cm	.025" (.64mm) 70cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x30x18</a>

### with extension

LOGICATH™ Set 1 Lumen 18G 16cm	0.8x1.3mm		20G (0.89mm) 4cm	.025" (.64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x16x18V</a>
LOGICATH™ Set 1 Lumen 18G 20cm	0.8x1.3mm		20G (0.89mm) 4cm	.025" (.64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x20x18V</a>
LOGICATH™ Set 1 Lumen 18G 30cm	0.8x1.3mm		20G (0.89mm) 4cm	.025" (.64mm) 70cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA231x30x18V</a>

### 2 LUMEN LOGICATH™ 4F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 2 Lumen 4F 5cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x05x40</a>
LOGICATH™ Set 2 Lumen 4F 5cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J-tip 1.5mm	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x05x40J</a>
LOGICATH™ Set 2 Lumen 4F 8cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x08x40</a>
LOGICATH™ Set 2 Lumen 4F 8cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J-tip 1.5mm	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x08x40J</a>
LOGICATH™ Set 2 Lumen 4F 10cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x10x40</a>
LOGICATH™ Set 2 Lumen 4F 10cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J-tip 1.5mm	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x10x40J</a>
LOGICATH™ Set 2 Lumen 4F 12cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x12x40</a>
LOGICATH™ Set 2 Lumen 4F 16cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x16x40</a>
LOGICATH™ Set 2 Lumen 4F 16cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J-tip 1.5mm	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x16x40J</a>
LOGICATH™ Set 2 Lumen 4F 20cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm straight	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x20x40</a>
LOGICATH™ Set 2 Lumen 4F 20cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J-tip 1.5mm	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x20x40J</a>

### 2 LUMEN LOGICATH™ Kit 4F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 2 Lumen 4F 8cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J tip 1.5mm	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x08x40-E</a>
LOGICATH™ Kit 2 Lumen 4F 12cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J tip 1.5mm	5.5F (1.8mm) 6cm	10	<a href="#">MXA232x12x40-E</a>

2 LUMEN LOGICATH™ 5.5F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 2 Lumen 5.5F 8cm	1.83mm	18G/18G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA232x08x55</a>
LOGICATH™ Set 2 Lumen 5.5F 10cm	1.83mm	18G/18G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA232x10x55</a>
LOGICATH™ Set 2 Lumen 5.5F 13cm	1.83mm	18G/18G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA232x13x55</a>
LOGICATH™ Set 2 Lumen 5.5F 16cm	1.83mm	18G/18G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA232x16x55</a>
LOGICATH™ Set 2 Lumen 5.5F 20cm	1.83mm	18G/18G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA232x20x55</a>

2 LUMEN LOGICATH™ Kit 5.5F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 2 Lumen 5.5F 16cm	1.83mm	18G/18G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA232x16x55-E</a>

3 LUMEN LOGICATH™ 4F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 3 Lumen 4F 5cm	1.35mm	20G/24G/24G	20G ( 0.89mm) 4cm	.018" ( .46mm) 40cm straight	5.5F ( 1.8mm) 6cm	10	<a href="#">MXA233x05x40</a>
LOGICATH™ Set 3 Lumen 4F 8cm	1.35mm	20G/24G/24G	20G ( 0.89mm) 4cm	.018" ( .46mm) 40cm straight	5.5F ( 1.8mm) 6cm	10	<a href="#">MXA233x08x40</a>
LOGICATH™ Set 3 Lumen 4F 10cm	1.35mm	20G/24G/24G	20G ( 0.89mm) 4cm	.018" ( .46mm) 40cm straight	5.5F ( 1.8mm) 6cm	10	<a href="#">MXA233x10x40</a>
LOGICATH™ Set 3 Lumen 4F 12cm	1.35mm	20G/24G/24G	20G ( 0.89mm) 4cm	.018" ( .46mm) 40cm straight	5.5F ( 1.8mm) 6cm	10	<a href="#">MXA233x12x40</a>
LOGICATH™ Set 3 Lumen 4F 16cm	1.35mm	20G/24G/24G	20G ( 0.89mm) 4cm	.018" ( .46mm) 40cm straight	5.5F ( 1.8mm) 6cm	10	<a href="#">MXA233x16x40</a>
LOGICATH™ Set 3 Lumen 4F 20cm	1.35mm	20G/24G/24G	20G ( 0.89mm) 4cm	.018" ( .46mm) 40cm straight	5.5F ( 1.8mm) 6cm	10	<a href="#">MXA233x20x40</a>

3 LUMEN LOGICATH™ Kit 4F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 3 Lumen 4F 8cm	1.35mm	20G/24G/24G	20G ( 0.89mm) 4cm	.018" ( .46mm) 40cm J tip 1.5mm	5.5F ( 1.8mm) 6cm	10	<a href="#">MXA233x08x40-E</a>

3 LUMEN LOGICATH™ 5.5F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Set 3 Lumen 5.5F 8cm	1.83mm	18G/22G/22G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA233x08x55</a>
LOGICATH™ Set 3 Lumen 5.5F 10cm	1.83mm	18G/22G/22G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA233x10x55</a>
LOGICATH™ Set 3 Lumen 5.5F 13cm	1.83mm	18G/22G/22G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA233x13x55</a>
LOGICATH™ Set 3 Lumen 5.5F 16cm	1.83mm	18G/22G/22G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA233x16x55</a>
LOGICATH™ Set 3 Lumen 5.5F 20cm	1.83mm	18G/22G/22G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA233x20x55</a>
LOGICATH™ Set 3 Lumen 5.5F 30cm	1.83mm	18G/22G/22G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA233x30x55</a>

3 LUMEN LOGICATH™ Kit 5.5F

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ Kit 3 Lumen 5.5F 8cm	1.83mm	18G/22G/22G	20G ( 0.89mm) 4cm	.025" ( .64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	<a href="#">MXA233x08x55-E</a>



LOGICATH™ *AgTive*

Smith Medical's innovative technology offers easy, safe and cost effective use of central venous catheters.

#### **Technology makes the difference**

Polyurethanes are hygroscopic and strongly attract water. The unique technology used in LOGICATH™ *AgTive* incorporates nano-sized particles of silver into the polyurethane, allowing the release of bactericidal concentrations of silver ions from both the inner and outer surface of the catheter when in contact with fluids.

LOGICATH™ *AgTive* utilises this technology within the catheter tubing and the delta to maximise infection prevention.

LOGICATH™ *AgTive* Adult Catheters

1 LUMEN LOGICATH™ *AgTive* antimicrobial

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Set 1 Lumen 14G 20cm	1.4x2.1mm		18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	7.5F ( 2.4mm) 8cm	10	<a href="#">MXA431x20x14V</a>
LOGICATH™ <i>AgTive</i> Set 1 Lumen 14G 30cm	1.4x2.1mm		18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	7.5F ( 2.4mm) 8cm	10	<a href="#">MXA431x30x14V</a>
LOGICATH™ <i>AgTive</i> Set 1 Lumen 16G 20cm	1.1x1.7mm		18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	6.5F ( 2.1mm) 6cm	10	<a href="#">MXA431X20X16V</a>
LOGICATH™ <i>AgTive</i> Set 1 Lumen 16G 30cm	1.1x1.7mm		18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	6.5F ( 2.1mm) 6cm	10	<a href="#">MXA431X30X16V</a>

1 LUMEN LOGICATH™ *AgTive* Kit antimicrobial

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 1 Lumen 14G 20cm	1.4x2.1mm		18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	7.5F ( 2.4mm) 8cm	10	<a href="#">MXA431x20x14V-E</a>
LOGICATH™ <i>AgTive</i> Kit 1 Lumen 16G 20cm	1.1x1.7mm		18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	6.5F ( 2.1mm) 6cm	10	<a href="#">MXA431x20x16V-E</a>

2 LUMEN LOGICATH™ *AgTive* 7F antimicrobial

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Set 2 Lumen 7F 16cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA432x16x70</a>
LOGICATH™ <i>AgTive</i> Set 2 Lumen 7F 20cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA432x20x70</a>
LOGICATH™ <i>AgTive</i> Set 2 Lumen 7F 30cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA432x30x70</a>

2 LUMEN LOGICATH™ *AgTive* Kit 7F antimicrobial

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 2 Lumen 7F 16cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA432x16x70-E</a>
LOGICATH™ <i>AgTive</i> Kit 2 Lumen 7F 20cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA432x20x70-E</a>
LOGICATH™ <i>AgTive</i> Kit 2 Lumen 7F 30cm	2.34mm	14G/18G	18G ( 1.27mm) 7cm	.035" ( .89mm) 70cm J-tip 3mm	8.5F ( 2.7mm) 8cm	10	<a href="#">MXA432x30x70-E</a>

2 LUMEN LOGICATH™ *AgTive* Kit 8.5F antimicrobial

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 2 Lumen 8.5F 20cm HF	2.77mm	14G/14G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	9.5F ( 3.0mm) 8cm	10	<a href="#">MXA432x20x85-E</a>

2 LUMEN LOGICATH™ *AgTive* 12F antimicrobial

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Set 12F 2 Lumen 16cm HF	4,00mm	12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F ( 4mm) 8cm	10	<a href="#">MXA432x16x12</a>
LOGICATH™ <i>AgTive</i> Set 12F 2 Lumen 20cm HF	4,00mm	12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F ( 4mm) 8cm	10	<a href="#">MXA432x20x12</a>

2 LUMEN LOGICATH™ *AgTive* Kit 12F antimicrobial

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 12F 2 Lumen 16cm HF	4,00mm	12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F ( 4mm) 8cm	10	<a href="#">MXA432x16x12-E</a>
LOGICATH™ <i>AgTive</i> Kit 12F 2 Lumen 20cm HF	4,00mm	12G/12G	18G ( 1.27mm) 7cm	.035" ( .89mm) 60cm J-tip 3mm	12F ( 4mm) 8cm	10	<a href="#">MXA432x20x12-E</a>

## 3 LUMEN LOGICATH™ *AgTive* 7F antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Set 3 Lumen 7F 16cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA433x16x70</a>
LOGICATH™ <i>AgTive</i> Set 3 Lumen 7F 20cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA433x20x70</a>
LOGICATH™ <i>AgTive</i> Set 3 Lumen 7F 30cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA433x30x70</a>

## 3 LUMEN LOGICATH™ *AgTive* Kit 7F antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 3 Lumen 7F 16cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA433x16x70-E</a>
LOGICATH™ <i>AgTive</i> Kit 3 Lumen 7F 20cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA433x20x70-E</a>
LOGICATH™ <i>AgTive</i> Kit 3 Lumen 7F 30cm	2.34mm	16G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	8.5F (2.7mm) 8cm	10	<a href="#">MXA433x30x70-E</a>

## 3 LUMEN LOGICATH™ *AgTive* Kit 8.5F antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 3 Lumen 8.5F 20cm	2.77mm	14G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA433x20x85-E</a>

## 4 LUMEN LOGICATH™ *AgTive* 8.5F antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Set 8.5F 4 Lumen 16cm F	2.77mm	16G/18G/ 18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA434x16x85</a>
LOGICATH™ <i>AgTive</i> Set 4 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA434x20x85</a>
LOGICATH™ <i>AgTive</i> Set 4 Lumen 8.5F 25cm	2.77mm	16G/18G/ 18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA434x25x85</a>
LOGICATH™ <i>AgTive</i> Set 4 Lumen 8.5F 30cm	2.77mm	16G/18G/ 18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA434x30x85</a>

## 4 LUMEN LOGICATH™ *AgTive* Kit 8.5F antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 4 Lumen 8.5F 16cm	2.77mm	16G/18G/ 18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA434x16x85-E</a>
LOGICATH™ <i>AgTive</i> Kit 4 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA434x20x85-E</a>
LOGICATH™ <i>AgTive</i> Kit 4 Lumen 8.5F 25cm	2.77mm	16G/18G/ 18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA434x25x85-E</a>
LOGICATH™ <i>AgTive</i> Kit 4 Lumen 8.5F 30cm	2.77mm	16G/18G/ 18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA434x30x85-E</a>

## 5 LUMEN LOGICATH™ *AgTive* 8.5F antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Set 5 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA435x20x85</a>
LOGICATH™ <i>AgTive</i> Set 5 Lumen 8.5F 25cm	2.77mm	16G/18G/ 18G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA435x25x85</a>
LOGICATH™ <i>AgTive</i> Set 5 Lumen 8.5F 30cm	2.77mm	16G/18G/ 18G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 70cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA435x30x85</a>

## 5 LUMEN LOGICATH™ *AgTive* Kit 8.5F antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 5 Lumen 8.5F 16cm	2.77mm	16G/18G/ 18G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA435x16x85-E</a>
LOGICATH™ <i>AgTive</i> Kit 5 Lumen 8.5F 20cm	2.77mm	16G/18G/ 18G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA435x20x85-E</a>
LOGICATH™ <i>AgTive</i> Kit 5 Lumen 8.5F 25cm	2.77mm	16G/18G/ 18G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA435x25x85-E</a>
LOGICATH™ <i>AgTive</i> Kit 5 Lumen 8.5F 30cm	2.77mm	16G/18G/ 18G/18G/18G	18G (1.27mm) 7cm	.035" (.89mm) 60cm J-tip 3mm	9.5F (3.0mm) 8cm	10	<a href="#">MXA435x30x85-E</a>

# LOGICATH™ AND LOGICATH™ AGTIVE ACCESSORIES

## LOGICATH™ *AgTive* Paediatric Catheters

### 1 LUMEN LOGICATH™ *AgTive* Antimicrobial

Description	I.D / O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Set 1 Lumen 20G 10cm	0.6x0.9mm		20G (0.89mm) 4cm	.018" (.46mm) 60cm	4.5F (1.5mm) 6cm	10	MXA431x10x20V
LOGICATH™ <i>AgTive</i> Kit 20G 1 Lumen 20cm	0.6x0.9mm		20G (0.89mm) 4cm			10	MXA431x20x20VS-E

### 2 LUMEN LOGICATH™ *AgTive* Kit 4F Antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Kit 4F 2 Lumen 8cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J tip 1.5mm	5.5F (1.8mm) 6cm	10	MXA432x08x40-E
LOGICATH™ <i>AgTive</i> Kit 4F 2 Lumen 12cm	1.35mm	18G/22G	20G (0.89mm) 4cm	.018" (.46mm) 40cm J tip 1.5mm	5.5F (1.8mm) 6cm	10	MXA432x12x40-E

### 2 LUMEN LOGICATH™ *AgTive* Kit 5.5F Antimicrobial

Description	O.D	Lumen	Introducer	Guidewire	Dilator	Case Qty	Product Codes
LOGICATH™ <i>AgTive</i> Set 5.5F2 Lumen 10cm	1.83mm	18G/18G	20G (0.89mm) 4cm	.025" (.64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	MXA432X10X55
LOGICATH™ <i>AgTive</i> Set 5.5F 2 Lumen 16cm	1.83mm	18G/18G	20G (0.89mm) 4cm	.025" (.64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	MXA432x16x55
LOGICATH™ <i>AgTive</i> Set 5.5F 2 Lumen 20cm	1.83mm	18G/18G	20G (0.89mm) 4cm	.025" (.64mm) 60cm J-tip 3mm	6.5F (2.1mm) 6cm	10	MXA432x20x55

## Accessories

### LOGICATH™ and LOGICATH™ *AgTive* Dilators

Description	Size and Length	Qty	Product Codes
LOGICATH™ Dilator	4.5F / 1.6mm, 6cm	25	MXA999x600x45
LOGICATH™ Dilator	5.5F / 1.8mm, 6cm	25	MXA999x600x55
LOGICATH™ Dilator	6.5F / 2.1mm, 6cm	25	MXA999x600x65
LOGICATH™ Dilator	7.5F / 2.4mm, 8cm	25	MXA999x800x75
LOGICATH™ Dilator	8.5F / 2.7mm, 8cm	25	MXA999x800x85
LOGICATH™ Dilator	9.5F / 3.0mm, 8cm	25	MXA999x800x95
LOGICATH™ Dilator	12F / 4.0mm, 8cm	25	MXA999x800x12

### LOGICATH™ and LOGICATH™ *AgTive* Introducer Needle

Description	Size and Length	Qty	Product Codes
Introducer Needle	18G / 1.27mm, 7cm	25	MXA999x7x18
Introducer Needle	18G / 1.27mm, 5cm	25	MXA999x5x18

### LOGICATH™ and LOGICATH™ *AgTive* Y-Connector

Description	Qty	Product Codes
Y-connector	50	MXA999-2020

### LOGICATH™ and LOGICATH™ *AgTive* Fixation Devices

Description	Qty	Product Codes
CVC Fixation Device 18G and 4F catheter	10	MXA999x40x00
CVC Fixation Device 16G and 5.5F catheter	10	MXA999x55x00
CVC Fixation Device 14G catheter	10	MXA999x14x00
CVC Fixation Device 7F catheter	10	MXA999x70x00
CVC Fixation Device 8.5F catheter	10	MXA999x85x00

### LOGICATH™ and LOGICATH™ *AgTive* Guidewires

Description	Size and Length	Qty	Product Codes
Guidewire .018"	0.46mm, 40cm straight	20	MXA999x40x18
Guidewire .018"	0.46mm, 60cm straight	20	MXA999x60x18
Guidewire .025"	0.64mm, 60cm J-tip 3mm	20	MXA999x60x25
Guidewire .025"	0.64mm, 70cm J-tip 3mm	20	MXA999x70x25
Guidewire .035"	0.89mm, 60cm J-tip 3mm	20	MXA999x60x35
Guidewire .035"	0.89mm, 70cm J-tip 3mm	20	MXA999x70x35



## LOGICATH™/LOGICATH™ *AgTive* Catheter Specifications

### 1 LUMEN

Gauge (G) Distal	Length (cm) Distal	Flow Rate (ml/hour) Distal	Priming Volume (ml)	Product Codes
14G	16	7460	0.54	<a href="#">MXA231X16X14V</a>
14G	20	6870	0.64	<a href="#">MXA231X20X14V</a>
14G	20	7850	0.41	<a href="#">MXA231X20X14</a>
16G	16	4480	0.45	<a href="#">MXA231X16X16V</a>
16G	20	3820	0.49	<a href="#">MXA231X20X16V</a>
16G	30	2960	0.59	<a href="#">MXA231X30X16V</a>
16G	20	4250	0.29	<a href="#">MXA231X20X16</a>
18G	16	1760	0.37	<a href="#">MXA231X16X18V</a>
18G	20	1510	0.38	<a href="#">MXA231X20X18V</a>
18G	20	1530	0.21	<a href="#">MXA231X20X18</a>
20G	16	560	0.33	<a href="#">MXA231X16X20V</a>
20G	20	500	0.36	<a href="#">MXA231X20X20V</a>

### 2 LUMEN

Gauge (G)		Length (cm)	French (Fr)	Flow Rate (ml/hour)		Priming Volume (ml)		Product Codes
Distal	Proximal			Distal	Proximal	Distal	Proximal	
12G	12G	16	12	16720	18540	1.62	1.45	<a href="#">MXA232X16X12</a>
12G	12G	20	12	15870	16950	1.69	1.62	<a href="#">MXA232X20X12</a>
12G	12G	25	12	15440	16950	1.97	1.77	<a href="#">MXA232X25X12</a>
12G	12G	16	12	16440	16470	1.55	1.43	<a href="#">MXA232X16X12S</a>
12G	12G	20	12	16080	17760	1.68	1.57	<a href="#">MXA232X20X12S</a>
14G	18G	16	7	5080	2380	0.63	0.70	<a href="#">MXA232X16X70</a>
14G	18G	20	7	4710	1930	0.66	0.69	<a href="#">MXA232X20X70</a>
14G	18G	30	7	3950	1370	0.81	0.77	<a href="#">MXA232X30X70</a>
14G	14G	20	8.5	5160	3800	0.70	0.69	<a href="#">MXA232X20X85</a>
14G	14G	25	8.5	4730	3140	0.74	0.77	<a href="#">MXA232X25X85</a>
18G	22G	8	4	1020	400	0.36	0.43	<a href="#">MXA232X08X40</a>
18G	22G	10	4	850	330	0.36	0.43	<a href="#">MXA232X10X40</a>
18G	22G	12	4	790	260	0.36	0.43	<a href="#">MXA232X12X40</a>
18G	22G	16	4	615	190	0.37	0.43	<a href="#">MXA232X16X40</a>
18G	22G	20	4	530	160	0.38	0.44	<a href="#">MXA232X20X40</a>
18G	18G	8	5.5	2370	1500	0.31	0.38	<a href="#">MXA232X08X55</a>
18G	18G	10	5.5	2200	1400	0.32	0.39	<a href="#">MXA232X10X55</a>
18G	18G	16	5.5	1750	960	0.34	0.43	<a href="#">MXA232X16X55</a>
18G	18G	20	5.5	1570	740	0.38	0.43	<a href="#">MXA232X20X55</a>

3 LUMEN

Gauge (G)			Length (cm)	French (Fr)	Flow Rate (ml/hour)			Priming Volume (ml)			Product Codes
Distal	Proximal	Middle			Distal	Proximal	Middle	Distal	Proximal	Middle	
14G	18G	18G	20	8.5	5360	1310	1250	0.69	0.65	0.59	<a href="#">MXA233X20X85</a>
14G	18G	18G	25	8.5	5520	1370	1320	0.91	0.65	0.58	<a href="#">MXA233X25X85</a>
16G	18G	18G	16	7	3650	1590	1440	0.48	0.45	0.43	<a href="#">MXA233X16X70</a>
16G	18G	18G	20	7	3270	1350	1260	0.54	0.48	0.43	<a href="#">MXA233X20X70</a>
16G	18G	18G	30	7	2600	910	850	0.62	0.51	0.45	<a href="#">MXA233X30X70</a>
16G	12G	12G	16	12	4190	13110	12610	0.58	1.26	1.28	<a href="#">MXA233X16X12</a>
16G	12G	12G	20	12	3780	11360	10970	0.64	1.33	1.46	<a href="#">MXA233X20X12</a>
16G	12G	12G	16	12	4540	13600	12995	0.54	1.25	1.31	<a href="#">MXA233X16X12S</a>
16G	12G	12G	20	12	4180	13210	12680	0.63	1.34	1.41	<a href="#">MXA233X20X12S</a>
18G	22G	22G	10	5.5	1980	380	350	0.36	0.38	0.33	<a href="#">MXA233X10X55</a>
18G	22G	22G	16	5.5	1640	240	230	0.42	0.38	0.33	<a href="#">MXA233X16X55</a>
20G	24G	24G	8	4	1020	170	165	0.36	0.39	0.34	<a href="#">MXA233X08X40</a>
20G	24G	24G	10	4	980	145	120	0.37	0.40	0.35	<a href="#">MXA233X10X40</a>
20G	24G	24G	12	4	850	110	100	0.37	0.40	0.35	<a href="#">MXA233X12X40</a>
20G	24G	24G	16	4	730	85	80	0.39	0.37	0.36	<a href="#">MXA233X16X40</a>
20G	24G	24G	20	4	650	65	60	0.37	0.39	0.35	<a href="#">MXA233X20X40</a>

4 LUMEN

Gauge (G)				Length (cm)	French (Fr)	Flow Rate (ml/hour)				Priming Volume (ml)				Product Codes
Distal	Proximal	Middle	Middle Distal			Distal	Proximal	Middle	Middle Distal	Distal	Proximal	Middle	Middle Distal	
16G	18G	18G	18G	16	8.5	4030	1840	1680	1580	0.50	0.45	0.44	0.41	<a href="#">MXA234X16X85</a>
16G	18G	18G	18G	20	8.5	3590	1440	1300	1160	0.51	0.46	0.43	0.45	<a href="#">MXA234X20X85</a>
16G	18G	18G	18G	25	8.5	3080	1170	1100	1010	0.56	0.54	0.49	0.45	<a href="#">MXA234X25X85</a>
16G	18G	18G	18G	30	8.5	3130	920	830	800	0.58	0.57	0.54	0.55	<a href="#">MXA234X30X85</a>

5 LUMEN

Gauge (G)					Length (cm)	French (Fr)	Flow Rate (ml/hour)					Priming Volume (ml)					Product Codes
Distal	Proximal	Middle	Middle Distal	Middle Proximal			Distal	Proximal	Middle	Middle Distal	Middle Proximal	Distal	Proximal	Middle	Middle Distal	Middle Proximal	
16G	18G	18G	18G	18G	16	8.5	4030	2080	1640	1560	1800	0.45	0.46	0.41	0.38	0.45	<a href="#">MXA235X16X85</a>
16G	18G	18G	18G	18G	20	8.5	3390	1440	1180	1080	1280	0.50	0.47	0.48	0.50	0.49	<a href="#">MXA235X20X85</a>
16G	18G	18G	18G	18G	25	8.5	3010	1120	950	890	1010	0.59	0.48	0.48	0.50	0.47	<a href="#">MXA235X25X85</a>
16G	18G	18G	18G	18G	30	8.5	2770	960	830	800	880	0.7	0.56	0.55	0.58	0.50	<a href="#">MXA235X30X85</a>



Single lumen standard

Smith Medical's Implantable Ports are available in a wide range of types and sizes to meet the needs of every therapy and every patient size. From titanium to plastic and from silicone to polyurethane, we have a port for every need.

#### **PORT-A-CATH® Implantable Access Systems**

- The gold standard in ports.

#### **PORT-A-CATH® II Implantable Access Systems**


- The power of two materials in one.

#### **ProPort® Implantable Access Systems**

- Cost effective and high quality.

#### **P.A.S. PORT® Implantable Peripheral Access Systems**

- The first system designed for arm placement.



## With Smiths Medical ports your choice continues with our range of

- Ports specifically for venous, arterial, intraspinal, epidural, peritoneal and peripheral placement.
- Catheter materials in Silicone or PolyFlow® Polyurethane.
- Single or dual lumen.
- Titanium, polysulfone or plastic.
- CATH-SHIELD® for easy and reliable catheter connections.
- Preassembled and unassembled Kits and Trays.

Smiths Medical portal systems and components are latex free and PVC free. Portals are MRI compatible up to 1.5 tesla.<sup>1</sup>

PolyFlow® Polyurethane features a smaller outer diameter than silicone catheters but the same internal diameter. High tensile and shear strength designed to withstand high pressure. PolyFlow® is firm at room temperature for easier advancement into the vein with silicone like softness at body temperature.

## Kit Content:

Kits include portal, catheter, PORT-A-CATH® straight needle, blunt needle (except preassembled systems) and vein pick.

## Tray Content:

Trays include portal, catheter, GRIPPER PLUS® Safety Needle, vein pick, stainless steel tunneling tool, PORT-A-CATH® straight needle, blunt needle, 18G extra thin wall introducer needle, J tip guidewire, vessel dilator/sheath introducer assembly, 90 degree PORT-A-CATH® needle and 2 syringes and POINT-LOK®.

## References

1. Shedlock, FG.MR Imaging and Implantable Vascular Access Devices: A practical Review, Monograph, published by SIMS Deltec, Inc. 1994.

## Standard System Single Lumen

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Titanium	28mm	14.6mm	16g	11.4mm	0.5ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4024-24 Tray: 21-4025-24
Titanium	28mm	14.6mm	16g	11.4mm	0.5ml	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	1.5ml	Kit: 21-4022-24 Tray: 21-4023-24
Titanium	28mm	14.6mm	16g	11.4mm	0.5ml	Silicone	2.8mm	1.0mm	76cm	9F	0.6ml	Kit: 21-4008-24 Tray: 21-4009-24
Titanium	25.4mm	13.5mm	16g	11.4mm	0.5ml	Silicone	2.8mm	1.0mm	76cm	9F	0.6ml	Kit: 21-4000-24 Tray: 21-4003-24

## Standard System Single Lumen Preassembled

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Titanium	28mm	13.5mm	16g	11.4mm	1.1ml per system	Silicone	2.8mm	1.0mm	76cm	9F	1.1ml per system	Kit: 21-4010-24 Tray: 21-4004-24

## Standard System Single Lumen Low Profile™

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Titanium	25mm	11.5mm	9.5g	9.5mm	0.3ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4036-24 Tray: 21-4037-24
Titanium	25mm	11.5mm	9.5g	9.5mm	0.3ml	Silicone	2.8mm	1.0mm	76cm	9F	0.6ml	Kit: 21-4034-24 Tray: 21-4035-24

## Standard System Double Lumen

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Titanium	46.7mm long 26.5mm wide	14.4mm	34g	11.4mm each	0.5ml per reservoir	Silicone	3.4mm per lumen	1.1mm per lumen	76cm per lumen	11F	0.7ml per lumen	Kit: 21-8010-24 Tray: 21-8011-24



● Single lumen low profile



● Single lumen standard



● Dual lumen standard

## PORT-A-CATH®

The titanium construction provides long term durability and gouge resistance, withstanding repeated needle access. The system's titanium housing features a contoured shape for increased patient comfort and ease of portal palpation. The distinct rounded septum ring is designed to assist in septum location and the high compression SECURE SITE® septum allows greater needle retention and stability.

PORT-A-CATH® II Single Lumen

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Polysulfone Titanium	30.5mm	14.7mm	10g	11.4mm	0.5ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4052-24 Tray: 21-4053-24
Polysulfone Titanium	30.5mm	14.7mm	10g	11.4mm	0.5ml	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	1.5ml	Kit: 21-4054-24 Tray: 21-4055-24
Polysulfone Titanium	30.5mm	14.7mm	10g	11.4mm	0.5ml	Silicone	2.8mm	1.0mm	76cm	9F	0.6ml	Kit: 21-4050-24 Tray: 21-4051-24

PORT-A-CATH® II Single Lumen Low Profile™

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Polysulfone Titanium	25mm	11.5mm	4.8g	9.5mm	0.3ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4082-24 Tray: 21-4083-24
Polysulfone Titanium	25mm	11.5mm	4.8g	9.5mm	0.3ml	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	1.5ml	Kit: 21-4070-24 Tray: 21-4071-24

PORT-A-CATH® II Single Lumen Preassembled

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Polysulfone Titanium	30.5mm	14.7mm	10g	11.4mm	0.7ml or 1.1ml system	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.5ml 1.1ml system	Tray: 21-4063-24
Polysulfone Titanium	30.5mm	14.7mm	10g	11.4mm	0.5ml or 2.0ml system	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	1.5ml 2.0ml system	Tray: 21-4065-24



● Single lumen low profile



● Single lumen standard

Patient Size Guide

- Small Patients
- Medium Patients
- Large Patients

PORT-A-CATH® II

Discover the powerful combination of two materials, polysulfone and titanium, together in one portal.

All the reliability and quality of PORT-A-CATH® in a lighter weight portal design. By engineering the portal with polysulfone on the outside and titanium on the inside we have created an access system designed to meet clinical needs for lighter weight and long term durability.

## PORT-A-CATH® II Single Lumen Low Profile™ Preassembled

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Polysulfone Titanium	25mm	11.5mm	4.8g	9.5mm	0.3ml or 1.0ml system	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml or 1.0ml system	<a href="#">Kit: 21-4084-24</a> <a href="#">Tray: 21-4085-24</a>
Polysulfone Titanium	25mm	11.5mm	4.8g	9.5mm	0.3ml or 1.8ml system	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	0.7ml or 1.5ml system	<a href="#">Tray: 21-4073-24</a>

## PORT-A-CATH® II Double Lumen

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Polysulfone Titanium	50mm long x30mm wide	16mm	24g	11.4mm each septum	0.5ml per reservoir	Silicone	3.4mm per lumen	1.1mm per lumen	76cm	11F	0.7ml per lumen	<a href="#">Kit: 21-8050-24</a> <a href="#">Tray: 21-8052-24</a>

## PORT-A-CATH® II Double Lumen Preassembled

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Polysulfone Titanium	50mm long x30mm wide	16mm	24g	11.4mm each septum	0.5ml	Silicone per reservoir 2.4ml system	3.4mm	1.1mm per lumen	76cm per lumen	11F	0.7ml per lumen 2.4ml system	<a href="#">Tray: 21-8053-24</a>

## PORT-A-CATH® II Double Lumen Low Profile™

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Polysulfone Titanium	38.7mm long x 23.5mm wide	11mm	11g	9.5mm	0.3ml per lumen	PolyFlow® Polyurethane	2.2mm per lumen	1.0mm per lumen	76cm	7F	0.6ml per lumen	<a href="#">Kit: 21-8065-24</a> <a href="#">Tray: 21-8066-24</a>
Polysulfone Titanium	38.7mm long x 23.5mm wide	11mm	11g	9.5mm	0.3ml per lumen	PolyFlow® Polyurethane	3.2mm per lumen	1.4mm per lumen	76cm	10F	1.2ml per lumen	<a href="#">Kit: 21-8067-24</a> <a href="#">Tray: 21-8068-24</a>

PORT-A-CATH® II systems are also available for Fluoro-Free® placement with the CATH-FINDER® catheter tracking system



● ● Dual lumen low profile



● ● Dual lumen standard

### Patient Size Guide

- Small Patients
- Medium Patients
- Large Patients

Standard System Single Lumen

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic	30.5mm	14.0	9g	12.7mm	0.85ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4152-24 Tray: 21-4153-24
Plastic	30.5mm	14.0	9g	12.7mm	0.85ml	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	1.5ml	Kit: 21-4154-24 Tray: 21-4155-24

Standard System Single Lumen Preassembled

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic	30.5mm	14.0	9g	12.7mm	0.92ml system	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml 0.9ml system	Kit: 21-4162-24 Tray: 21-4163-24
Plastic	30.5mm	14.0	9g	12.7mm	2.35ml system	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	1.5ml 2.35ml system	Kit: 21-4164-24 Tray: 21-4165-24

Standard System Single Lumen Low Profile™

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic	25mm	10.0mm	4g	10.5mm	0.45ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4182-24 Tray: 21-4183-24
Plastic	25mm	10.0mm	4g	10.5mm	0.45ml	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	1.5ml	Kit: 21-4170-24 Tray: 21-4171-24

Standard System Single Lumen Low Profile™ Preassembled

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic	25mm	10.0mm	4g	10.5mm	0.45ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4184-24 Tray: 21-4185-24
Plastic	25mm	10.0mm	4g	10.5mm	0.45ml	PolyFlow® Polyurethane	2.6mm	1.6mm	76cm	8.5F	1.5ml	Kit: 21-4172-24 Tray: 21-4173-24

**ProPort®**

ProPort® is our range of lightweight plastic ports. Like the PORT-A-CATH®, PORT-A-CATH® II, P.A.S. PORT® T2 ranges ProPort® also has a SECURE SITE® septum. It also has a unique SlideLock™ Connector which is designed to provide easy and secure port to catheter locking. The low profile ProPort® option, with its contour shape provides a more discreet alternative helping to maximise comfort.



Standard System Single Lumen

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic	30.5mm	14.0mm	9g	12.7mm	0.85ml	Silicone	2.8mm	1.0mm	76cm	9F	0.7ml	Kit: 21-4150-24 Tray: 21-4151-24

Standard System Single Lumen Preassembled

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic	30.5mm	14.0mm	9g	12.7mm	0.85ml	Silicone	2.8mm	1.0mm	76cm	9F	0.7ml or 1.55ml system	Kit: 21-4160-24 Tray: 21-4161-24

Standard System Single Lumen Low Profile™

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic	25mm	10.0mm	4g	10.5mm	0.45ml	Silicone	2.8mm	1.0mm	76cm	9F	0.7ml	Kit: 21-4186-24 Tray: 21-4187-24

Standard System Single Lumen Low Profile™ Preassembled

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic	25mm	10.0mm	4g	10.5mm	0.45ml or 1.15ml system	Silicone	2.8mm	1.0mm	76cm	9F	0.7ml or 1.15ml system	Kit: 21-4188-24 Tray: 21-4189-24



● Single lumen low profile



●● Single lumen standard

Patient Size Guide

- Small Patients
- Medium Patients
- Large Patients

Standard System Single Lumen P.A.S. PORT® T2

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Titanium	24.5mm long x 18.2mm wide	11.5mm	8.3g	9.5mm	0.3ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4572-24 Tray: 21-4573-24

Standard System Single Lumen P.A.S.PORT® Elite

Material	Base	Height	Weight	Septum Volume	Priming Volume	Catheter Material	O.D.	I.D.	Length	Introducer Size	Priming Volume	Product Codes
Plastic Titanium	24.9mm long x 17.1mm wide	9.8mm	2.8g	9.7mm	0.3ml	PolyFlow® Polyurethane	1.9mm	1.0mm	76cm	6F	0.7ml	Kit: 21-4590-24 Tray: 21-4591-24

P.A.S.PORT® systems are also available for Fluoro-Free® placement with the CATH-FINDER catheter tracking system



● P.A.S.PORT® Elite



●● P.A.S.PORT® T2

Accessories

Ports Introducer Sets

Size	Qty	Product Codes
6F	5 per box	21-2320-24
8.5F	5 per box	21-2470-24
9F	5 per box	21-2340-24
11F	5 per box	21-2386-24

Sets include 18G thin wall introducer needle, J-tip guidewire, vessel dilator sheaths introducer assembly, 90 degree PORT-A-CATH® needle, and 2 syringes.

Ports Subcutaneous Tunneling Device

0.D. x 24cm length, solid malleable stainless steel, disposable

Size	Qty	Product Codes
3.2mm	5 per box	221-0125-24

P.A.S PORT®

For those of your patients requiring peripheral port placement, Smiths Medical has designed the P.A.S. PORT® range. P.A.S. PORT® T2 is a titanium port suitable for upper or lower arm placement. It also has a SECURE SITE® septum designed for needle retention and stability. A lighter weight, low profile option is P.A.S. PORT® Elite, which is made of an acetal (plastic) housing exterior and a titanium base. It has a unique WIDESITE® septum designed for a large needle puncture area.

# GRIPPER® AND GRIPPER PLUS®

## GRIPPER® Needles

Gauge	Length Inches	Length	With Y site	Without Y site	Qty
22G	1 ¼	32mm	21-2941-24	21-2715-24	12 per box
22G	1	25mm	21-2940-24	21-2714-24	12 per box
22G	¾	19mm	21-2939-24	21-2733-24	12 per box
22G	½	16mm	21-2938-24	21-2737-24	12 per box
20G	1 ¼	32mm	21-2949-24	21-2718-24	12 per box
20G	1	25mm	21-2948-24	21-2717-24	12 per box
20G	¾	19mm	21-2947-24	21-2734-24	12 per box
20G	½	16mm	21-2946-24	21-2736-24	12 per box
19G	1 ¼	32mm	21-2957-24	21-2721-24	12 per box
19G	1	25mm	21-2956-24	21-2720-24	12 per box
19G	¾	19mm	21-2955-24	21-2735-24	12 per box

## GRIPPER PLUS® Safety Needles

Gauge	Length Inches	Length	With Y site	Without Y site	Qty
22G	1 ¼	32mm	21-2863-24	21-2761-24	12 per box
22G	1	25mm	21-2862-24	21-2760-24	12 per box
22G	¾	19mm	21-2861-24	21-2766-24	12 per box
22G	½	16mm	21-2860-24	21-2770-24	12 per box
20G	1 ¼	32mm	21-2867-24	21-2763-24	12 per box
20G	1	25mm	21-2866-24	21-2762-24	12 per box
20G	¾	19mm	21-2865-24	21-2767-24	12 per box
20G	½	16mm	21-2864-24	21-2769-24	12 per box
19G	1 ¼	32mm	21-2870-24	21-2765-24	12 per box
19G	1	25mm	21-2869-24	21-2764-24	12 per box
19G	¾	19mm	21-2868-24	21-2768-24	12 per box



GRIPPER PLUS®

## GRIPPER®

GRIPPER® needles designed to provide more confident needle placement for implanted access systems.

### GRIPPER® provides

- Natural rubber latex free to reduce potential allergic reactions.
- Colour coded clamps for needle gauge identification.
- Removable contoured grip for controlled needle placement.
- Cushioned needle platform for greater patient comfort, needle stabilisation and access site protection.
- Variety of gauges and lengths to accommodate most portals.
- Available with or without y-site.

## GRIPPER PLUS® (Safety)

Effective needlestick protection is only part of the GRIPPER PLUS® safety needle story, our exciting new needle was created with your input, and the new GRIPPER PLUS® Safety Needle delivers on all counts.

### GRIPPER PLUS® provides

- Audible and visible confirmation of safety mechanism.
- Minimal change to current technique.
- Kink- resistant TOTM tubing available with or without y-site.
- Colour coded clamp for ease of needle gauge identification.
- Fully threaded Luer lock connectors fits all standard accessories.
- Removable slit septum cap offers needleless access flexibility.
- Cushioned needle platform for greater patient comfort.
- Non-coring, non siliconised needle for secure needle retention.

## Plastic Huber Needles straight

Gauge	Length Inches	Length	With Y site	Qty
24G	1	25mm	21-2003-24	12 per box
22G	1½	38mm	21-2007-24	12 per box
22G	1	25mm	21-2005-24	12 per box
20G	1½	38mm	21-2015-24	12 per box
20G	1	25mm	21-2013-24	12 per box
19G	1	25mm	21-2021-24	12 per box

## Plastic Huber Needles 90° Bend

Gauge	Length Inches	Length	With Y site	Qty
22G	1 1/2	38mm	21-2008-24	12 per box
22G	1	25mm	21-2006-24	12 per box
22G	1/2	13mm	21-2316-24	12 per box
20G	1 1/2	38mm	21-2016-24	12 per box
20G	1	25mm	21-2014-24	12 per box
20G	3/4	19mm	21-2287-24	12 per box
20G	1/2	13mm	21-2317-24	12 per box



## Traditional Huber Needles

- Available with plastic hub.
- Straight or 90 degree bent.
- Non coring, non siliconised needles.