


# Formula Cutters & Burs



**Resector**   **Aggressive Plus**   **Tomcat**   **Scalloped**   **Cougar End**   **End Cutter**   **True End**   **Double Bite**   **Slotted Whisker**   **Whisker**

## Formula Series Cutters

Stryker's new advances in cutter technology combine increased sharpness, enhanced suction, and the tightest tolerances for optimal cutting blades.

<b>Resector Cutters</b> (Quantity: box of 5)	375-532-000	3.5 mm	<b>Scalloped Cutter</b> (Quantity: box of 5)	375-546-000	4.0 mm	
	375-542-000	4.0 mm		375-566-000	5.5 mm	
	375-552-000	5.0 mm		<b>Cougar End Cutters</b> (Quantity: box of 5)	375-531-000	3.5 mm
	375-562-000	5.5 mm			375-541-000	4.0 mm
	380-542-150	4.5 mm Angled			375-551-000	5.0 mm
	380-552-100	5.0 mm Angled			<b>End Cutters</b> (Quantity: box of 5)	375-737-000
<b>Aggressive Plus Cutters</b> (Quantity: box of 5)	375-534-000	3.5 mm	375-747-000	4.0 mm		
	375-544-000	4.0 mm	375-757-000	5.0 mm		
	375-554-000	5.0 mm	<b>True End Cutter</b> (Quantity: box of 5)	375-884-000	4.0 mm	
375-564-000	5.5 mm	<b>Double Bite Cutter</b> (Quantity: box of 5)		375-543-000	4.0 mm	
380-544-100	4.0 mm Angled		380-543-150	4.5 mm Angled		
<b>Tomcat Cutters</b> (Quantity: box of 5)	375-535-000	3.5 mm	<b>Slotted Whisker Cutters</b> (Quantity: box of 5)	375-548-000	4.0 mm	
	375-545-000	4.0 mm		375-558-000	5.0 mm	
	375-555-000	5.0 mm	<b>Whisker Cutter</b> (Quantity: box of 5)	375-549-000	4.0 mm	
	375-565-000	5.5 mm				
	380-545-100	4.0 mm Angled				
	380-545-150	4.5 mm Angled				
380-555-100	5.0 mm Angled					

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.



**12 Flute Barrel Bur**



**6 Flute Barrel Bur**



**Unhooded Barrel Bur**



**12 Flute Round Bur**



**6 Flute Round Bur**



**Unhooded Round Bur**



**Egg Bur**

## Formula Series Burs

<b>Standard 12 Flute Barrel Burs</b> (Quantity: box of 5)	<b>375-941-012</b> 4.0 mm	<b>Standard 12 Flute Round Burs</b> (Quantity: box of 5)	<b>375-940-012</b> 4.0 mm
	<b>375-951-112</b> 5.0 mm		<b>375-950-112</b> 5.0 mm
	<b>375-951-012</b> 5.5 mm		<b>375-950-012</b> 5.5 mm
<b>Aggressive 6 Flute Barrel Burs</b> (Quantity: box of 5)	<b>375-941-000</b> 4.0 mm	<b>Aggressive 6 Flute Round Burs</b> (Quantity: box of 5)	<b>375-930-000</b> 3.5 mm
	<b>375-951-100</b> 5.0 mm		<b>375-940-000</b> 4.0 mm
	<b>375-951-000</b> 5.5 mm		<b>375-950-100</b> 5.0 mm
			<b>375-950-000</b> 5.5 mm
<b>Unhooded 6 Flute Barrel Burs</b> (Quantity: box of 5)	<b>375-941-200</b> 4.0 mm	<b>Unhooded 6 Flute Round Burs</b> (Quantity: box of 5)	<b>375-940-200</b> 4.0 mm
	<b>375-951-200</b> 5.5 mm		<b>375-950-200</b> 5.5 mm
		<b>Egg Bur 6 Flute</b> (Quantity: box of 5)	<b>375-952-000</b> 5.5 mm

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.

### Specialty Burs



6 Flute SLAP Bur



Auger Bur



Subchondral Drill



6 Flute Left Helix Barrel Bur



12 Flute Left Helix Barrel Bur

#### Formula Series Specialty Burs

(Quantity: box of 5)

Protected bur design allows safe access to the glenoid without compromising the condition of the labrum

375-941-500 4.0 mm 6 Flute SLAP Bur

Unique cork screw design enables the removal of both bone and soft tissue

375-450-500 5.0 mm Auger Bur

For the initial treatment of articular cartilage defects

375-832-000 4.0 mm Subchondral Drill (1.5 mm Drill Tip)

Eliminates the bur head from skipping off of the bone and enables surgeons to use the same technique on right and left joint

375-953-000 5.5 mm 6 Flute Left Helix Barrel Bur  
375-953-012 5.5 mm 12 Flute Left Helix Barrel Bur

### Angled Cutters



Resector



Aggressive Plus



Tomcat



Double Bite



Proprietary "puzzle piece" Laser Cut Design shaver inner-housing

#### Formula Series Angled Cutters

(Quantity: box of 5)

380-542-150 4.5 mm Angled Resector Cutter

380-552-100 5.0 mm Angled Resector Cutter

380-544-100 4.0 mm Angled Aggressive Plus Cutter

380-544-150 4.5 mm Angled Aggressive Plus Cutter

380-554-100 5.0 mm Angled Aggressive Plus Cutter

380-545-100 4.0 mm Angled Tomcat Cutter

380-545-150 4.5 mm Angled Tomcat Cutter

380-555-100 5.0 mm Angled Tomcat Cutter

380-543-150 4.5 mm Angled Double Bite Cutter

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.