

Urological Products PRODUCT CATALOG



Care. Comfort. Choice.™

BAIRD | **CARE**
has joined BD







TABLE OF CONTENTS

INTERMITTENT CATHETERS	5
SILICONE INTERMITTENT CATHETERS	
Magic ³ Go [®] Female Intermittent Catheter	6
Magic ³ Go [®] Pediatric Intermittent Catheter	7
Magic ³ Go [®] Male Intermittent Catheter	8
Magic ^{3®} Intermittent Catheter	9
Magic ^{3®} Hydrophilic Intermittent Catheter	10
Personal Catheter [®] Intermittent Catheter	12
Personal Catheter [®] Hydrophilic Intermittent Catheter	13
INTERMITTENT CATHETERS	
Clean-Cath [®] Intermittent Catheter	14
Red Rubber Intermittent Catheter	16
CLOSED SYSTEMS AND URETHRAL TRAYS	19
Touchless [®] Plus Intermittent Catheter System	20
Bardia [®] Urethral Trays	22
INCONTINENCE SOLUTIONS	25
PureWick [™] Urine Collection System	26
PureWick [™] Female External Catheter	27
MALE EXTERNAL CATHETERS	
Spirit [®] Male External Catheter Style 1	28
Spirit [®] Male External Catheter Style 2	28
Spirit [®] Male External Catheter Style 3	28
UltraFlex [®] Male External Catheter	29
Pop-On [®] Male External Catheter	29
WideBand [®] Male External Catheter	30
Natural [™] Male External Catheter	30
LEG BAGS, TUBING AND ACCESSORIES	
Leg Bags and Leg Bag Systems	31
Straps	32
Tubing	32
Drainage Bags	33
FOLEY CATHETERS AND FOLEY TRAYS	35
2-Way Latex Foley Catheters	36
2-Way All-Silicone Foley Catheters	37
Foley Catheter Trays	38
Closed System Foley Trays	39
Piston Syringe	39
Index & HCPCS Codes	40





Intermittent Catheters

Summary of Indications and Warnings

INTERMITTENT CATHETERS

Intermittent catheters are intended for draining urine from the bladder. Reuse of a single use device may create a risk of patient or user infection. Please consult product labels and inserts for more detailed safety information and instructions for use.

Please consult product label and insert for any indications, contraindications, hazards, warnings, cautions and directions for use.

Magic³ Go[®] Female


HYDROPHILIC INTERMITTENT CATHETER

with CLEANGLIDE™ Technology

- Hydrophilic, female length catheter designed to glide easily and comfortably through the urinary tract
- New, unique CLEANGLIDE™ technology that maintains the catheter in a hydrated state so that it is ready to go out of the package
- Three-layered silicone catheter with an ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation
- Comfort-designed tapered tip with four comfort-sized drainage eyes
- Robust handle for better control that is ergonomically designed for a secure grip that can easily be connected to a collection device
- Designed to minimize the risk of contamination by enabling the user to hold the catheter without touching the catheter shaft
- Discreet and compact consumer package that is easy to carry, easy to open, and easy to use right out of the package without the need of adding water or lubricant
- Environmentally conscious design - does not contain DEHP and/or phthalates



MAGIC³ GO™ catheter with three layers of silicone



CLEANGLIDE™ Technology

Hydrophilic coating

Reorder Code	Length	French Size	Qty/Case
51808	6"	8 Fr.	30/Case
51810	6"	10 Fr.	
51812	6"	12 Fr.	
51814	6"	14 Fr.	
51816	6"	16 Fr.	
51818	6"	18 Fr.	

Magic³ Go[®] Pediatric

HYDROPHILIC INTERMITTENT CATHETER

with CLEANGLIDE[™] Technology

- Hydrophilic, pediatric length catheter designed to glide easily and comfortably through the urinary tract
- New, unique CLEANGLIDE[™] technology that maintains the catheter in a hydrated state so that it is ready to go out of the package.
- Three-layered silicone catheter with an ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation
- Comfort-designed tapered tip with four comfort-sized drainage eyes.
- Robust SureGrip[™] insertion sleeve for better control that is ergonomically designed for a secure grip that can easily be connected to a collection device.
- Designed to minimize the risk of contamination by enabling the user to hold the catheter without touching the catheter shaft
- Discreet and compact consumer package that is easy to carry, easy to open, and easy to use right out of the package without the need of adding water or lubricant.
- Environmentally conscious design - does not contain DEHP and/or phthalates

Reorder Code	Length	French Size	Tip Style	Qty/Case
52806G	10"	6 Fr.	Straight	30/Case
52808G	10"	8 Fr.	Straight	
52810G	10"	10 Fr.	Straight	



100% SILICONE INTERMITTENT CATHETERS

Magic³ Go[®] Male

HYDROPHILIC INTERMITTENT CATHETER

with CLEANGLIDE™ Technology

- Hydrophilic, male length catheter designed to glide easily and comfortably through the urinary tract
- New, unique CLEANGLIDE™ technology that maintains the catheter in a hydrated state so that it is ready to go out of the package.
- Three-layered silicone catheter with an ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation
- Comfort-designed tapered tip with four comfort-sized drainage eyes.
- Robust SureGrip™ insertion sleeve for better control.
- Designed to minimize the risk of contamination by enabling the user to hold the catheter without touching the catheter shaft
- Discreet and compact consumer package that is easy to carry, easy to open, and easy to use right out of the package without the need of adding water or lubricant.
- Environmentally conscious design - does not contain DEHP and/or phthalates



Reorder Code	Length	French Size	Tip Style	Qty/Case
53810G	16"	10 Fr.	Straight	30/Case
53812G	16"	12 Fr.	Straight	
53814G	16"	14 Fr.	Straight	
53816G	16"	16 Fr.	Straight	
53818G	16"	18 Fr.	Straight	
53820G	16"	20 Fr.	Straight	
50810G	16"	10 Fr.	Coudé	
50812G	16"	12 Fr.	Coudé	
50814G	16"	14 Fr.	Coudé	
50816G	16"	16 Fr.	Coudé	
50818G	16"	18 Fr.	Coudé	




Magic³[®]

INTERMITTENT CATHETER

100% SILICONE INTERMITTENT CATHETERS

- Designed to be gentle on delicate urethral tissue¹
- Ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation
- Comfort-designed tapered tip with four comfort-sized drainage eyes
- Not made with natural rubber latex
- Phthalate-free, eliminating the risk of health concerns with phthalate plasticizers
- PVC-free, eliminating concerns on the disposal of intermittent catheters



Reorder Code	Length	French Size	Qty/Case
51310	6"	10 Fr.	30/Case
51312	6"	12 Fr.	
51314	6"	14 Fr.	
51316	6"	16 Fr.	
51318	6"	18 Fr.	
51320	6"	20 Fr.	
52306	10"	6 Fr.	
52308	10"	8 Fr.	
52310	10"	10 Fr.	
53310	16"	10 Fr.	
53312	16"	12 Fr.	
53314	16"	14 Fr.	
53316	16"	16 Fr.	
53318	16"	18 Fr.	
53320	16"	20 Fr.	

1. Data on file. Preclinical testing may not correlate with outcomes in humans

Magic³[®]

HYDROPHILIC INTERMITTENT CATHETER

- Packaged with its own sterile water packet, immediate hydration creates virtually friction-free surface¹
- Ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation¹
- Available with SURE-GRIP™ insertion sleeve, designed to ease handling and eliminate the need to touch the catheter surface²
- Comfort-designed tapered tip with four comfort-sized drainage eyes¹
- Not made with natural rubber latex
- Phthalate-free, eliminating the risk of health concerns with phthalate plasticizers¹
- PVC-free, eliminating concerns on the disposal of intermittent catheters¹



Reorder Code	Length	French Size	SureGrip™	Qty/Case
51610	6"	10 Fr.	No	30/Case
51612	6"	12 Fr.	No	
51614	6"	14 Fr.	No	
51616	6"	16 Fr.	No	
51618	6"	18 Fr.	No	
51620	6"	20 Fr.	No	
52606	10"	6 Fr.	No	
52608	10"	8 Fr.	No	
52610	10"	10 Fr.	No	
53610	16"	10 Fr.	No	
53612	16"	12 Fr.	No	
53614	16"	14 Fr.	No	
53616	16"	16 Fr.	No	
53618	16"	18 Fr.	No	
53620	16"	20 Fr.	No	

1. Data on file. Preclinical testing may not correlate with outcomes in humans

Magic³[®]

HYDROPHILIC INTERMITTENT CATHETER



100% SILICONE INTERMITTENT CATHETERS

With SUREGRIP™ Insertion Sleeve + Coudé Tip

Reorder Code	Length	French Size	SureGrip™	Qty/Case
50610	16"	10 Fr.	Yes	30/Case
50612	16"	12 Fr.	Yes	
50614	16"	14 Fr.	Yes	
50616	16"	16 Fr.	Yes	
50618	16"	18 Fr.	Yes	

With SUREGRIP™ Insertion Sleeve

Reorder Code	Length	French Size	SureGrip™	Qty/Case
52606G	10"	6 Fr.	Yes	30/Case
52608G	10"	8 Fr.	Yes	
52610G	10"	10 Fr.	Yes	
53610G	16"	10 Fr.	Yes	
53612G	16"	12 Fr.	Yes	
53614G	16"	14 Fr.	Yes	
53616G	16"	16 Fr.	Yes	
53618G	16"	18 Fr.	Yes	
53620G	16"	20 Fr.	Yes	

With Insertion Supply Kit, which includes Nitrile gloves, BZK wipes, underpad, and collection bag. 536XXGS codes include SUREGRIP™ insertion sleeve

Reorder Code	Length	French Size	SureGrip™	Qty/Case
51614S	6"	14 Fr.	Yes	50/Case
53614GS	16"	14 Fr.	Yes	
53616GS	16"	16 Fr.	Yes	



Personal Catheter®

SILICONE INTERMITTENT CATHETER

The 100% silicone intermittent catheter features a soft silicone construction designed to provide comfortable, easy advancement, while minimizing the potential for trauma¹

Reorder Code	Length	French Size	Qty/Case
61310	6"	10 Fr.	30/Case
61312	6"	12 Fr.	
61314	6"	14 Fr.	
61316	6"	16 Fr.	
61318	6"	18 Fr.	
62308	10"	8 Fr.	
62310	10"	10 Fr.	
63310	16"	10 Fr.	
63312	16"	12 Fr.	
63314	16"	14 Fr.	
63316	16"	16 Fr.	
63318	16"	18 Fr.	



1. Data on file. Preclinical testing may not correlate with outcomes in humans

Personal Catheter®

HYDROPHILIC SILICONE INTERMITTENT CATHETER

100% SILICONE INTERMITTENT CATHETERS

Packaged with its own sterile water packet. Immediate hydration creates virtually friction-free surface¹

Reorder Code	Length	French Size	Qty/Case
61610	6"	10 Fr.	30/Case
61612	6"	12 Fr.	
61614	6"	14 Fr.	
61616	6"	16 Fr.	
61618	6"	18 Fr.	
62608	10"	8 Fr.	
62610	10"	10 Fr.	
63610	16"	10 Fr.	
63612	16"	12 Fr.	
63614	16"	14 Fr.	
63616	16"	16 Fr.	
63618	16"	18 Fr.	




1. Data on file. Preclinical testing may not correlate with outcomes in humans

Clean-Cath[®]

PVC INTERMITTENT CATHETER

- Finger grips for ease of use
- Flexible funnel
- Not made with natural rubber latex
- Individually packaged in easy to open peel-strip package
- Smooth round tip for easy insertion

Reorder Code	Length	French Size	Qty/Case
420706	6"	6 Fr.	50/Case
420708	6"	8 Fr.	
420710	6"	10 Fr.	
420712	6"	12 Fr.	
420714	6"	14 Fr.	
420716	6"	16 Fr.	
420718	6"	18 Fr.	
421108	10"	8 Fr.	
421110	10"	10 Fr.	
421708	16"	8 Fr.	
421710	16"	10 Fr.	
421712	16"	12 Fr.	
421714	16"	14 Fr.	
421716	16"	16 Fr.	
421718	16"	18 Fr.	



Clean-Cath[®]

PVC INTERMITTENT CATHETER

- Olive coudé tip helpful for patients with insertion difficulties

With Coude Tip

Reorder Code	Length	French Size	Qty/Case
423710	16"	10 Fr.	30/Case
423712	16"	12 Fr.	
423714	16"	14 Fr.	
423716	16"	16 Fr.	
423718	16"	18 Fr.	



Straight Tip Without Funnel

Reorder Code	Length	French Size	Qty/Case
420614	6"	14 Fr.	50/Case
421008	10"	8 Fr.	
421010	10"	10 Fr.	
421610	16"	10 Fr.	
421612	16"	12 Fr.	
421614	16"	14 Fr.	

Bardia®

RED RUBBER URETHRAL CATHETERS

- Radiopaque red rubber catheter with a round, hollow tip
- Can be used as a Robinson or a Nelaton catheter
- Two opposing drainage eyes

Red Rubber All-Purpose Urethral Catheters

Reorder Code	Length	French Size	Qty/Case
802408	16"	8 Fr.	100/Case
802410	16"	10 Fr.	
802412	16"	12 Fr.	
802414	16"	14 Fr.	
802416	16"	16 Fr.	
802418	16"	18 Fr.	
802420	16"	20 Fr.	
802422	16"	22 Fr.	



- Radiopaque red rubber catheter with a round, hollow tip
- Single drainage eye

Red Rubber Coudé Olive Tip Urethral Catheters

Reorder Code	Length	French Size	Qty/Case
802512	16"	12 Fr.	100/Case
802514	16"	14 Fr.	
802516	16"	16 Fr.	
802518	16"	18 Fr.	
802520	16"	20 Fr.	
802522	16"	22 Fr.	



CAUTION: This product contains natural rubber latex which may cause allergic reactions.





Closed Systems and Urethral Trays

Summary of Indications and Warnings

INTERMITTENT CATHETERS

Intermittent catheters are intended for draining urine from the bladder. Reuse of a single use device may create a risk of patient or user infection. Please consult product labels and inserts for more detailed safety information and instructions for use.

Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Please consult product label and insert for any indications, contraindications, hazards, warnings, cautions and directions for use.

Touchless® Plus

INTERMITTENT CATHETER SYSTEM

- Pre-lubricated vinyl catheter with finger guide
- 1100cc collection bag with built-in sample/drainage port
- The catheter cap has pull ring for easy removal - especially for people with limited dexterity
- Accessory codes include 2 gloves not made with natural rubber latex, povidone iodine swabs, and underpad

Vinyl

Reorder Code	French Size	Tip	Accessories	Qty/Case
4A5108	8 Fr.	Straight	Yes	50/case
4A5110	10 Fr.	Straight	Yes	
4A5142	12 Fr.	Straight	Yes	
4A5144	14 Fr.	Straight	Yes	
4A5146	16 Fr.	Straight	Yes	
4A6110	10 Fr.	Straight	No	
4A6142	12 Fr.	Straight	No	
4A6144	14 Fr.	Straight	No	
4A6146	16 Fr.	Straight	No	
4A7114	14 Fr.	Coudé	Yes	
4A7116	16 Fr.	Coudé	Yes	
4A8114	14 Fr.	Coudé	No	
4A8116	16 Fr.	Coudé	No	



Touchless® Plus

INTERMITTENT CATHETER SYSTEM

- Pre-lubricated red rubber catheter with finger guide
- 550cc or 1100cc collection bag with built-in sample/drainage port
- The catheter cap has pull ring for easy removal - especially for people with limited dexterity
- Accessory codes include 2 gloves not made with natural rubber latex, povidone iodine swabs, and underpad

Red Rubber

Reorder Code	French Size	Sex	Bag Size	Accessories	Qty/Case
4A3054	12 Fr.	Female	550cc	Yes	50/Case
4A3055	14 Fr.	Female	550cc	Yes	
4A2056	14 Fr.	Female	1100cc	Yes	
4A3052	12 Fr.	Male	550cc	No	
4A3053	14 Fr.	Male	550cc	No	
4A2044	14 Fr.	Male	1100cc	Yes	
4A2045	14 Fr.	Male	1100cc	Yes	
4A2046	12 Fr.	Male	1100cc	Yes	
4A5042	12 Fr.	Unisex	1100cc	Yes	
4A5044	14 Fr.	Unisex	1100cc	Yes	
4A6044	14 Fr.	Unisex	1100cc	No	
4A7044 (coudé)	14 Fr.	Unisex	1100cc	Yes	



CAUTION: The red rubber product codes contain natural rubber latex which may cause allergic reactions.


Bardia®


URETHRAL TRAYS

- Waterproof underpad, fenestrated drape, 3 skin cleansing swab sticks, 5g lubricating jelly, 2 gloves not made with natural rubber latex, specimen container with label
- Peel-top tray

Reorder Code	Catheter	Qty/Case
802115	With 15 Fr. red rubber urethral catheter	20/Case
802214	With 14 Fr. plastic urethral catheter	
802100	Without catheter	

CAUTION: The red rubber product contains natural rubber latex which may cause allergic reactions.





Incontinence Solutions

Summary of Indications and Warnings

PUREWICK™ URINE COLLECTION SYSTEM

The PureWick™ Urine Collection System is to be used with PureWick™ External Catheters which are intended for non-invasive urine output management. Do not use on individuals with urinary retention. Discontinue use if an allergic reaction occurs. Not recommended for users who are experiencing skin irritation or skin breakdown in device contact areas.

The PureWick™ Urine Collection System and accessories may be a potential biohazard. Clean and disinfect the PureWick™ Urine Collection System and accessories before disposal. Handle in accordance with accepted medical practice and applicable local and federal regulations. The PureWick™ Urine Collection System (PW100 and PW200) contains electrical and/or electronic equipment, and the PW200 product contains a lithium-ion battery that must be disposed of per local laws and regulations.

Please consult the PureWick™ Urine Collection System user guide for more detailed safety information and instructions for use.

PUREWICK™ FEMALE EXTERNAL CATHETER

The PureWick™ Female External Catheters is intended for non-invasive urine output management in female patients. To avoid potential skin injury, never push or pull the PureWick™ Female External Catheter against the skin during placement or removal. Never insert into the vagina, anal canal, or other body cavities. Discontinue use if an allergic reaction occurs. Not recommended for patients who are: agitated, combative, or uncooperative and might remove the PureWick™ Female External Catheter.

After use, this product may be a potential biohazard. Dispose of in accordance with applicable local, state and federal laws and regulations.

Please consult the PureWick™ Female External Catheter insert for more detailed safety information and instructions for use.

MALE EXTERNAL CATHETERS

Male external catheters are intended for use by male patients for the management of male urinary incontinence. Reuse of a single use device may create a risk of patient or user infection. Use of a single device for longer periods than 24 hours may increase the risk of complications. Do not use on irritated or compromised skin. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Please consult product label and insert for any indications, contraindications, hazards, warnings, cautions and directions for use.

PureWick™

URINE COLLECTION SYSTEM

- The PureWick™ Urine Collection System, the system to be used with the PureWick™ Female External Catheter, is an innovative device that is intended for non-invasive urine output
- Designed to draw urine away, helping to keep skin dry^{1,2}



*PureWick™ Urine
Collection System*

*PureWick™ Female
External Catheter*

PureWick™ System

Order#	Description	Qty.
PW100	PureWick™ Urine Collection System (<i>non-battery</i>)	1
PW200	PureWick™ Urine Collection System (<i>lithium-ion battery</i>)	1
PWKIT03	Accessory kit (<i>Replacement pump tubing, collector tubing and 2000cc collection canister</i>)	1

1. Data on file.
2. Ritter, J: Effectiveness and Functionality of a Female External Catheter. 2017

PureWick™

FEMALE EXTERNAL CATHETER

- Non-invasive way to manage urinary incontinence
- May reduce the risk of nighttime falls, by staying in bed¹
- Can make better sleep possible¹
- Can be placed in a few steps

Order#	Description	Qty.
PWFX30	PureWick™ Female External Catheter <i>(Not made with natural rubber latex).</i>	30

Attaches to the PureWick™ Urine Collection System or existing hospital wall suction

Aligns with pubic bone

Soft material wicks urine down to the base

Continuous low pressure suction at the base removes urine



¹ Based on survey responses (01/19) by 141 caregivers and users who purchased the PureWick™ System and a subsequent supply of PureWick™ Female External Catheters.

Spirit®

HYDROCOLLOID ADHESIVE SHEATH

- First and only self-adhering male sheath to combine the excellent moisture-wicking properties of a hydrocolloid adhesive with the exceptional breathability of silicone^{1,2}
- A broad line of male external catheters with a diverse range of coverage
- Hydrocolloid adhesive wicks away moisture to help protect skin integrity, minimize skin maceration and maintain adhesion^{1,2}
- Made with 100% silicone




Style 1



Offers built-in adhesive in a length and placement found with traditional male external catheters.

Size		30/Box	100/Box
Small	25mm	35301	35101
Medium	29mm	35302	35102
Intermediate	32mm	35303	35103
Large	36mm	35304	35104
X-Large	41mm	35305	35105

Style 2



The 1.5" sheath with forward placement of adhesive designed to minimize unintended detachment.

Size		30/Box	100/Box
Small	25mm	37301	37101
Medium	29mm	37302	37102
Intermediate	32mm	37303	37103
Large	36mm	37304	37104
X-Large	41mm	37305	37105

Style 3



Our most popular style provides additional adhesive area compared to similar silicone self-adhering catheters. Combined with forward placement of adhesive designed to minimize unintended detachment.²

Size		30/Box	100/Box
Small	25mm	39301	39101
Medium	29mm	39302	39102
Intermediate	32mm	39303	39103
Large	36mm	39304	39104
X-Large	41mm	39305	39105

1. Data on file. Preclinical testing may not correlate with outcomes in humans
 2. Voegali, David, British Journal of Community Nursing, Vol. 18, No. 1, Jan. 2013

UltraFlex®

SELF-ADHERING MALE EXTERNAL CATHETER

MALE EXTERNAL CATHETERS




- Offers built-in adhesive in a length and placement found with traditional male external catheters
- Made with 100% silicone

Size		30/Box	100/Box
Small	25mm	33301	33101
Medium	29mm	33302	33102
Intermediate	32mm	33303	33103
Large	36mm	33304	33104
X-Large	41mm	33305	33105



Pop-On®

SELF-ADHERING MALE EXTERNAL CATHETER



- The 1.5" sheath with forward placement of adhesive designed to minimize unintended detachment
- Made with 100% silicone

Size		30/Box	100/Box
Small	25mm	32301	32101
Medium	29mm	32302	32102
Intermediate	32mm	32303	32103
Large	36mm	32304	32104
X-Large	41mm	32305	32105



WideBand®

SELF-ADHERING MALE EXTERNAL CATHETER



- Provides additional adhesive area compared to similar silicone self-adhering catheters. Combined with forward placement of adhesive designed to minimize unintended detachment
- Made with 100% silicone

Size		30/Box	100/Box
Small	25mm	36301	36101
Medium	29mm	36302	36102
Intermediate	32mm	36303	36103
Large	36mm	36304	36104
X-Large	41mm	36305	36105



Natural™

NON-ADHESIVE, ALL-SILICONE MALE EXTERNAL CATHETER



- Non-adhesive male external catheter. Each box comes with three reusable foam straps for securement
- Made with 100% silicone

Size		30 Catheters and 3 Straps/Box
Small	25mm	38301
Medium	29mm	38302
Intermediate	32mm	38303
Large	36mm	38304
X-Large	41mm	38305



Dispoz-A-Bag®

LEG BAGS

- Made of heavy duty vinyl for long wear time
- FLIP-FLO™ drainage valve opens and closes with a simple flip
- Reusable, sterile



Leg Bag 4-Pack

Reorder Code	Description	Qty/Case
150419, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (Bonus pack 4 leg bags, 1 pair fabric leg straps)	12/Case

Leg Bag Systems

Reorder Code	Description	Qty/Case
150819, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (1 pair fabric leg straps, 18" PVC extension tubing)	50/Case
150919, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (18" PVC extension tubing, 4 1/2" drainage tube, no straps)	
151919, 32*	19, 32 oz. Leg Bag with FLIP-FLO™ valve (18" PVC extension tubing, 4 1/2" drainage tube, 1 pair latex straps)	
650819, 32*	19, 32 oz. Leg Bag with FLIP-FLO™ valve (1 pair fabric leg straps, 18" latex extension tubing)	
651919, 32*	19, 32 oz. Leg Bag with FLIP-FLO™ valve (18" latex extension tubing, 4 1/2" drainage tube, no straps)	

Leg Bag with Fabric Straps

Reorder Code	Description	Qty/Case
150719, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (1 pair fabric straps)	50/Case

***CAUTION:** Components of this product contain natural rubber latex which may cause allergic reactions.

Leg Bag with Latex Straps

Reorder Code	Description	Qty/Case
150106, 107, 108*	9, 19, 32 oz. Leg Bag with rubber cap valve (1 pair latex straps)	48/Case
150319, 32*	19, 32 oz. Leg Bag with FLIP-FLO™ valve (1 pair latex straps)	50/Case
4A4173, 74*	19, 32 oz. Leg Bag with drain tube and clamp (1 pair latex straps)	24/Case

Leg Bag with No Straps

Reorder Code	Description	Qty/Case
150619, 32	19, 32 oz. Leg Bag with FLIP-FLO™ valve (no straps)	50/Case

Straps

Reorder Code	Description	Qty/Case
150507	Deluxe fabric leg straps - 24" long, 3/4" wide	24 pair/Case
162110 (Small) 162210 (Medium) 162310 (Large)	2" wide leg straps with fabric backing and Velcro™ closure (9"-13" long, Small) (13"-20" long, Medium) (20"-27" long, Large)	6 pair/Case

Tubing

Reorder Code	Description	Qty/Case
000538*	8 1/2" tubing and adapter to connect leg bag to male external or Foley catheters	50/Case
6A4194*	18" Latex Extension tubing with connector, Sterile	

***CAUTION:** Components of this product contain natural rubber latex which may cause allergic reactions.

Bardia®

DRAINAGE BAGS

- Available with either swivel hanger or hook & loop cloth strap which facilitate attachment to bed frame or wheelchair
- 2000ml capacity with anti-reflux chamber
- Urine sampling area
- Easy-to-drain outlet tube with clamp and tube holder
- Plastic bag, single-use, sterile fluid path
- Not made with natural rubber latex



Reorder Code	Hanger Type	Qty/Case
802001	Swivel hanger	20/Case
802002	Hook and loop cloth strap	



Swivel hanger



Hook and loop cloth strap





Foley Catheters and Trays

Summary of Indications and Warnings

FOLEY CATHETERS

Foley catheters are intended for use in the drainage and/or collection and/or measurement of urine.

Please consult product label and insert for any indications, contraindications, hazards, warnings, cautions and directions for use.

Coated 2-Way Latex Foley Catheters

- Silicone-elastomer coating
- Available with either a 5cc or 30cc balloon
- Peel-pouch

Reorder Code	French Size	Balloon Size	Qty/Case
123512A	12 Fr.	5cc	12/Case
123514A	14 Fr.	5cc	
123516A	16 Fr.	5cc	
123518A	18 Fr.	5cc	
123520A	20 Fr.	5cc	
123522A	22 Fr.	5cc	
123524A	24 Fr.	5cc	
123526A	26 Fr.	5cc	
123528A	28 Fr.	5cc	
123530A	30 Fr.	5cc	
123612A	12 Fr.	30cc	
123614A	14 Fr.	30cc	
123616A	16 Fr.	30cc	
123618A	18 Fr.	30cc	
123620A	20 Fr.	30cc	
123622A	22 Fr.	30cc	
123624A	24 Fr.	30cc	
123626A	26 Fr.	30cc	
123628A	28 Fr.	30cc	
123630A	30 Fr.	30cc	

CAUTION: This product contains natural rubber latex which may cause allergic reactions.



2-Way All-Silicone Foley Catheters

- Available with either a 5cc or 30cc balloon
- Peel-pouch
- Not made with natural rubber latex

Reorder Code	French Size	Balloon Size	Qty/Case
806508	8 Fr.	3cc	12/Case
806510	10 Fr.	3cc	
806512	12 Fr.	5cc	
806514	14 Fr.	5cc	
806516	16 Fr.	5cc	
806518	18 Fr.	5cc	
806520	20 Fr.	5cc	
806522	22 Fr.	5cc	
806524	24 Fr.	5cc	
806316	16 Fr.	30cc	
806318	18 Fr.	30cc	
806320	20 Fr.	30cc	
806322	22 Fr.	30cc	
806324	24 Fr.	30cc	



Bardia®

ADD-A-FOLEY TRAY

Drainage bag and prepping components only:

- 2000ml center-entry urine drainage bag
- Waterproof underpad, fenestrated drape, 3 povidone iodine swabsticks, 5g lubricating jelly, 2 gloves not made with natural rubber latex, specimen container with label
- Pre-filled inflation syringe for 5cc and 30cc catheters
- CSR-wrap covered tray

Reorder Code	Description	Qty/Case
802015	For use with 5cc balloon catheters	10/Case
802035	For use with 30cc balloon catheters	



Bardia®

FOLEY INSERTION TRAYS

Foley catheter prepping components only:

- Waterproof underpad, fenestrated drape, 5g lubricating jelly, 2 gloves not made with natural rubber latex, specimen container with label and/or graduated collection container
- For 5cc or 30cc balloon Foley catheters
- 3 povidone iodine (PVI) or 3 skin cleansing swab sticks
- CSR-wrap package or Peel-top package

Reorder Code	Description	Qty/Case
Peel-Top Package		20/Case
802010	Includes PVI swabs and 10cc syringe	
802030	Includes PVI swabs and 30cc syringe	
802110	Includes BZK swabs and 10cc syringe	
802130	Includes BZK swabs and 30cc syringe	
CSR Wrap		
802011	Includes PVI swabs and 10cc syringe	
802031	Includes PVI swabs and 30cc syringe	

Bardia®

CLOSED SYSTEM FOLEY TRAYS

- 5cc silicone-elastomer coated latex Foley catheter pre-connected to a 2000ml urine drainage bag
- Waterproof underpad, fenestrated drape, 3 povidone iodine swabsticks, pre-filled 10cc inflation syringe, 5g lubricating jelly, 2 gloves not made with natural rubber latex, specimen container with label
- CSR-wrap covered tray

Reorder Code	French Size	Qty/Case
802016	16 Fr.	10/Case
802018	18 Fr.	

CAUTION: This product contains natural rubber latex which may cause allergic reactions.



Piston Syringe

Reorder Code	Description	Qty/Case
750379	BARD® Enteral 60cc Piston syringe in polyethylene bag, non-sterile	30/Case
802060	Irrigation Systems, BARDIA®, Piston Irrigation Tray	20/case



Reorder Code	Page Number	HCPCS Code	Reorder Code	Page Number	HCPCS Code	Reorder Code	Page Number	HCPCS Code
000538	31	A4331	35104	27	A4349	423714	15	A4352
123512A	36	A4338	35105	27	A4349	423716	15	A4352
123514A	36	A4338	35301	27	A4349	423718	15	A4352
123516A	36	A4338	35302	27	A4349	4A4173	31	A4358
123518A	36	A4338	35303	27	A4349	4A4174	31	A4358
123520A	36	A4338	35304	27	A4349	4A4194	31	A4331
123522A	36	A4338	35305	27	A4349	4A2044	21	A4353
123524A	36	A4338	36101	29	A4349	4A2045	21	A4353
123526A	36	A4338	36102	29	A4349	4A2046	21	A4353
123528A	36	A4338	36103	29	A4349	4A2056	21	A4353
123530A	36	A4338	36104	29	A4349	4A3052	21	A4353
123612A	36	A4338	36105	29	A4349	4A3053	21	A4353
123614A	36	A4338	36301	29	A4349	4A3054	21	A4353
123616A	36	A4338	36302	29	A4349	4A3055	21	A4353
123618A	36	A4338	36303	29	A4349	4A5042	21	A4353
123620A	36	A4338	36304	29	A4349	4A5044	21	A4353
123622A	36	A4338	36305	29	A4349	4A5108	20	A4353
123624A	36	A4338	37101	27	A4349	4A5110	20	A4353
123626A	36	A4338	37102	27	A4349	4A5142	20	A4353
123628A	36	A4338	37103	27	A4349	4A5144	20	A4353
123630A	36	A4338	37104	27	A4349	4A5146	20	A4353
150106	31	A4358	37105	27	A4349	4A6044	21	A4353
150107	31	A4358	37301	27	A4349	4A6110	20	A4353
150108	31	A4358	37302	27	A4349	4A6142	20	A4353
150319	31	A4358	37303	27	A4349	4A6144	20	A4353
150332	31	A4358	37304	27	A4349	4A6146	20	A4353
150419	30	A4358	37305	27	A4349	4A7114	20	A4353
150432	30	A4358	38301	29	A4349	4A7116	20	A4353
150507	31	A5114	38302	29	A4349	4A8114	20	A4353
150619	31	A4335	38303	29	A4349	4A8116	20	A4353
150632	31	A4335	38304	29	A4349	4A7044	21	A4353
150719	30	A4358	38305	29	A4349	50610	11	A4352
150732	30	A4358	39101	27	A4349	50612	11	A4352
150819	30	A4358	39102	27	A4349	50614	11	A4352
150832	30	A4358	39103	27	A4349	50616	11	A4352
150919	30	A4335	39104	27	A4349	50618	11	A4352
150932	30	A4335	39105	27	A4349	50810G	7	A4352
151919	30	A4358	39301	27	A4349	50812G	7	A4352
151932	30	A4358	39302	27	A4349	50814G	7	A4352
162110	31	A5114	39303	27	A4349	50816G	7	A4352
162210	31	A5114	39304	27	A4349	50818G	7	A4352
162310	31	A5114	39305	27	A4349	51310	9	A4351
32101	28	A4349	420614	15	A4351	51312	9	A4351
32102	28	A4349	420706	14	A4351	51314	9	A4351
32103	28	A4349	420708	14	A4351	51316	9	A4351
32104	28	A4349	420710	14	A4351	51318	9	A4351
32105	28	A4349	420712	14	A4351	51320	9	A4351
32301	28	A4349	420714	14	A4351	51610	10	A4351
32302	28	A4349	420716	14	A4351	51612	10	A4351
32303	28	A4349	420718	14	A4351	51614	10	A4351
32304	28	A4349	421008	15	A4351	51616	10	A4351
32305	28	A4349	421010	15	A4351	51618	10	A4351
33101	28	A4349	421108	14	A4351	51620	10	A4351
33102	28	A4349	421110	14	A4351	51614S	11	A4353
33103	28	A4349	421610	15	A4351	51808	6	A4351
33104	28	A4349	421612	15	A4351	51810	6	A4351
33105	28	A4349	421614	15	A4351	51812	6	A4351
33301	28	A4349	421708	14	A4351	51814	6	A4351
33302	28	A4349	421710	14	A4351	51816	6	A4351
33303	28	A4349	421712	14	A4351	51818	6	A4351
33304	28	A4349	421714	14	A4351	52306	9	A4351
33305	28	A4349	421716	14	A4351	52308	9	A4351
35101	27	A4349	421718	14	A4351	52310	9	A4351
35102	27	A4349	423710	15	A4352	52606	10	A4351
35103	27	A4349	423712	15	A4352	52606G	11	A4351

Reorder Code	Page Number	HCPCS Code	Reorder Code	Page Number	HCPCS Code
52608	10	A4351	802010	38	A4310
52608G	11	A4351	802011	38	A4310
52610	10	A4351	802015	38	A4354
52610G	11	A4351	802016	39	A4314
52806G	8	A4351	802018	39	A4314
52808G	8	A4351	802030	38	A4310
52810G	8	A4351	802031	38	A4310
53310	9	A4351	802035	38	A4354
53312	9	A4351	802060	39	A4322
53314	9	A4351	802100	22	A4310
53316	9	A4351	802110	38	A4310
53318	9	A4351	802115	22	A4353
53320	9	A4351	802130	38	A4310
53610	10	A4351	802214	22	A4353
53610G	11	A4351	802408	16	A4351
53612	10	A4351	802410	16	A4351
53612G	11	A4351	802412	16	A4351
53614	10	A4351	802414	16	A4351
53614G	11	A4351	802416	16	A4351
53614GS	11	A4353	802418	16	A4351
53616	10	A4351	802420	16	A4351
53616G	11	A4351	802422	16	A4351
53616GS	11	A4353	802512	16	A4352
53618	10	A4351	802514	16	A4352
53618G	11	A4351	802516	16	A4352
53620	10	A4351	802518	16	A4352
53620G	11	A4351	802520	16	A4352
53810G	7	A4351	802522	16	A4352
53812G	7	A4351	806316	37	A4344
53814G	7	A4351	806318	37	A4344
53816G	7	A4351	806320	37	A4344
53818G	7	A4351	806322	37	A4344
53820G	7	A4351	806324	37	A4344
61310	12	A4351	806508	37	A4344
61312	12	A4351	806510	37	A4344
61314	12	A4351	806512	37	A4344
61316	12	A4351	806514	37	A4344
61318	12	A4351	806516	37	A4344
61610	13	A4351	806518	37	A4344
61612	13	A4351	806520	37	A4344
61614	13	A4351	806522	37	A4344
61616	13	A4351	806524	37	A4344
61618	13	A4351	PWFX30	26	N/A
62308	12	A4351	PW100	26	N/A
62310	12	A4351	PW200	26	N/A
62608	13	A4351	PWKIT03	26	N/A
62610	13	A4351			
63310	12	A4351			
63312	12	A4351			
63314	12	A4351			
63316	12	A4351			
63318	12	A4351			
63610	13	A4351			
63612	13	A4351			
63614	13	A4351			
63616	13	A4351			
63618	13	A4351			
650819	30	A4335			
650832	30	A4335			
651919	30	A4335			
651932	30	A4335			
6A4194	31	A4331			
750379	39	A4322			
802001	32	A4357			
802002	32	A4357			

*NOTE: Reimbursement policies vary widely from insurer and reimbursement policies of the same insurer may vary in different sections of the U.S. As reimbursement policies are subject to change, the information contained herein may not be accurate at the time of use. Prior to submission of a claim for reimbursement, the user should contact the insurer (i.e., Medicare, Medicaid, or private payer) to verify applicable codes and reimbursement levels. BD does not guarantee that the purchases of any similar item is reimbursable in whole or in part, by any public or private health care insurer. Source for HCPCS codes: <http://www.hcfa.gov/statspufiles.htm>

Care. Comfort. Choice.™

For more information on any of our products, contact
Bard Care at 800.526.4455 or visit www.bardcare.com.

8195 Industrial Blvd
Covington, GA 30014
800 526.4455
www.bardcare.com

BD, the BD Logo and Bard, Bardia, Care, Comfort, Choice., Clean-Cath, CleanGlide, Dispoz-a-Bag, Flip-Flo, Magic³, Magic³ Go, Magic³ Touchless, Natural, Personal Catheter, Pop-On, Spirit, StatLock, Sure-Grip, Touchless, UltraFlex and WideBand are trademarks of Becton, Dickinson and Company and its affiliates. © 2019 BD. All rights reserved.
2002-05 BD-XXXXX

BARD | **CARE**
has joined BD

